

FEB 14 2000

Page 1 of 2

<genX> international, Inc. 510(K) Premarket Notification

510(K) Summary

K993881

a) **DEVICE NAME**

Proprietary Name: <genX> Culture dish

Classification Name: Assisted Reproduction Labware

b) **Submitted by:**

<genX> international, inc.

170 Fort Path Rd.

Madison, CT 06443

ESTABLISHMENT REGISTRATION No.: 9003605

Telephone number: 203-245-4901

Fax number: 203-245-4994

E-mail: genXintl@aol.com

Contact individual: Michael D. Cecchi

President

c) **CLASSIFICATION: Class II**

Assisted Reproductive Labware

Product: <genX> Culture dish

Procode: 85 MQK

CFR#: 884.6160

d) **PERFORMANCE STANDARDS**

Performance Standards under Section 514 of the ACT have not been developed for this device. However, Special controls have been identified by the FDA to provide reasonable assurance of safety and effectiveness of the device in assisted reproductive procedures.

e) **PROPOSED LABELS, LABELING, AND ADVERTISING**

The proposed labeling and instruction material is included in this package in the appropriate section.

f) **DESCRIPTION AND INTENDED USE**

The <genX> Culture Dish is intended to be used when culturing tissues and cells. It may be used for the use with sperm and for the culturing of embryos.

The catalog number is GDEC-100

The basic design is that of a traditional culture dish. Approx 60mm round with ½ inch sidewalls. This dish does have defined areas for the cell to be confined within. These areas are squares of picket-fenced walls, which will allow fluids, nutrients and byproducts to freely flow among the culturing cell, embryos.

g) **SUBSTANTIALLY EQUIVALENT**

The <genX> Culture Dish is substantially equivalent to several products currently marketed in the US. See the appropriate section for the comparisons. (Section 2)

1. Same intended use
2. Same design
3. Same materials configured the same for use.
4. Perform similar procedures, by method.
5. Same performance criteria
6. Same sterility methods and testing
7. Same intended population

h) **STERILIZATION PROCEDURES AND FACILITIES**

All dishes manufactured, will be Gamma Irradiated by and according to the guidelines in place.

Dosage: 2.5 Mrad

Sterility Assurance Level (SAL): of (10 to the -6)

After sterilization the contractor will issue a "Certificate of Sterilization".

Validation assurance is in accordance with ANSI / AAMI / ISO 11137-1994 Standards.

See Section 11.0



FEB 14 2000

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850Mr. Michael D. Cecchi
President
<genX> international, Inc.
170 Fort Path Road
Madison, CT 06443Re: K993881
<genX> Tissue Culture Dish
Dated: November 15, 1999
Received: November 16, 1999
Regulatory Class: II
21 CFR §884.6110/Procode: 85 MQK

Dear Mr. Cecchi:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4613. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597, or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

CAPT Daniel G. Schultz, M.D.
Acting Director, Division of Reproductive,
Abdominal, Ear, Nose and Throat,
and Radiological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATIONS FOR USE

510 (k) Number (if known) K99 3881

Device Names: <genX> Tissue Culture Dish

Indication for Use:


"To culture tissues and cells"

May be used to hold sperm and culture embryos

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X or Over-the Counter Use _____


(Division Sign-Off)
Division of Reproductive, Abdominal, ENT,
and Radiological Devices
510(k) Number K993881