

**510(k) Summary**  
**Special 510K**

APR 25 2014

**Precision Medical, Inc. Easy Go Aspirator.**

**Submitter Information**

|                         |  |
|-------------------------|--|
| Submitter               | Precision Medical, Inc.<br>300 Held Drive<br>Northampton, Pa.<br>18067 |
| Facility Registration # | 2523148  |
| Contact                 | James Parker<br>Quality Assurance Manager                              |
| Tel:                    | (610)-262-6090 Extensions 228  |
| Fax:                    | (610)-262-6080   |
| Preparation Date:       | January 20, 2014   |

**Device Name**

|                              |                           |
|------------------------------|---------------------------|
| Proprietary Name:            | Easy Go Vac               |
| Model number                 | PM66                      |
| Common Name:                 | Aspirator                 |
| Classification Name:         | Pump, Portable,Aspiration |
| Classification Product Code: | BTA                       |
| Regulation number:           | 878.4780                  |
| Regulatory Status            | Class II                  |

**Predicate Device Equivalence**

Precision Medical, Inc. is claiming substantial equivalence to the Precision Medical Inc  
 Portable aspirator  
 Model # pm 65  
 510Knumber K971749  
 Manufactured by:  
 Precision Medical Inc.  
 300 Held Drive  
 Northampton Pa 18067

**Device Descriptions ;**

Predicate Device (PM65)

Small DC motor that powers a suction pump monitored by a gauge, vacuum is set by the care provider up to 21 inches of mercury. Uses a internal 12 volt lead acid battery, with external charger/power supply.

## New Device PM66 (3 models)

### **General description;**

Small DC motor that powers a suction pump monitored by a gauge, vacuum is set by the care provider up to 21 inches of mercury. Uses an internal battery, with internal charger/power supply, on two models, and AC only on the third model.

### **PM66LI model**

Small DC motor that powers a suction pump monitored by a gauge, vacuum is set by the care provider up to 21 inches of mercury. Uses an internal Lithium battery, with internal charger/power supply. (charge status indicator along with a battery level indicator, disposable or reusable waste container)

### **PM66AC**

Small DC motor that powers a suction pump monitored by a gauge, vacuum is set by the care provider up to 21 inches of mercury. Uses AC 100-240 volts power only.

### **PM66S**

Small DC motor that powers a suction pump monitored by a gauge, vacuum is set by the care provider up to 21 inches of mercury. Uses an internal lead acid battery, with internal charger/power supply. (charge status indicator along with a battery level indicator, disposable or reusable waste container)

### Differences from the new devices to the Predicate device

- Battery charger is internal on new devices for PM66LI and PM66S
- Lithium battery on the PM66LI model
- AC only Model
- Addition of the charge indicator
- Addition of the battery level indicator
- Choice of reusable or disposable waste container

### **Conclusion:**

The new devices are adding features that enhance the operation,

The lithium battery extends battery run time.

The AC model simplifies the device, for operators that do not need a portable device but want a simple small, easy to use device.

The lead acid battery model is the same as the predicate, with the added features of optional containers, a battery status indicator and the charge indicator.

The safety and effectiveness of the new devices compared to the predicate are substantially equivalent, with added features that ease, and enhance the operation.

**Intended Use**

The Easy Go Aspirator provides a portable, AC/DC powered medical vacuum source. It is intended for use in the homecare / health care environments.

**Applied Standards;**

**ISO 15223-1:2012**

**BS EN 1041:2008**

**EN ISO 10079-1:2009**

**IEC 60601-1:2005 ed3.1 Consol. with am1**

**IEC 60601-1-2:2007**

**IEC 60601-1-6:2007**

**IEC 60601-1-11:2010**

**EN ISO 14971:2012**

**ISO 13485:2012**



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

April 25, 2014

Precision Medical Incorporated  
Mr. James Parker  
Quality Assurance Manager  
300 Held Drive  
Northampton, Pennsylvania 18067

Re: K140179  
Trade/Device Name: Easy Go Vac  
Regulation Number: 21 CFR 878.4780  
Regulation Name: Powered suction pump  
Regulatory Class: Class II  
Product Code: BTA  
Dated: April 11, 2014  
Received: April 14, 2014

Dear Mr. Parker:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

Page 2 – Mr. James Parker

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**David Krause -S**

for Binita S. Ashar, M.D., M.B.A., F.A.C.S.  
Acting Director  
Division of Surgical Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indications for Use**

510(k) Number (if known)  
K140179

Device Name  
Precision Medical Easy Go Vac (PM66)

Indications for Use (Describe)

The Easy Go Vac Aspirator provides a portable, AC/DC powered medical vacuum source. It is intended for use in the homecare / healthcare environments.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.**

**FOR FDA USE ONLY**

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

**Joshua C. Nipper -S**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRASStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*