

AUG 01 2002

Section 5-1

K021503



Bisco, Inc.
1100 W. Irving Park Road, Schaumburg, IL. 60193 U.S.A.
Telephone: (847) 534-6000 or 1-800-BIS-DENT Fax: (847) 534-6396
WEB SITE <http://www.bisco.com>

Bringing Science to the Art of Dentistry™

Contact: Stephen D. Smith

SUMMARY

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 190 and 21 CFR par 807.92

Trade Name: **UNINHIBITED**
Common Name: **Surface Sealant**
Classification name: Material, Tooth Shade, Resin
Class II per 21 CFR 872.3690

Description of Applicant Device:

UNINHIBITED contains Hexafunctional Urethane Acrylate (30.6%), Polyethylene Glycol (400) Diacrylate (61.4%), phosphine oxide, diphenyl (Photoinitiator) 8.0% 4-methyl methacrylate (Inhibitor) 500 ppm. It is supplied in a bottle. It is applied with a brush to seal the surface of a composite restoration.

Intended uses of Applicant Device:

The principle use of **UNINHIBITED** is to seal the surface of:
Direct Composite Restorations, Indirect Composite Restorations, Provisional Restorations, Amalgam Restorations, Enamel Surfaces around Orthodontic Brackets, Processed Acrylic Prostheses-Removable Partial & Complete Dentures and Maxillofacial Prostheses, Processed Acrylic Appliances-Orthodontic Retainers, Lingual Retainers and Extracoronary Splints

Predicate Devices: Fortify K882915 cleared by FDA on 08/12/88

Significant Performance Characteristics:

	UNINHIBITED	FORTIFY
Intended Use	Resin sealant.	Resin sealant.
Product Description	Clear solution.	Clear solution.
Delivery System	Brush	Brush

Side by side comparisons of **UNINHIBITED** to the predicate device **FORTIFY** clearly demonstrates that the applicant device is substantially equivalent to the legally marketed devices. The ingredients of **UNINHIBITED** were tested for biocompatibility and were found to be non-toxic.

It is concluded that the information supplied in this submission has proven the safety and efficacy of **UNINHIBITED**.



AUG 01 2002

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Stephen Smith
Manager of Regulatory Affairs
Bisco, Incorporated
1100 West Irving Park Road
Schaumburg, Illinois 60193

Re: K021503
Trade/Device Name: Uninhibited
Regulation Number: 872.3310
Regulation Name: Coating Material for Resin Fillings
Regulatory Class: II
Product Code: EBD
Dated: May 07, 2002
Received: May 09, 2002

Dear Mr. Smith:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.


Page 2 – Mr. Smith

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4613. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Timothy A. Ulatowski
Director

Division of Anesthesiology, General Hospital,
Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K021503

Device Name: UNINHIBITED

Indications For Use:

TO SEAL THE SURFACE OF:

**DIRECT COMPOSITE RESTORATIONS
INDIRECT COMPOSITE RESTORATIONS
AMALGAM RESTORATIONS
INTRACORONAL SPLINTS
EXTRACORONAL SPLINTS
LINGUAL RETAINERS
ENAMEL SURFACES AROUND ORTHODONTIC BRACKETS
PROVISIONAL RESTORATIONS
PROCESSED ACRYLIC PROSTHESES-REMOVABLE PARTIAL & COMPLETE DENTURES
AND MAXILLOFACIAL PROSTHESES
PROCESSED ACRYLIC APPLIANCES-ORTHODONTIC RETAINERS**

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Optional Format 3-10-98)

Susan Purpura

(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number K021503