

Food and Drug Administration 2098 Gaither Road Rockville MD 20850

## JUL 1 6 2003

Ms. Cynthia C. Knapp Director Lab Services TREK Diagnostic System 982 Keynote Circle, Suite 6 Cleveland, Ohio 44131

Re:

k031545

Trade/Device Name: ESBL Confirmatory Test Plate

Ceftazimdime, Ceftazidime/clavulanic acid, Cefotaxime,

Cefotaxime/ clavulanic acid

Regulation Number: 21 CFR 866.1640

Regulation Name: Antimicrobial Susceptibility Test

Regulatory Class: Class II Product Code: JWY Dated: May 16, 2003 Received: May 20, 2003

Dear Ms. Knapp:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

## Page 2 -

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <a href="http://www.fda.gov/cdrh/dsma/dsmamain.html">http://www.fda.gov/cdrh/dsma/dsmamain.html</a>.

Sincerely yours,

Steven I. Gutman, M.D., M.B.A.

Director

Office of In Vitro Diagnostic Device

Steven Butman

**Evaluation and Safety** 

Center for Devices and

Radiological Health

Enclosure

510 (k) Number (If known):

K031545

Device Name: ESBL Confirmatory Test Plate

Indications For Use:

The Sensititre ESBL confirmatory test plate is an *in vitro* diagnostic product for detection of ESBL's in clinical isolates of *Klebsiella pneumoniae*, *Klebsiella oxytoca* and *Escherichia coli*. Test results can be read either manually or automatically on a Sensititre® ARIS or Autoreader.

This 510(K) is for ESBL Confirmatory Test Plate for <u>Klebsiella pneumoniae</u>, <u>Klebsiella oxytoca</u>, and <u>Escherichia coli</u> containing the following antimicrobials as per the NCCLS reference standard (M7):

Ceftazimdime: 0.06-128

Ceftazidime/clavulanic acid: 0.06/4-128/4

Cefotaxime: 0.06-64

Cefotaxime/clavulanic acid: 0.06/4-64/4

(PLEASE DO NOT WE	RITE BELOW	THIS LINE- CONTINU	E ON ANOTHER	(PAGE IF NEEDED)	)
	Concurrence of CDRH, Office of Device Evaluation (ODE)  Division Sign-Off  Office of In Vitro Diagnostic Device Evaluation and Safety				
	510(k)	K031242	·	,	
Prescription Use_ (Per 21 CFR 801.109	<b></b>	OR	Over-The	e-Counter Use	