OCT 1 4 2003

Section 4. 510(k) Summary

| General Provisions | Submitter's Name and Address | EKOS Corporation 22030 20 th Ave. SE Suite 101 Bothell, WA 98021 |
|--|--|--|
| | Contact Person | Jocelyn Kersten 425-482-1108 425-482-1109 (fax) jkersten@EKOSCORP.com |
| | Classification Name | Catheter, Continuous Flush (KRA) |
| | Common or Usual Name | Continuous Flush Catheter |
| | Proprietary Name | EKOS Peripheral Infusion System |
| Name of Predicate Device | Predicate Device EKOS Peripheral Infusion System | 510(k) Reference No. K030637 |
| Device Description | The system consists of a disposable infusion/ultrasound catheter and an instrument that generates and controls the delivery of energy to the catheter. The catheter contains a single ultrasound transducer, located at the distal tip, a thermal sensor and a distal end hole for placement over a guide wire and fluid infusion. | |
| Intended Use | The EKOS Peripheral Infusion System is intended for the controlled and selective infusion of physician-specified fluids, including thrombolytics, into the peripheral vasculature. | |
| Summary of Technological Characteristics | The proposed EKOS Peripheral Infusion System is similar in construction and materials to the previously cleared EKOS Peripheral Infusion System. | |
| Test Summary | The proposed EKOS Peripheral Infusion System is considered to be substantially equivalent to the currently marketed EKOS Peripheral Infusion System based on a comparison of the intended uses and designs and results of the testing and evaluations performed. | |
| | | |



OCT 1 4 2003

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

EKOS Corp. c/o Ms. Jocelyn Kersten 22030 20th Ave SE, Suite 101 Bothell, WA 98021

Re: K033214

EKOS Peripheral Infusion System Regulation Number: 21 CFR 870.1210

Regulation Name: Catheter, Continuous Flush

Regulatory Class: Class II

Product Code: KRA Dated: October 2, 2003 Received: October 3, 2003

Dear Ms. Kersten:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act

Page 2 – Ms. Jocelyn Kersten

or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4648. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address

http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours,

Bram D/Zuckerman, l

Director

Division of Cardiovascular Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indication for Use Section 3. 510(k) Number (033214 **Device Name** The EKOS Peripheral Infusion System is intended for the controlled and **Indications for** selective infusion of physician-specified fluids, including thrombolytics, Use into the peripheral vasculature. (PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED) Concurrence of CDRH, Office of Device Evaluation (ODE) Prescription Use_ OR Over The Counter Use____ (Per 21 CFR 801.109) (Optional Format 1-2-96) (Division Sign-Off) Division of Cardiovascular Devices 510(k) Number