

JUL - 8 2005

K051530
1/2

2.6

510(k) Summary

Submitter: American Medical Systems
10700 Bren Road West
Minnetonka, MN 55343
Phone: 952-933-4666
Fax: 952-930-6496

Contact Person: Denise Thompson

Date Summary Prepared: June 8, 2005

Device Common Name: Urethral Sling, Surgical Mesh

Device Trade Name: Monarc™, Monarc™ +, and Monarc™ C
Subfascial Hammocks / BioArc TO™,
BioArc™ TO +, and BioArc TO - C
Subfascial Hammocks

Device Classification Name: Surgical Mesh, polymeric

Predicate Device: Monarc™ Subfascial Hammock, K023516
BioArc™ TO Subfascial Hammock, K040538

Device Description:

The Monarc and BioArc TO Subfascial Hammocks are suburethral sling procedure that uses a transobturator surgical approach to treat stress urinary incontinence. They are sterile, single use procedure kits consisting of two stainless steel curved needle passers and a mesh or mesh and graft sling assembly.

Indications for Use:

The Monarc and BioArc TO Subfascial Hammocks are intended for the placement of a suburethral sling for the treatment of female stress urinary incontinence (SUI) resulting from urethral hypermobility and / or intrinsic sphincter deficiency.

Comparison to Predicate Device:

The Monarc and BioArc TO + and C needle passers offer physicians alternative needle options to place the suburethral sling. The needle passers are all designed for a transobturator approach.

The Indications for Use, fundamental scientific technology, surgical approach, sling placement, and materials are all the same as the predicates.

Supporting Information:

The risk analysis and the verification / validation activities reported in this Special 510(k) application substantiate equivalence to the predicate devices and did not raise any new questions of safety or efficacy.

DNC
7/8/05

K051530

Conclusion:

The Monarc and BioArc TO + and C Subfascial Hammock versions are substantially equivalent to their predicates with respect to intended use, technological characteristics, and performance.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room –WO66-G609
Silver Spring, MD 20993-0002

Ms. Denise Thompson
Regulatory Specialist
American Medical Systems
10700 Bren Road West
MINNETONKA MN 55343

SEP 28 2012

Re: K051530
Trade/Device Name: Monarc, Monarc +, and Monarc C Subfascial Hammocks
and the BioArc TO, BioArc To +, and BioArc TO-C
Subfascial Hammocks
Regulation Number: 21 CFR 878.3300
Regulation Name: Surgical mesh
Regulatory Class: II
Product Code: OTN
Dated: June 8, 2005
Received: June 9, 2005

Dear Ms. Thompson:

This letter corrects our substantially equivalent letter of July 8, 2005.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must

comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Benjamin R. Fisher". The signature is fluid and cursive, with a large initial "B" and "F".

Benjamin R. Fisher, Ph.D.
Director
Division of Reproductive, Gastro-Renal,
and Urological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

