

MAR 20 2007

K063556

1/1

### 510(k) Summary of Safety & Effectiveness

Submitter: PM Devices Inc.  
13700 Mayfield Place, Unit 2135  
Richmond, British Columbia V6V 2E4  
Canada  
Phone: (604) 270-4344  
Fax: (604) 270-4384

Contact Person: Tim Verspagen

Device Name: PeriPatch Aegis

Common Name: Surgical Mesh

Classification Name: Surgical Mesh (Product Code: FTM)

Predicate Devices: 1) PeriPatch Sleeve (K033985)  
2) Peri-Strips Staple Line Reinforcement - Sleeve Configuration (K992537)

Device Description: The PeriPatch Aegis is intended be used as a surgical patch material to reinforce surgical staples during stapling procedures. The PeriPatch Aegis device assembly consists of a strip of buttress material, a cartridge which mates the device to the endoscopic stapler, and finally a loader which assists the user to place the device on to the stapler.

Intended Use: The PeriPatch Endo-Sleeve is intended for reinforcing staples and suture lines during a number of resection techniques - surgical repair of tissue deficiencies such as:

- Lung resections
- Biopses
- Lobectomies
- Gastric banding
- Abdominal and thoracic wall repair

Technological Comparison & Performance Data: Tissue functional testing was carried out in accordance with the *Guidance for the Preparation of a Premarket Notification Application for a Surgical Mesh*. The test results demonstrate suitable material properties for the indications of the device.

The PeriPatch Aegis is substantially equivalent to the predicate devices as it has similar indications for use, similar materials and similar technological characteristics.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

PM Devices, Inc.  
% Mr. Tim Verspagen  
Regulatory Affairs Director  
13700 Mayfield Place, Unit 2135  
Richmond, British Columbia  
V6V 2E4 Canada

Re: K063556

MAR 20 2007

Trade/Device Name: PeriPatch Aegis  
Regulation Number: 21 CFR 878.3300  
Regulation Name: Surgical mesh  
Regulatory Class: II  
Product Code: FTM  
Dated: February 23, 2007  
Received: February 26, 2007

Dear Mr. Verspagen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set


Page 2 – Mr. Tim Verspagen

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a stylized flourish at the end.

Mark N. Melkerson

Director

Division of General, Restorative  
and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

## Statement of Indications for Use

---

510(k) Number (if known): K063556

Device Name: PeriPatch Aegis

### Indications for Use:

The PeriPatch Aegis is intended for reinforcing staples and suture lines during a number of resection techniques - surgical repair of tissue deficiencies such as:


- Lung resections
- Biopses
- Lobectomies
- Gastric banding
- Abdominal and thoracic wall repair

Prescription Use   X   AND/OR Over-The-Counter Use             
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER  
PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)



**(Division Sign-Off)**  
**Division of General, Restorative,**  
**and Neurological Devices**

510(k) Number   K063556