

3.0 510(k) Summary**K083025**Page 1 of 1

Sponsor: Synthes (USA)
1301 Goshen Parkway
West Chester, PA 19380
(610) 719-6940

Contact: Sheri L. Musgnung
Synthes (USA)
1301 Goshen Parkway
West Chester, PA 19380
(610) 719-6940
FAX (610) 484-356-9682

JAN 14 2009

Device Name: Synthes 4.5 mm VA-LCP Curved Condylar System

Classification: Class II, §888.3030 – Single/multiple component metallic bone fixation appliances and accessories

Predicate Device: Synthes 4.5 mm Curved Broad LCP Plates
Synthes Locking Condylar Plates
Synthes Condylar Buttress Plate – DCP
Synthes 4.5 mm Broad DCP Plates

Device Description: The Synthes 4.5 mm VA-LCP Curved Condylar System consists of plates and screws that offer a screw to plate variable angle locking construct. The 4.5 mm VA-LCP Curved Condylar System enhances fracture fixation by providing the flexibility to lock screws in trajectories that can diverge from the central axis of the plate hole.

Intended Use: The Synthes 4.5 mm VA-LCP Curved Condylar System is intended for buttressing multi-fragmentary distal femur fractures, including supracondylar, intra-articular and extra-articular condylar fractures, periprosthetic fractures, fractures in normal or osteopenic bone, and non-unions and malunions.

Substantial Equivalence: Information presented supports substantial equivalence.



JAN 14 2009

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Synthes® (USA)
% Ms. Sheri L. Musgnung
Regulatory Affairs Manager
1301 Goshen Parkway
West Chester, Pennsylvania 19380

Re: K083025

Trade/Device Name: Synthes® 4.5mm VA-LCP Curved Condylar System
Regulation Number: 21 CFR 888.3030
Regulation Name: Single/multiple component metallic bone fixation appliances and accessories
Regulatory Class: II
Product Code: HRS, JDP, HWC
Dated: December 10, 2008
Received: December 11, 2008

Dear Ms. Musgnung:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Sheri L. Musgnung

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at toll-free number (800) 638-2041 or (240) 276-3150 or the Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure



2.0

Indications for Use

510(k) Number (if known):

K083025

Device Name:

Synthes 4.5 mm VA-LCP Curved Condylar System

Indications for Use:

The Synthes 4.5 mm VA-LCP Curved Condylar System is intended for buttressing multi-fragmentary distal femur fractures including: supracondylar, intra-articular and extra-articular condylar fractures, periprosthetic fractures, fractures in normal or osteopenic bone, and non-unions and malunions.

Prescription Use X
(Per 21 CFR 801.109)

AND/OR

Over-The-Counter Use
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of General, Restorative,
and Neurological Devices

510(k) Number K083025

000004