



SEP 24 2009

K092059  
Page 1 of 2

Page 1-11

## 510(k) Summary

### Submitter Information

Submitter's Name: Riverpoint Medical  
Address: 825 NE 25<sup>th</sup> Ave.  
Portland, OR 97232  
Phone Number: (503) 517-8001 or 866 445-4923  
Fax Number: (503) 517-8002  
Registration Number: 3006981798  
Contact Person: Douglas Rowley  
(503) 517-8001  
Date of Preparation: July 1st, 2009

### Device Name

Trade Name: RP Cutting Needle  
Common Name: Biopsy Needle (Disposable)  
Classification Name: Instrument, Biopsy

### Device Classification

FDA Class: 2  
Product Classification: 876.1075(a) Gastroenterology-Urology Biopsy Instrument, Needle  
Code: KNW  
Classification Panel: Gastroenterology-Urology

### Predicate Devices

*C.R. Bard®:* Biopsy-Cut® Needle; Magnum® Needle  
*Remington Medical:* MLL Needle; M-Needle  
*Cook® Medical:* Cook Biopsy Needle

510(k) Summary - RP Cutting Needle



K092059  
Page 2 of 2

### Special Controls

Not applicable for this device

### Device Description

The RP Cutting Needle is a sterile, disposable biopsy instrument composed of a stainless steel cannula with a molded plastic hub, and a stainless steel stylet with a notch removed for specimen collection and a molded plastic hub. The cannula and stylet hubs are provided with a spacer attached which is designed to allow for easy insertion of both hubs into a Bard® Magnum® reusable biopsy instrument. The needle is covered with a LDPE sheath for safety. Needles are available in 14-20 gauge, with lengths between 10 and 25 cm. The predicate devices listed contain identical or substantially equivalent materials and characteristics as described above.

### Intended Use

The RP Cutting Needle is intended to be used by medical professionals with the Bard® Magnum® reusable biopsy instrument for biopsies of soft tissues, such as the liver, lung, kidney or prostate. The RP Cutting Needle is not intended for use in bone. The RP Cutting Needle is provided sterile as a single use device.

### Safety and Effectiveness

The RP Cutting needle is designed and manufactured to be substantially equivalent to the predicate devices listed below for safety and effectiveness. All materials used were selected based on known biocompatibility and established histories of use in the medical device industry for non-implantable devices, and are identical or substantially equivalent to the materials used in the predicate devices listed. The RP Cutting Needle has been designed to have the same general shape, size and method of function as the predicate devices. Dimensional, functional and quality verifications, along with material certification has confirmed that the RP Cutting Needle is substantially equivalent to the predicate devices in size, shape, and function.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room – WO66-G609  
Silver Spring, MD 20993-0002

SEP 24 2009

Mr. Doug Rowley  
RA/QA Manager  
Riverpoint Medical  
825 NE 25<sup>th</sup> Avenue  
PORTLAND OR 97232

Re: K092059

Trade/Device Name: RP Cutting Needle  
Regulation Number: 21 CFR 876.1075  
Regulation Name: Gastroenterology-urology biopsy instrument  
Regulatory Class: II  
Product Code: KNW  
Dated: July 1, 2009  
Received: July 7, 2009

Dear Mr. Rowley:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

Page 2 –

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.


If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Janine M. Morris  
Acting Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure



### Indications for Use

510(k) Number: K092059

Device Name: RP Cutting Needle

**Indications for Use:**

The RP Cutting Needle is intended to be used by medical professionals with the Bard® Magnum® reusable biopsy instrument for biopsies of soft tissues, such as the liver, lung, kidney or prostate. The RP Cutting Needle is not intended for use in bone. The RP Cutting Needle is provided sterile as a single use device.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1

(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number K092059