

SPECIAL 510(k): Device Modification  
ODE Review Memorandum (Decision Making Document is Attached)

To: THE FILE

RE: DOCUMENT NUMBER k081715

---

This 510(k) submission contains information/data on modifications made to the SUBMITTER'S own Class II, Class III or Class I devices requiring 510(k). The following items are present and acceptable (delete/add items as necessary):

1. The name and 510(k) number of the SUBMITTER'S previously cleared device. (For a preamendments device, a statement to this effect has been provided.) k072784
2. Submitter's statement that the **INDICATION/INTENDED USE** of the modified device as described in its labeling **HAS NOT CHANGED** along with the proposed labeling which includes instructions for use, package labeling, and, if available, advertisements or promotional materials (labeling changes are permitted as long as they do not affect the intended use).
3. A description of the device **MODIFICATION(S)**, including clearly labeled diagrams, engineering drawings, photographs, user's and/or service manuals in sufficient detail to demonstrate that the **FUNDAMENTAL SCIENTIFIC TECHNOLOGY** of the modified device **has not changed**.  
**This change was for adding code number selection feature for the test strips that differs from predicate device with no code function for the test strips. In comparison to the predicate, there are minor changes in the physical appearance of the meter.**
4. **Comparison Information** (similarities and differences) to applicant's legally marketed predicate device including, labeling, intended use, physical characteristics, reagents, test steps, and performance characteristics.  
**The only differences are, adding code number selection feature, and physical appearance of size and weight of the device.**
5. A **Design Control Activities Summary** which includes:
  - a) Identification of Risk Analysis method(s) used to assess the impact of the modification on the device and its components, and the results of the analysis
  - b) Based on the Risk Analysis, an identification of the verification and/or validation activities required, including methods or tests used and acceptance criteria to be applied
  - c) A declaration of conformity with design controls. The declaration of conformity included:
    - i) A statement signed by the individual responsible, that, as required by the risk analysis, all verification and validation activities were performed by the designated individual(s) and the results demonstrated that the predetermined acceptance criteria were met, and
    - ii) A statement signed by the individual responsible, that the manufacturing facility is in conformance with design control procedure requirements as specified in 21 CFR 820.30 and the records are available for review.
6. A **Truthful and Accurate Statement, a 510(k) Summary or Statement and the Indications for Use Enclosure (and Class III Summary for Class III devices).**

The labeling for this modified subject device has been reviewed to verify that the indication/intended use for the device is unaffected by the modification. In addition, the submitter's description of the particular modification(s) and the comparative information between the modified and unmodified devices demonstrate that the fundamental scientific technology has not changed. The submitter has provided the design control information as specified in The New 510(k) Paradigm and on this basis, I recommend the device be determined substantially equivalent to the previously cleared (or their preamendment) device.

---

(Reviewer's Signature)

(Date)

Comments:

See the attached email correspondence with the requests to the sponsor for changes in the labeling, making control solutions available to the user, and addition of a statement to the indications for use statement on speaker function.

revised:8/1/03

### "SUBSTANTIAL EQUIVALENCE" (SE) DECISION MAKING DOCUMENTATION

	Yes	No	
1. Same Indication Statement?	X		If <b>YES</b> = Go To 3
2. Do Differences Alter The Effect Or Raise New Issues of Safety Or Effectiveness?		X	If <b>YES</b> = Stop <b>NSE</b>
3. Same Technological Characteristics?	X		If <b>YES</b> = Go To 5
4. Could The New Characteristics Affect Safety Or Effectiveness?			If <b>YES</b> = Go To 6
5. Descriptive Characteristics Precise Enough?		X	If <b>NO</b> = Go To 8 If <b>YES</b> = Stop <b>SE</b>
6. New Types Of Safety Or Effectiveness Questions?			If <b>YES</b> = Stop <b>NSE</b>
7. Accepted Scientific Methods Exist?			If <b>NO</b> = Stop <b>NSE</b>
8. Performance Data Available?	X		If <b>NO</b> = Request Data
9. Data Demonstrate Equivalence?	X		Final Decision: SE

Note: See

[http://erom.fda.gov/eRoomReg/Files/CDRH3/CDRHPremarketNotification510kProgram/0\\_4148/FLOWCHART%20DECISION%20TREE%20.DOC](http://erom.fda.gov/eRoomReg/Files/CDRH3/CDRHPremarketNotification510kProgram/0_4148/FLOWCHART%20DECISION%20TREE%20.DOC) for Flowchart to assist in decision-making process. Please complete the following table and answer the corresponding questions. "Yes" responses to questions 2, 4, 6, and 9, and every "no" response requires an explanation.

1. Explain how the new indication differs from the predicate device's indication:
2. Explain why there is or is not a new effect or safety or effectiveness issue: *Adding coding feature does not bring new safety concerns*
3. Describe the new technological characteristics:
4. Explain how new characteristics could or could not affect safety or effectiveness:
5. Explain how descriptive characteristics are not precise enough: New device contains Code selection feature that was not in the predicate device, that the user has to select the specific code for each new strip vial.
6. Explain new types of safety or effectiveness question(s) raised or why the question(s) are not new:
7. Explain why existing scientific methods can not be used:
8. Explain what performance data is needed:
9. Explain how the performance data demonstrates that the device is or is not substantially equivalent: