

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION  
DECISION SUMMARY**

**I. Background Information**

**A. 510(k) Number**

K260235

**B. Applicant**

Datar Cancer Genetics Pvt Ltd

**C. Proprietary and Established Names**

CellDx - Tissue

**D. Regulatory Information**

Product Code:	PZM
Device Class:	Class II
Classification Regulation:	21 CFR 866.6080 – Next Generation Sequencing Based Tumor Profiling Assay
Classification Panel:	Pathology

**II. Submission/Device Overview**

**A. Purpose for submission**

New Device

**B. Measurand**

Somatic single nucleotide variants (SNVs), insertions and deletions (Indels), Copy Number Variant (CNV), and select gene fusions, in genomic DNA and RNA isolated from formalin-fixed paraffine-embedded (FFPE) tumor tissue. A complete list of genes included in the assay can be found in Appendix A.

**C. Type of Test**

Next-generation sequencing tumor profiling test

**III. Intended Use/Indications for Use**

**A. Intended Use(s)**

The CellDx-Tissue is a qualitative in vitro diagnostic (IVD) test that uses next-generation sequencing of DNA and RNA isolated from formalin-fixed paraffine-embedded (FFPE) tumor tissue from patients previously diagnosed with solid malignant neoplasm to detect tumor gene

alterations in a broad multi-gene panel. The test is intended to provide tumor mutation profiling information on somatic mutations, including single nucleotide variants (SNVs), insertions and deletions (indels), one gene amplification, and three fusions.

Information provided by the CellDx-Tissue test is intended for use by qualified healthcare professionals in accordance with professional guidelines in oncology. Results from CellDx-Tissue are not intended to be conclusive or prescriptive for the labeled use of any specific therapeutic product. CellDx-Tissue is a single-site assay performed at Datar Cancer Genetics (DCG).

#### **B. Indication(s) for Use:**

Same as above

#### **C. Special Instrument Requirements**

Rx – For prescription use  
For *in vitro* diagnostic use

#### **D. Special Instrument Requirements**

Thermo Fisher Ion GeneStudio S5 Prime NGS Systems (qualified by Datar Cancer Genetics)

### **IV. Device/System Characteristics**

#### **A. Device Description**

The CellDx-Tissue utilizes an amplicon-based sequencing approach to detect a broad range of clinically significant somatic alterations, including SNVs, InDels, and a gene amplification in tumor tissue DNA (ttDNA), and three fusions in tumor tissue RNA (ttRNA). The test requires FFPE tumor tissue specimen having a surface area of 125 mm<sup>2</sup> or more, a thickness of at least 3 mm, with at least 20% tumor content. Specimens with lower tumor content may be enriched by macro-dissection. A description of required equipment, software, reagents, vendors, and storage conditions is provided in the product labeling (Datar Cancer Genetics Technical Information). The CellDx-Tissue system includes a sequencing instrument, reagents (nucleic acid extraction, library preparation, and sequencing) software (operation of the sequencing instrument and variant calling), and instructions for use (IFU) of the system. All instruments and reagents are qualified by Datar Cancer Genetics.

##### **1. Sample Preparation**

The CellDx – Tissue requires genomic DNA and RNA isolated from formalin-fixed paraffin-embedded (FFPE) tissue specimens. The tumor volume and minimum tumor content needed to obtain sufficient DNA and RNA for testing to achieve the stated performance are listed in Table 1. The tumor content of the FFPE specimen is determined by hematoxylin-eosin (HE) staining. If the histopathologist assessed tumor content is < 20%, the tissue should be macro-dissected to select as much viable tumor as possible and minimize the amount of adjacent non-tumor tissue.

**Table 1. Specimen Handling and Processing for Validated Specimen Type**

<b>Tissue Type</b>	Formalin-fixed paraffin-embedded (FFPE) tumor tissue blocks
<b>Specimen Volume</b>	≥ 3 mm thick; >125 mm <sup>2</sup> total surface area
<b>Tumor Content</b>	≥ 20%
<b>Macrodissection Requirements</b>	If the minimum tumor content is less than 20%, macrodissection will be done to enrich neoplastic content
<b>Limitations</b>	Low (<20%) tumor content, Inadequate tissue volume, Long-term storage degradation.
<b>Storage</b>	Room Temperature

**2. Nucleic Acid Extraction**

DNA extraction is performed using the QIAamp DNA FFPE Advanced UNG Kit (qualified by DCG) as per the manufacturer’s instructions for use. The procedure includes Uracil-N-Glycosylase (UNG) treatment to remove deaminated cytosine residues and an RNase A digestion to enhance the purity of the extracted DNA. The yield of extracted DNA is evaluated by real-time PCR. At least 20 ng DNA is required for library synthesis. RNA extraction procedure includes a DNase digestion step to eliminate contaminating DNA. The yield of extracted RNA is determined using the fluorometric assay, based on a standard curve generated from calibration standards. At least 20 ng RNA is required for reverse transcription (cDNA synthesis) and library synthesis. The assay has been validated with extracted DNA and RNA stored at -20°C and -80°C, respectively for up to 28 days.

**3. Library Synthesis**

DNA libraries are synthesized using the Oncomine DNA kit (qualified by DCG). The assay employs an amplicon-based targeted enrichment workflow, in which predefined genomic regions are amplified and converted into sequencing-ready libraries through adapter ligation and barcode incorporation. This approach enables multiplexed sequencing across the 517-gene panel. Library preparation entails the targeted amplification of genomic regions, partial digestion of primer sequences, ligation of barcode adapters, and subsequent clean-up and amplification to produce sequencing-competent libraries.

RNA libraries are synthesized using the Oncomine RNA Kit (qualified by DCG). The assay employs an amplicon-based targeted RNA sequencing workflow, in which input RNA is reverse-transcribed to cDNA and targeted regions are selectively amplified to generate sequencing-ready libraries. The workflow includes cDNA synthesis, targeted amplification, adapter ligation with unique barcodes, and subsequent clean-up steps to prepare multiplexed libraries suitable for NGS. DNA and RNA libraries are quantified using the kit which employs a real-time PCR-based method to measure the concentration of adapter-ligated libraries prior to sequencing. Library concentration results are reviewed against predefined quality acceptance criteria, and all quantitation outputs are retained as part of routine quality assurance and traceability.

#### **4. Template Preparation and Chip Loading**

Template preparation and chip loading are performed on the Ion Chef™ System (qualified by DCG) using the manufacturer's automated workflow for emulsion PCR (ePCR), enrichment, and chip loading. Quantified libraries meeting the predefined quality criteria are pooled and loaded onto the instrument along with the required reagents and consumables supplied in the template preparation kit. The Ion Chef System (qualified by DCG) performs clonal amplification of library molecules on beads through ePCR, followed by automated enrichment to isolate beads containing amplified templates. The instrument then loads the enriched, template-positive beads onto an Ion 550™ sequencing chip (qualified by DCG), producing a sequencing-ready chip for downstream processing on the compatible semiconductor sequencing platform qualified by Datar Cancer Genetics. Instrument checks and onboard controls ensure workflow integrity, and template-prepared chips proceed directly to sequencing.

#### **5. Sequencing**

Sequencing is performed on the Ion GeneStudio™ S5 Prime System (qualified by DCG) using the manufacturer's semiconductor sequencing workflow. The sequencing chip generated during template preparation is loaded into the instrument along with the required sequencing reagents once all automated system checks confirm acceptable setup. A predefined run plan specifies the assay configuration, sample assignments, and chip parameters. After sequencing is initiated, the system performs automated nucleotide flows and signal detection to generate base-level data across all loaded libraries. Sequencing data are accepted only when run-level QC metrics and control results meet validated acceptance criteria.

#### **6. Sequencing Data Analysis**

##### ***a. Data Management***

Sample tracking, data processing, and archival of sequencing data for the CellDx-Tissue assay are managed through a Laboratory Information Management System (LIMS) integrated with DCG's validated bioinformatics pipeline operating on a high-performance server environment. The system tracks and archives run-associated metadata including barcode identifiers, sequencing run identifiers, sample accession numbers, specimen source, library batch identifiers, and assay modality (DNA or RNA). Key functions include tracking sample status across all stages of data analysis, logging analysis iterations applied to each sample, recording software, algorithm, and database versions used for analysis, and archiving pipeline output files (BAM and VCF) together with sequencing run statistics (e.g., total reads generated, mean depth of coverage, uniformity, and on-target rates). All annotation and population databases are maintained as access-restricted, version-locked static copies in a controlled repository. Any update to software, databases, or internal scripts is managed under a formal change-control process requiring documented impact assessment, regression testing using predefined golden datasets, QA/RA approval, and release prior to clinical use.

***b. Signal process, Base Calling, Read Alignment, BAM Generation, and Coverage Statistics***

Primary data processing is performed using Torrent Suite Software (qualified by DCG). During sequencing, raw semiconductor flow signals are converted into nucleotide bases through signal processing and base calling algorithms. Reads are assigned to individual samples through barcode classification, followed by trimming of barcode and adapter sequences and quality filtering. Filtered reads are aligned to the human reference genome (hg19 / GRCh37.p5) using platform-validated alignment algorithms. Aligned reads are written to Binary Alignment Map (BAM) files, which serve as the input for downstream variant calling and copy number analysis. The software also generates coverage statistics including mean depth of coverage, base coverage uniformity, and percentage of reads on target for each sample.

***c. Read Alignment Quality Check***

Read alignment metrics are used to assess sequencing and library quality. Reads are aligned to hg19, and mismatches arising from biological variation (true variants) or technical sequencing errors are recorded. Alignment performance is summarized using the Alignment Quality (AQ) score derived from a Phred-scaled  $-10\log_{10}$  transformation. The manufacturer default setting of AQ17, corresponding to a base accuracy of 98% (2% error rate), is applied. AQ17 allows longer effective read lengths while accommodating mismatches expected in variant-rich tumor samples. Decreases in AQ are expected in samples with higher mutational burden and are evaluated in the context of overall run and sample QC metrics.

***d. Run - Level and Sample - Level Quality Control Checks***

***i. Run - Level QC***

A sequencing run is considered valid only if all of the following criteria are met:

- (1) The total sequencing output is  $\geq 40$  million reads per run.
- (2) Acceptable alignment quality (AQ17).

Runs failing any of the above criteria are invalid, investigated, and repeated following troubleshooting and corrective actions.

***ii. Sample – Level QC***

The sequencing performance is evaluated using the following criteria:

- (1) For DNA libraries, the mean target depth is  $\geq 500x$  and  $\geq 90\%$  of target amplification with  $\geq 100$  reads.
- (2) For RNA Libraries, the total fusion-mapped reads per run is  $\geq 500,000$  and reads per RNA pool is  $\geq 100,000$ .
- (3) The Contamination Score is  $< 0.120$

Samples that fail QC due to specimen quality, extraction, or library preparation may be re-sequenced, re-prepared, or re-extracted according to the predefined decision rules.

***e. Mutation Calling: Single Nucleotide Variants (SNVs) and Insertions / Deletions (indels)***

**i. Analysis of Positive and Negative Controls**

Positive and negative controls are analyzed using the same coverage requirements as patient samples (mean depth  $\geq 500\times$  and  $\geq 100\times$  coverage across  $\geq 90\%$  of target regions). Expected variants must be detected in positive controls within validated variant allele frequency (VAF) ranges, and no reportable variants may be detected in negative controls.

**ii. Filters on Sample Coverage**

Tumor samples must achieve a minimum mean sequencing depth of  $\geq 500\times$  to be eligible for variant calling and reporting.

**iii. Filtering for High-Confidence Mutations**

Raw SNV and indel calls are subjected to locked filtering thresholds to ensure only high-confidence somatic variants are reported. Variants are classified as:

- (1) Hotspot variants: VAF  $\geq 2\%$ , total depth  $\geq 40$ , mutant reads  $\geq 8$ , strand bias  $\leq 0.9$
- (2) Non-hotspot variants: VAF  $\geq 5\%$ , total depth  $\geq 40$ , mutant reads  $\geq 10$ , strand bias  $\leq 0.85$

Variants are annotated using AMP/ASCO/CAP somatic guidelines, HGVS nomenclature, and major public databases (e.g., ClinVar, gnomAD, dbSNP, 1000 Genomes). Population variants with a frequency  $\geq 5\%$  in any gnomAD or 1000 Genomes subpopulation are excluded. Filtering thresholds are fixed for clinical use, and any modification requires formal change-control and re-validation.

***f. Mutation Annotation***

Variants with functional consequences including missense, nonsense, frameshift, in-frame insertions/deletions, splice-site, and splice-region alterations are retained. Annotation incorporates gene context, predicted functional impact, and clinical relevance. Common polymorphisms present in population databases or internal normal datasets are excluded from somatic reporting.

***g. Tumor Purity Analysis***

Tumor purity of the sequenced specimen is estimated bioinformatically using allele frequencies of informative variants and sequencing read distributions. This estimate may differ from the pathologist-reported tumor percentage due to macrodissection or sampling effects. When discrepancies or software limitations are identified, tumor purity is manually reviewed using pathologist estimates and driver mutation VAFs to ensure accurate interpretation of variant calls.

#### ***h. Copy Number Analysis***

Copy number variation analysis is performed for *ERBB2* using a variability-corrected informatics baseline derived from normal samples. Amplicon coverage log<sub>2</sub> ratios are normalized and adjusted for tumor cellularity to estimate copy number and confidence intervals. *ERBB2* amplification is reported as Positive when the observed copy number is  $\geq 8.5$ . Copy numbers  $\geq 4$  and  $< 8.5$  are reported as Indeterminate / Low-Level Amplification (below the LoD).

#### ***i. Fusion Analysis***

RNA-based fusion detection is validated for *ALK*, *RET*, and *ROSI*. Fusions are identified using split-read and discordant read-pair evidence. Clinically reportable fusions are reported only when  $\geq 500$  fusion-supporting reads are detected and validated against internal positive controls. Fusion calls rely exclusively on RNA-level evidence.

### **7. Controls**

- a. Positive Control (PC):*** An external control that consist of a mixture of well-characterized cell line DNA with known somatic variants at defined variant allele frequency (VAF), and RNA with defined fusions (Table 2) and is included in each assay run (batch level).
- b. Negative Control (NC):*** An external control that is derived from genomic DNA and RNA of a known wild-type (WT) material (Table 2) and is included in each assay run (batch level).

**Table 2. Positive and Negative Controls Used in the CellDx-Tissue Assay**

<b>Control Type</b>	<b>Qualification</b>
Positive Control DNA	Qualified by DCG
Positive Control RNA	
Negative Control DNA	
Negative Control RNA	

- c. No-Template Control (NTC):*** DNase, RNase free water is used as the no template control and will be included in each assay run (batch level).

### **8. Results Reporting**

The CellDx-Tissue assay report provides structured information on detected genomic alterations, level of clinical significance and their clinical implications. The detected genomic alterations are categorized into Level 2 (Variants with Evidence of Clinical Significance) or Level 3 (Variants with Potential Clinical Significance). The CellDx-Tissue does not report mutations in 217 regions among 13473 interrogated target regions due to low coverage and high GC content. A list of all 517 genes is provided in Appendix A; and a list of excluded regions in the genes or genes with variants that are excluded due

to challenging regions (e.g., low complexity / repeats) is provided in Appendix B and Appendix C, respectively.

## 9. Quality Metrics

Quality metrics are evaluated across the following categories:

- Run Level: Metrics that evaluate the overall sequencing performance, including instrument functions, reagent performance, data yield, and controls.
- Sample Level: Metrics that evaluate the sequencing quality, library performance, and coverage sufficiency.
- Variant (Analyte) Level: Metrics that are quantified for locus-level coverage and read support. Variants passing analyte-level QC are reported.

**Table 3. Quality Control Metrics and Assay Cutoff**

Quality Metric	Frequency	Acceptance Criteria
Specimen	Sample	Labels intact and Legible, Test Requisition Form(TRF) completely filled, consent signed, No loss of or damage to specimen, Number of FFPE blocks matches TRF, FFPE specimen identifiers match TRF and, LIMS entry
Specimen	Sample	>125mm <sup>2</sup> surface area; ≥3mm thick Tumor content ≥20% (Macro-dissect if <20%)
DNA Yield	Sample	>20 ng
RNA Yield	Sample	>20 ng
Library Yield	Sample	>100 pM
Total Data Output	Run	≥ 40 million reads
Positive Control (PC)	Run (Batch)	DNA: All variants detected within expected VAF ranges  RNA: All 6 fusions detected.

Quality Metric	Frequency	Acceptance Criteria
Negative Control (NC)	Run (Batch)	DNA: None of the 890 hotspot mutations should be detected RNA: No fusions detected.
Base Quality (AQ Score)	Run	AQ17 (98% accuracy)
Mean Target Coverage (DNA)	Sample	$\geq 500x$ .
Coverage Uniformity	Sample	$\geq 90\%$ of target regions with $\geq 100x$ coverage
Contamination Score	Sample	$\leq 0.120$
Mean Target Coverage (RNA)	Sample	Total mapped fusion reads $\geq 500,000$ Pool 1, Pool 2 (each): $>100,000$ reads
SNVs and Indels Calling	Variant	Mutation coverage (total depth) $\geq 40$ Number of mutant reads (variant read count) $\geq 8$ hotspot variants / $\geq 10$ non-hotspot variants Variant Allele Fraction (VAF) $\geq 0.02$ hotspot variants / $\geq 0.05$ non-hotspot variants Strand bias $< 0.9$ hotspot variants / $< 0.85$ non-hotspot variants
<i>ERBB2</i> Amplification	Variant	CN $\geq 8.5$ is Positive CN $\geq 4$ to $< 8.5$ is Indeterminate / Low Level Amplification CN $< 4$ is Negative
<i>ALK</i> , <i>RET</i> , and <i>ROS1</i> Fusions	Variant	Positive where supporting read counts are $\geq 500$

Quality Metric	Frequency	Acceptance Criteria
Test Failure Criteria	Sample	DNA mean coverage < 500x, or Coverage uniformity < 90%, or Contamination score >0.120, or Fusion mapped reads <500,000, or Pool reads <100,000
Reprocessed Specimen	Sample	All above QC criteria met
Report	Sample	Patient ID matches LIMS entry  Accreditation Logos match standard template  Signing authority signatures present

## B. Principles of Operation

The CellDx-Tissue utilizes an amplicon-based sequencing approach to detect a broad range of clinically significant somatic alterations, including small nucleotide variants (SNVs), insertions and deletions (Indels), and *ERBB2* gene amplification in tumor tissue DNA (ttDNA); and *ALK*, *RET*, and *ROS1* fusions in tumor tissue RNA (ttRNA). The test requires FFPE tumor tissue specimen having a surface area of more than 125 mm<sup>2</sup>, a thickness of at least 3 mm, with at least 20% tumor content. Specimens with lower tumor content may be enriched by macro-dissection. The assay involves target enrichment and deep sequencing of specific regions across a comprehensive panel of 517 genes, including oncogenes, tumor suppressor genes, and other clinically relevant genes. DNA and RNA primers are designed to target all coding exons, selected introns and targeted RNA sequences including chimeric fusion reads. Sequence libraries are prepared using DNA and RNA through a multiplex polymerase chain reaction (PCR) amplification step to enrich the target sequences. The target sequences are tagged with unique barcode oligonucleotides for individual sample identification and adaptor oligonucleotides that facilitate anchoring to the sequencing platform. These target sequences are then clonally amplified on microscopic beads using emulsion PCR before sequencing. Multiple barcoded sequence libraries are subsequently pooled and loaded onto a semiconductor sequencing chip. The resulting sequence reads are then aligned to the reference human genome (hg19), to identify gene variants.

## C. Determination of Assay Thresholds


### 1. Requirements on Exon Coverage

A power analysis was conducted to estimate the minimum sequencing depth (total number of reads) needed to detect a mutation with a true underlying variant allele fraction (VAF) of 0.02 or greater, for varying levels of statistical power (0.8 to 0.99), assuming a fixed alpha (Type I error rate) of 0.05. The 95% confidence interval (CI) ranges of observed VAF as a function of sequencing depth were also calculated. This study showed that when

a mutation is present at 0.1 VAF, the 95% CI with a sequencing depth of 500X is expected to fall between 0.075 and 0.13. When the overall coverage is 100X, the 95% CI for a mutation at 10% VAF is estimated to fall between 0.05 and 0.176.

To confirm these estimates, empirical data was obtained to measure the range of observed VAF for expected VAF, utilizing DNA from 20 formalin-fixed, paraffin-embedded (FFPE) normal tissue specimens from unrelated individuals. Equimolar parts of the DNA from these specimens were mixed to create five secondary specimen pools with a range of genetic variants (SNPs) having expected frequencies as low as 2%. A total of 863 common SNPs were considered for this study.

The observed variant frequencies for these common SNPs genotyped in the pooled normal samples, binned by their true underlying variant frequency is shown in Figure 1. The empirical data demonstrated a strong correlation between expected and observed VAF (Pearson's  $r = 0.99$ ), with a slope of 0.99 and an intercept of 0.013. For SNPs with a true underlying variant fraction of 0.1, the observed variant fraction ranged from 0.069 to 0.138 when the mean coverage of the specimen was 1151X. This range is consistent with the theoretical statistical assessment for a depth of 500X (0.075 to 0.13). This data supports using 0.05 as the lower limit for reporting mutations detected with a true underlying frequency of 10%.





**Figure 1. Observed vs. Expected Variant Allele Frequency for the CellDx-Tissue Assay**

## 2. Requirements on Sample Coverage

Twenty FFPE normal (diploid) specimens were profiled using CellDx-Tissue to generate summary statistics across all targeted exons. The mean coverage across all amplicons at optimal depth was 2087X (Range: 1420X-3369X, SD = 481X). The percentage of amplicons with >100X coverage was 99.3% (Range: 98.0-99.9%, SD = 0.5%). When mapped to the exon level, the mean coverage across all targeted exons was 2087X (Range: 1420X-3371X, SD = 482X). Summary statistics were also computed on coverage values per exon normalized by per-sample coverage. Exons with consistently low coverage (median normalized coverage < 0.05) were excluded from SNV/indel reporting, primarily due to high GC content (Appendix B), with the exception of two *TERT* promoter amplicons where variant calling parameters were relaxed with medical justification.


Sequence coverage was further evaluated to establish minimum criteria for the analysis and reporting of variants. Based on the power calculations, a minimum 100X coverage threshold per exon is required to call mutations with a true underlying mutation frequency of 10% or greater, with 95% power at an alpha level of 0.05 (Figure 2).





**Figure 2. Distribution of Mean and Medium Coverage Values for the Targeted Regions of the CellDx-Tissue Assay using High Coverage Samples.** The dashed blue line indicates 100X coverage.

A second set of 20 samples were evaluated at lower sequencing depth. The mean coverage across all amplicons was 471X (Range: 284X-663X, SD = 101X). The percentage of amplicons with >100X coverage was 92.6% (Range: 86.6-96.2%, SD = 2.4%). When mapped to the exon level, the mean coverage across all targeted exons was 471X (Range: 284X-662X, SD = 101X). With the same remaining exons across all genes, 94.1% (Range: 88.0-97.7%, SD = 2.42%) were sequenced to a depth of 100X or greater. The distribution of these mean and medium coverage values for the targeted exons is shown in Figure 3.



**Figure 3. Distribution of Mean and Medium Coverage Values for Targeted Regions of the CellDx – Tissue Assay with Low Coverage Samples.** The dashed blue line indicates 100X coverage.

Based on these two studies, the threshold of  $\geq 90\%$  of amplicons achieving 100X coverage was established as a critical quality metric for the CellDx-Tissue assay.

### **3. Requirements on Variant Calling Thresholds**

Variant calling parameters such as sequence coverage (Alternate Allele Observation, AO), variant coverage (Coverage Depth, DP) and strand bias (SB) were assessed as filters for specificity while maintaining the ability to detect true positive calls. Thresholds were established to ensure specificity is maintained at targeted VAF levels for reporting, particularly at 0.02 and 0.05 for hotspot and non-hotspot categories, respectively.

The cutoffs for AO, DP, and SB for hotspot variants and non-hotspot variants were established using two key development studies. The first study was designed to select the optimal parameter cutoff combination utilizing a dataset of 20 CellDx-Tissue runs using a Reference Standard. This involved evaluating multiple cutoffs for each of the three parameters, generating a total of 7,650 combinations. A list of 3,723 positive variants and 4,630 negative variants was used to evaluate the performance of each parameter combination, with the goal of achieving a Positive Predictive Value (PPV) greater than 96%. The relationship between these parameters and PPV, highlighting the selected cutoff combination: AO = 8, DP = 40, SB = 0.90 for hotspot variants; and AO = 10, DP = 40, SB = 0.85 for non-hotspot variants is illustrated in Figure 3.

The second study aimed to determine the optimal cutoff for Coverage Depth (DP) and minimum variant read count. This study used ten non-cancer FFPE samples and evaluated a total of 902 hotspot variants and 6,418 non-hotspot variants at different Coverage Depth cutoffs. The lowest Coverage Depth (DP) that effectively filtered out >99.7% of noise variants were selected, corresponding to a minimum of 8 reads for hotspot variants and 10 reads for non-hotspot variants.

### **4. Tumor-Only vs. Matched-Normal Germline Filtering**

A study was conducted to demonstrate analytical equivalence between the CellDx-Tissue tumor-only (T/O) germline filtering approach and a matched tumor-normal (T/N) reference method, as required under 21 CFR 866.6080(b)(1)(iv). The study utilized 52 prospectively selected matched peripheral blood and tumor specimens representing 16 solid tumor types (including bladder / urinary tract (n = 1), bowel (n = 10), breast (n = 13), cervix (n = 2), esophagus / stomach (n = 1), head and neck (n = 3), liver (n = 1), lung (n = 6), ovary / fallopian tube (n = 5), pancreas (n = 3), peritoneum (n = 1), prostate (n = 1), testis (n = 1), thymus (n = 1), uterus (n = 2), and unknown primary (n = 1), and diverse ethnic representation (65% Asian, 27% White/Caucasian, 6% African American and 2% unspecified ethnicity from the US).

The T/O pipeline employed an in-silico filtering strategy that excluded variants with population frequency  $\geq 5\%$  in gnomAD (v2.1.1) or 1000 genome (v Phase 3v5), while the reference method used variants called using torrent variant caller (TVC, v5.18) with total white blood cell (WBC) genomic DNA (gDNA) as the germline control. Clinically actionable variants were defined as those with ClinVar (v 20250323) pathogenic/likely pathogenic classification or OncoKB (v 4.27) Level 1/2 oncogenic annotation with population frequency  $\leq 0.1\%$ . A total of 345 clinically actionable variants were evaluated, including germline pathogenic alterations in *BRCA2*, *CHEK2*, *ATM*, *MSH2*, *PALB2*, and *RAD51D*.

The study demonstrated 100% PPA (345/345; 95% CI: 98.9% - 100%) and 100% NPA (2,745/2,745; 95% CI: 99.9% - 100%), with zero false negatives (FN) for actionable germline variants. Twenty-one (6.1%) actionable variants had VAFs between 0.02 -0.05, confirming assay sensitivity near the LoD. This analysis concluded the T/O filtering strategy is analytically equivalent to matched T/N subtraction.

**D. Substantial Equivalence Information**

**1. Predicate Device Name(s)**

MSK-IMPACT

**2. Predicate 510(k) Number(s)**

DEN170058

**3. Comparison with Predicate(s)**

Characteristics	Predicate Device: MSK-IMPACT	Subject Device: CellDx-Tissue
<b>Similarities</b>		
<b>Indications for Use</b>	The MSK-IMPACT assay is a qualitative in vitro diagnostic test that uses targeted next generation sequencing of formalin-fixed paraffin-embedded tumor tissue matched with normal specimens from patients with solid malignant neoplasms to detect tumor gene alterations in a broad multi gene panel.	The CellDx-Tissue is a qualitative in vitro diagnostic (IVD) test that uses next-generation sequencing of DNA and RNA isolated from formalin-fixed paraffine-embedded (FFPE) tumor tissue from patients previously diagnosed with solid malignant neoplasm to detect tumor gene alterations in a broad multi-gene panel. The test is intended to provide tumor mutation profiling information on somatic mutations, including single nucleotide variants (SNVs), insertions and deletions (indels), one gene amplification, and three fusions.

<b>Characteristics</b>	<b>Predicate Device: MSK-IMPACT</b>	<b>Subject Device: CellDx-Tissue</b>
<b>Similarities</b>		
<b>Indications for Use</b>	The test is intended to provide information on somatic mutations (point mutations and small insertions and deletions) and microsatellite instability for use by qualified health care professionals in accordance with professional guidelines and is not conclusive or prescriptive for labeled use of any specific therapeutic product. MSK-IMPACT is a single-site assay performed at Memorial Sloan Kettering Cancer Center.	Information provided by the CellDx-Tissue test is intended for use by qualified healthcare professionals in accordance with professional guidelines in oncology. Results from CellDx-Tissue are not intended to be conclusive or prescriptive for the labeled use of any specific therapeutic product. CellDx-Tissue is a single-site assay performed at Datar Cancer Genetics (DCG).
<b>Specimen Type</b>	Formalin-fixed, paraffin-embedded (FFPE) tumor tissue matched with normal specimens from patients with solid malignant neoplasms	Formalin-fixed, paraffin-embedded (FFPE) tumor tissue from patients with solid malignant neoplasms.
<b>Target Population</b>	Patients with solid malignant neoplasms	Same
<b>Test Environment</b>	Single Site Assay	Same
<b>Assay Cut-Off</b>	Does not report mutations below 2% for known hotspot mutations and 5% for non-hotspot mutations	Same
<b>Report Format</b>	Results are reported under one of the following categories: <ul style="list-style-type: none"> <li>• Variants with Evidence of Clinical Significance</li> <li>• Variants with Potential Clinical Significance</li> </ul>	Same

<b>Characteristics</b>	<b>Predicate Device: MSK-IMPACT</b>	<b>Subject Device: CellDx-Tissue</b>
<b>Similarities</b>		
<b>Clinical Evidence Curation</b>	<ul style="list-style-type: none"> <li>• Uses OncoKB, knowledge base</li> <li>• Classification criteria were developed by MSK to communicate the level of clinical evidence available for individual mutations in the test report.</li> </ul>	Variant calls are organized into Variant with Evidence of Clinical Significance or Variant with Potential Clinical Significance depending on the designated cancer type.
<b>Differences</b>		
<b>Genes on Panel</b>	468 (6,357 exons)	517 (6,597 exons)
<b>Variant Types</b>	Somatic Variants including point mutations and small insertions and deletions, and microsatellite instability	Somatic Variants including point mutations and small insertions and deletions, <i>ERBB2</i> amplification, and three fusions ( <i>ALK</i> , <i>RET</i> , and <i>ROS1</i> )
<b>Black List</b>	73 exons	157 exons
<b>Analyte</b>	DNA	DNA and RNA
<b>Minimum Tumor Content</b>	10% (20% preferred; 25% for MSI)	> 20%
<b>Nucleic Acid Input</b>	50 ng – 250 ng	DNA: 20 ng – 50 ng RNA: 20 ng – 50 ng
<b>Sequencing Instrument</b>	Illumina HiSeq 2500 Sequencer	Thermo Fisher Ion GeneStudio S5 Prime System
<b>Technology</b>	Hybrid Capture	Amplicon
<b>Sequencing Chemistry</b>	Sequencing by Synthesis (Fluorescent Detection)	Sequencing by Synthesis (Hydrogen Ion Detection)
<b>Controls</b>	<ul style="list-style-type: none"> <li>• Positive Control</li> <li>• Negative Control</li> <li>• No Template Control (NTC)</li> <li>• Matched Normal</li> </ul>	<ul style="list-style-type: none"> <li>• Positive Control</li> <li>• Negative Control</li> <li>• No Template Control (NTC)</li> </ul>

Characteristics	Predicate Device: MSK-IMPACT	Subject Device: CellDx-Tissue
<b>Differences</b>		
<b>Germline Filtering</b>	Matched normal specimen analysis	Computational, population database filtering
<b>Coverage Requirements</b>	≥ 200x mean; 100x for ≥ 98% target regions	≥ 500x; 100x for ≥ 90% target regions
<b>Variant Calling Thresholds</b>	<u>Hotspot:</u> <ul style="list-style-type: none"> <li>• Mutant reads (AD) ≥ 8</li> <li>• Mutation coverage (DP) ≥ 20</li> <li>• Mutation frequency (VAF) ≥ 0.02</li> </ul> <u>Non-hotspot</u> <ul style="list-style-type: none"> <li>• Mutant reads (AD) ≥ 10,</li> <li>• Mutation coverage (DP) ≥ 20</li> <li>• Mutation frequency (VAF) ≥ 0.05</li> </ul>	<u>Hotspot:</u> <ul style="list-style-type: none"> <li>• Mutant reads (AD) ≥ 8,</li> <li>• Mutation coverage (DP) ≥ 40</li> <li>• Strand bias (SB) &lt; 0.9</li> <li>• Mutation frequency (VAF) ≥ 0.02</li> </ul> <u>Non-hotspot</u> <ul style="list-style-type: none"> <li>• Mutant reads (AD) ≥ 10,</li> <li>• Mutation coverage (DP) ≥ 40</li> <li>• Strand bias (SB) &lt; 0.85</li> <li>• Mutation frequency (VAF) ≥ 0.05</li> </ul>
<b>Contamination QC</b>	Percent heterozygous sites at fingerprint SNPs < 55%; Average MAF at homozygous fingerprint SNPs < 2%	Estimated contamination score < 0.120, Contamination score estimated from the signal derived from reference reads at homozygous alternate sites.
<b>Criteria for Calling “Failed” Samples</b>	If a sample presents with mean coverage across all exons < 50x and no mutations are detected due to the low overall coverage, the test is deemed “failed” for the sample	If a sample presents with mean coverage across all exons < 500x for DNA and <500,000 reads for RNA, the test is deemed “failed” for the sample.

## E. Standards/Guidance Documents Referenced

The following FDA guidance documents and standards were consulted:

1. Medical Devices – Quality Management Systems: Requirements for Regulatory Purposes. ISO 13485: 2016
2. Medical Devices – Application of Risk Management to Medical Devices. ISO 14971: 2019
3. Medical Laboratories – Requirements for Quality and Competence. ISO 15189: 2022
4. Medical Devices – Symbols to be used with Information to be Supplies by the Manufacturer – Part 1: General Requirements. ISO 15223-1: 2021
5. Medical Device Software – Software Life Cycle Processes. ISO / IEC 62304: 2006
6. Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements. IEC 61010-1: 2010
7. Electrical Equipment for Measurement, Control and Laboratory Use – EMC Requirements – Part 1: General Requirements. EN 61326-1: 2013
8. Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline - CLSI EP25-A
9. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline - Second Edition. CLSI EP17-A2
10. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline - Third Edition. CLSI EP05-A3
11. User Protocol for Evaluation of Qualitative Test Performance; Approved Guideline - Second Edition. CLSI EP12-A2

## F. Performance Characteristics

### 1. Analytical Performance – General

The CellDx – Tissue assay is a targeted NGS panel with 517 genes. The targeted regions of interest in CellDx – Tissue are designed to detect SNVs, insertions up to 20 bp and deletions up to 30 bp in length of the targeted genes, as well as *ERBB2* amplifications, *ALK*, *RET*, and *ROSI* fusions. For SNVs and indels, a representative approach to validation of the targeted genes in the panel was submitted with data representing variant types for SNVs and indels, and at the gene levels for amplifications and fusions indicated with this assay. In addition, the assay was evaluated for performance regarding the panel-wide quality metrics.

#### a. *Invalid Rates*

Multiple factors can influence the overall robustness and performance of complex molecular tests, including pre-analytical factors and overall sample quality. If key in-process or automated data quality metrics are not met, the CellDx – Tissue assay supports repeat samples through the workflow once. Performance throughout the verification and validation of the device was tracked, and a summary of the rates for first pass (no repeat) and overall pass (allowing for a single repeat) are presented in

Table 4. Data were aggregated for clinical cases from >25 tumor types. Resulting pass rates by tumor type across the workflow are shown in Table 5.

**Table 4. Acceptability Rates of CellDx – Tissue Assay**

Clinical FFPE Samples	Acceptability Rate (n/N) (95% CI)
First Pass	81.8% (2190/2676) (80.33%–83.25%)
After Repeat Test	93.1% (2492/2676) (92.10%–94.02%)

**Table 5. Comparability of Tumor Invalid Rates for the CellDx – Tissue Assay**

Organ	Number of Samples Failed				Invalid Rate (n/N) (%)
	Pre-analytical	Pre-run	Post-Run		
	Low Tumor Purity <sup>1</sup>	Low Nucleic Acid Yield <sup>2</sup>	Library Failure <sup>3</sup>	Data QC Failure	
Adrenal Gland	0	0	0	0	(0/2) 0%
Ampulla of Vater	1	0	0	0	(1/33) 3%
Biliary Tract	0	0	0	0	(0/7) 0%
Bladder/Urinary Tract	0	0	0	2	(2/20) 10%
Bone	0	0	0	0	(0/2) 0%
Bowel	6	4	2	35	(47/479) 10%
Breast	3	5	2	3	(13/317) 4%
Cervix	0	0	0	0	(0/41) 0%
CNS/Brain	0	1	1	3	(5/67) 7%
Esophagus/Stomach	0	0	0	9	(9/115) 8%
Head and Neck	12	3	0	17	(32/220) 15%
Kidney	0	0	0	0	(0/18) 0%
Liver	0	0	0	0	(0/12) 0%
Lung	2	3	0	38	(43/695) 6%
Ovary/Fallopian Tube	0	1	0	3	(4/73) 5%
Ovary/Fallopian Tube, Uterus	0	0	0	0	(0/3) 0%
Pancreas	2	0	0	2	(4/27) 15%
Peritoneum	0	0	0	0	(0/2) 0%

Organ	Number of Samples Failed				Invalid Rate (n/N) (%)
	Pre-analytical	Pre-run	Post-Run		
	Low Tumor Purity <sup>1</sup>	Low Nucleic Acid Yield <sup>2</sup>	Library Failure <sup>3</sup>	Data QC Failure	
Pleura	0	1	0	1	(2/2) 100%
Prostate	0	0	0	6	(6/217) 3%
Skin	0	0	5	7	(12/51) 24%
Soft Tissue	1	0	0	0	(1/18) 6%
Testis	0	0	0	1	(1/2) 50%
Thymus	0	0	0	0	(0/4) 0%
Thyroid	0	0	0	0	(0/88) 0%
Unknown Primary	0	0	0	0	(0/2) 0%
Uterus	0	0	0	2	(2/155) 1%
Vulva/Vagina	0	0	0	0	(0/4) 0%
<b>Total</b>	<b>27</b>	<b>18</b>	<b>10</b>	<b>129</b>	<b>(184/2676) 6.9%</b>

<sup>1</sup> < 20 %

<sup>2</sup> < 20 ng

<sup>3</sup> No amplification

**b. Invalid Rates Observed in Method Comparison**

Invalid rates were infrequent across all tested biomarkers and were consistently within the acceptance criteria for the CellDx-Tissue assay workflow. A summary of the invalid rates observed during the method comparison study is shown in Table 6.

**Table 6. Invalid Rates Observed in the Method Comparison Study.**

Category	DNA	RNA	Total
1st Pass Valid	233	221	454
After Repeat Test Valid	10	23	33
Low Nucleic Acid Yield (<20 ng)	14	6	17
Library Failure	1	0	1
Invalid (Sequencing QC Failure)	7	5	12
Total Invalid	22	11	33
Total Attempts	265	252	517
<b>Total Invalid (%)</b>	<b>8.30%</b>	<b>4.4%</b>	<b>6.38%</b>

## 2. Precision

### *a. Precision Panel*

The within lab precision (between-run and within-run) of the CellDx – Tissue assay was assessed using 17 FFPE clinical samples originating from eight primary tumor types (Bowel (n = 5, 29.4%), Lung (n = 4, 23.5%), Prostate (n = 2, 11.8%), CNS (n = 1, 5.9%), Breast (n = 1, 5.9%), Ovary (n = 2, 11.8%), Esophagus (n=1, 5.9%) and Thyroid (n = 1, 5.9%) to represent different variants and a range of frequencies. Extracted DNA and RNA from each of the 17 samples was tested in duplicate by two different operators on multiple sequencing instruments across  $\geq 3$  non-consecutive days using two reagent lots, yielding 24-48 replicates per sample. For each replicate tested, all observed variants which passed the CellDx – Tissue assay QC metrics were reported and assessed for precision.

### *b. Panel-Wide Precision*

Precision was assessed for each variant across all replicates. The positive call rates (PCR) and negative call rates (NCR) were calculated along with the two-sided 95% confidence intervals.

Table 7 summarizes the positive call rate (PCR) and negative call rate (NCR) stratified by mutation type (SNV, insertions, and deletions) and variant allele frequency (VAF). An overall PCR of 98.69% across all samples and replicates (2871/2909; 95% CI: 98.21% – 99.05%; Average VAF range: 2.0 – 85.52), with an increase in PCR at higher VAFs observed, and an overall NCR of 100% (533539/533539; 95% CI: 99.99, 100.00).

The positive call rates for selected individual sequence mutations assessed in the precision study, along with the VAF range, mean SD, and percent CV per variant are presented in Appendix D. A total of 48 SNVs and 15 indels (3 insertions, 15 deletions) are provided.

**Table 7. Precision Positive (PCR) and Negative (NCR) Call Rates**

Variant Type	VAF Level (%)	Unique Mutations	PCR (%) (n/N)	NCR (%) (n/N)	AF Range (%)	Mean Allele Depth/FSR Range	Mean Loci Depth Rang
All	AF ≥ 0	63	98.69 (2871/2909)	100 (533539/533539)	2.00 - 85.52	46 - 1455	440 - 1984
	AF ≥ 2.0	63	98.69 (2871/2909)	100 (533539/533539)	2.00 - 85.52	46 - 1455	440 - 1984
	AF ≥ 5.0	60	98.99 (2758/2786)	100 (533539/533539)	3.20 - 85.52	46 - 1455	440 - 1978
	AF ≥ 10.0	53	99.23 (2451/2470)	100 (450769/450769)	7.43 - 85.52	82 - 1455	440 - 1978
	AF ≥ 15.0	48	99.15 (2211/2230)	100 (450769/450769)	15.44 - 85.52	142 - 1455	440 - 1978
Variants with Evidence of Clinical Significance	AF ≥ 0	10	99.56 (455/457)	100 (365377/365377)	2.26 - 79.57	74 - 1455	604 - 1978
	AF ≥ 2.0	10	99.56 (455/457)	100 (365377/365377)	2.26 - 79.57	74 - 1455	604 - 1978
	AF ≥ 5.0	9	100 (411/411)	100 (324483/324483)	20.07 - 79.57	290 - 1455	604 - 1978
	AF ≥ 10.0	9	100 (411/411)	100 (324483/324483)	20.07 - 79.57	290 - 1455	604 - 1978
	AF ≥ 15.0	9	100 (411/411)	100 (324483/324483)	20.07 - 79.57	290 - 1455	604 - 1978
Hotspot Variants	AF ≥ 0	24	97.99 (1073/1095)	100 (409829/409829)	2.00 - 60.77	46 - 1092	737 - 1984
	AF ≥ 2.0	24	97.99 (1073/1095)	100 (409829/409829)	2.00 - 60.77	46 - 1092	737 - 1984
	AF ≥ 5.0	21	98.77 (960/972)	100 (409829/409829)	3.20 - 60.77	46 - 1092	737 - 1978
	AF ≥ 10.0	17	99.25 (791/797)	100 (409829/409829)	15.44 - 60.77	214 - 1092	810 - 1978
	AF ≥ 15.0	17	99.25 (791/797)	100 (409829/409829)	15.44 - 60.77	214 - 1092	810 - 1978

N/A: Not Applicable; FSR: Fusion Supporting Reads

Variant Type	VAF Level (%)	Unique Mutations	PCR (%) (n/N)	NCR (%) (n/N)	AF Range (%)	Mean Allele Depth/FSR Range	Mean Loci Depth Rang
Non-hotspot Variants	AF ≥ 0	39	99.12 (1798/1814)	100 (407205/407205)	5.04 - 85.52	63 - 1455	440 - 1977
	AF ≥ 2.0	39	99.12 (1798/1814)	100 (407205/407205)	5.04 - 85.52	63 - 1455	440 - 1977
	AF ≥ 5.0	39	99.12 (1798/1814)	100 (407205/407205)	5.04 - 85.52	63 - 1455	440 - 1977
	AF ≥ 10.0	36	99.22 (1660/1673)	100 (324435/324435)	7.43 - 85.52	82 - 1455	440 - 1977
	AF ≥ 15.0	31	99.09 (1420/1433)	100 (324435/324435)	17.53 - 85.52	142 - 1455	440 - 1977
SNVs	AF ≥ 0	48	99.13 (2172/2191)	100 (449879/449879)	2.00 - 79.57	46 - 1455	449 - 1984
	AF ≥ 2.0	48	99.13 (2172/2191)	100 (449879/449879)	2.00 - 79.57	46 - 1455	449 - 1984
	AF ≥ 5.0	45	99.56 (2059/2068)	100 (449879/449879)	3.20 - 79.57	46 - 1455	449 - 1978
	AF ≥ 10.0	41	99.84 (1891/1894)	100 (408049/408049)	10.33 - 79.57	101 - 1455	449 - 1978
	AF ≥ 15.0	39	99.83 (1795/1798)	100 (408049/408049)	15.44 - 79.57	142 - 1455	449 - 1978
Insertions	AF ≥ 0	3	100 (144/144)	100 (85248/85248)	7.43 - 46.99	82 - 571	758 - 1430
	AF ≥ 2.0	3	100 (144/144)	100 (85248/85248)	7.43 - 46.99	82 - 571	758 - 1430
	AF ≥ 5.0	3	100 (144/144)	100 (85248/85248)	7.43 - 46.99	82 - 571	758 - 1430
	AF ≥ 10.0	3	100 (144/144)	100 (85248/85248)	7.43 - 46.99	82 - 571	758 - 1430
	AF ≥ 15.0	1	100 (48/48)	100 (42720/42720)	35.92 - 46.99	571 - 571	1375 - 1375
Deletions	AF ≥ 0	12	96.69 (555/574)	100 (254204/254204)	3.54 - 85.52	69 - 578	440 - 1699
	AF ≥ 2.0	12	96.69 (555/574)	100 (254204/254204)	3.54 - 85.52	69 - 578	440 - 1699
	AF ≥ 5.0	12	96.69 (555/574)	100 (254204/254204)	3.54 - 85.52	69 - 578	440 - 1699
	AF ≥ 10.0	9	96.30 (416/432)	100 (170736/170736)	9.28 - 85.52	108 - 578	440 - 1699
	AF ≥ 15.0	8	95.83 (368/384)	100 (170736/170736)	15.64 - 85.52	269 - 578	440 - 1699

N/A: Not Applicable; FSR: Fusion Supporting Reads

Variant Type	VAF Level (%)	Unique Mutations	PCR (%) (n/N)	NCR (%) (n/N)	AF Range (%)	Mean Allele Depth/FSR Range
<i>ERBB2</i> Amplification	1	100 (120/120)	-	N/A	N/A	N/A
<i>ALK</i> Fusion	1	100 (45/45)	-	N/A	7121-30568	N/A
<i>ROS1</i> Fusion	1	95.0 (38/40)	-	N/A	504-2328	N/A
<i>RET</i> Fusion	1	100 (46/46)	-	N/A	740-175619	N/A

N/A: Not Applicable; FSR: Fusion Supporting Reads

**c. Per-Specimen Precision for SNVs and Indels**

The precision was calculated for each individual specimen as shown in Table 8. Results from the precision studies were combined across all reportable genes for each specimen. The positive and negative call rates for sequenced mutations (SNVs and indels) in each sample were calculated based on the total number of mutations along with the two-sided 95% CI.

**Table 8. Positive and Negative Call Rates per Sample**

Sample ID	Cancer Type	Unique Mutations Detected*	PCR Per Mutation	PCR (%) (n/N) (95% CI)	NCR (%) (n/N) (95% CI)
Sample 1	Ovary/ Fallopian Tube	4	45/48 for 1 <sup>***</sup>	98.43	100
			48/48 for 3	189/192 (95.5, 99.5)	42672/42672 (99.99, 100)
Sample 2	Thyroid	1	48/48 for 1	100 48/48 (92.6, 100.0)	100 42672/42672 (99.99, 100)
Sample 3	Prostate	1	48/48 for 1	100 48/48 (92.6, 100.0)	100 42672/42672 (99.99, 100)
Sample 4	Bowel	6	25/31 for 1 <sup>**</sup>	94.09, 175/186, (88.0, 95.5)	100, 42672/42672, (99.99,100)
			26/31 for 1 <sup>**</sup>		
			31/31 for 4 <sup>**</sup>		

\* Across the 48 Replicates

\*\* Reduction in the number of replicates due to QC failure

\*\*\* Variant with VAF at LoD

Sample ID	Cancer Type	Unique Mutations Detected*	PCR Per Mutation	PCR (%) (n/N) (95% CI)	NCR (%) (n/N) (95% CI)
Sample 5	Lung	1	46/48 for 1 <sup>***</sup>	95.83, 46/48, (85.8, 99.5)	100, 42720/42720, (100, 100)
Sample 6	Bowel	16	42/48 for 2 <sup>***</sup>	98.04, 753/768, (96.8, 98.8)	99.99, 42522/42528, (99.97,99.99)
			45/48 for 1 <sup>***</sup>		
			48/48 for 13		
Sample 7	Lung	1	47/48 for 1	97.91, 47/48, (89.1, 99.6)	99.91, 42634/42672, (99.88,99.94)
Sample 8	Prostate	2	46/46 for 2 <sup>**</sup>	100, 92/92, (96.0, 100)	100, 40940/40940, (99.99,100)
Sample 9	Lung	6	48/48 for 6	100, 288/288, (98.7, 100)	100, 42623/42624, (99.99,100)
Sample 10	Bowel	12	45/48 for 1 <sup>***</sup>	99.3, 572/576, (98.7, 100)	100, 42432/42432, (99.99,100)
			47/48 for 1 <sup>***</sup>		
			48/48 for 10		
Sample 11	CNS/ Brain	4	48/48 for 4	100, 192/192, (98.0, 100)	99.94, 42646/42672, (99.91,99.96)
Sample 12	Bowel	3	44/46 for 1 <sup>**,***</sup>	97.1, 134/138, (98.0, 100)	100, 40801/40802, (99.99,100)
			44/46 for 1 <sup>**</sup>		
			46/46 for 1 <sup>**</sup>		
Sample 13	Breast	6	47/48 for 1 <sup>***</sup>	99.65, 287/288, (98.1, 99.9)	100, 42623/42624, (99.99, 100)
			48/48 for 5		

\* Across the 48 Replicates

\*\* Reduction in the number of replicates due to QC failure

\*\*\* Variant with VAF at LoD

**d. Precision for *ERBB2* Amplification**

Precision of *ERBB2* amplification was evaluated using 15 samples. Three out of the 15 samples were tested at two DNA input levels using two different operators on two instruments across three non-consecutive days using two kit lots, yielding 12 replicates per sample per input level. The remaining samples were tested in duplicate by two different operators on six instruments across four non-consecutive days using two kit lots, yielding 24 replicates per sample per input level. Of the 15 samples, four contained *ERBB2* amplifications and the remaining samples contained no *ERBB2* amplifications. The CellDx-Tissue reports an *ERBB2* amplification when the observed copy number for the gene is determined by the test to be 8.5 copies or more. The tumor purity, input, mean observed copy number (CN), coefficient of variation (%CV), positive and negative call rates along with the 95% CIs for each sample are summarized in Table 9. The call rates were 100% for both amplification and no amplification groups.

**Table 9. Summary of the *ERBB2* Gene Amplification Precision**

Sample ID	Cancer Type	<i>ERBB2</i> Status	Tumor Purity (%)	Input (ng)	Mean CN	CN CV (%)	PCR % (n/N) (95% CI)	NCR % (n/N) (95% CI)
Sample 1	Breast	Positive	58	10	16.3	11.0	100 (24/24) (85.7, 100)	-
				20	16.4	7.0	100 (24/24) (85.7, 100)	-
Sample 2	Bowel	Positive	21	10	8.06	6.3	100 (12/12) (75.8, 100)	-
				20	7.82	4.9	100 (12/12) (75.8, 100)	-
Sample 3	Esophagus	Positive	49	10	15.57	10.2	100 (12/12) (75.8, 100)	-
				20	15.34	7.7	100 (12/12) (75.8, 100)	-
Sample 4	Ovarian / Fallopian Tube	Positive	55	10	14.54	7.1	100 (12/12) (75.8, 100)	-
				20	14.02	7.0	100 (12/12) (75.8, 100)	-

Sample ID	Cancer Type	ERBB2 Status	Tumor Purity (%)	Input (ng)	Mean CN	CN CV (%)	PCR % (n/N) (95% CI)	NCR % (n/N) (95% CI)
Sample 5	CNS	Negative	92	10	1.98	5.1	-	100 (24/24) (86.2, 100)
				20	1.98	4.6	-	100 (24/24) (86.2, 100)
Sample 6	Bowel	Negative	37	10	1.34	15.2	-	100 (24/24) (86.2, 100)
				20	1.32	14.9	-	100 (24/24) (86.2, 100)
Sample 7	Lung	Negative	28	10	1.97	8.8	-	100 (24/24) (86.2, 100)
				20	1.93	8.8	-	100 (24/24) (86.2, 100)
Sample 8	Prostate	Negative	22	10	1.79	18.6	-	100 (24/24) (86.2, 100)
				20	1.75	13.7	-	100 (22/22) (85.1, 100)
Sample 9	Lung	Negative	20	10	1.28	15.1	-	100 (24/24) (86.2, 100)
				20	1.28	16.4	-	100 (24/24) (86.2, 100)
Sample 10	Bowel	Negative	29	10	1.54	5.7	-	100 (24/24) (86.2, 100)
				20	1.54	10.3	-	100 (24/24) (86.2, 100)
Sample 11	Lung	Negative	45	10	1.58	11.8	-	100 (24/24) (86.2, 100)
				20	1.62	9.2	-	100 (24/24) (86.2, 100)

Sample ID	Cancer Type	<i>ERBB2</i> Status	Tumor Purity (%)	Input (ng)	Mean CN	CN CV (%)	PCR % (n/N) (95% CI)	NCR % (n/N) (95% CI)
Sample 12	Bowel	Negative	39	10	2.87	24.2	-	100 (7/7) (59.0, 100)
				20	2.70	8.0	-	100 (24/24) (86.2, 100)
Sample 13	Prostate	Negative	36	10	1.98	6.6	-	100 (24/24) (86.2, 100)
				20	2.08	5.7	-	100 (24/24) (86.2, 100)
Sample 14	Thyroid	Negative	22	10	1.99	2.9	-	100 (24/24) (86.2, 100)
				20	1.99	3.3	-	100 (24/24) (86.2, 100)
Sample 15	Ovary	Negative	56	10	1.52	13.7	-	100 (24/24) (86.2, 100)
				20	1.58	9.1	-	100 (24/24) (86.2, 100)

*e. Precision for ALK, RET, and ROS1 Fusions*

The precision of RNA fusions was evaluated across different samples with varying tumor purities. The results are summarized in Table 10.

**Table 10. Summary of *ALK*, *RET*, and *ROS1* Fusions**

Fusion	Tumor Purity (%)	LoD Level	Mean Supporting Reads (Range)	PCR (%) n/N (95% CI)	NCR (%) n/N (95% CI)
<i>ALK</i>	20	1x*	15062 (7121-30568)	100 45/45 (92.1, 100)	100 609/609 (99.4, 100)

<b>Fusion</b>	<b>Tumor Purity (%)</b>	<b>LoD Level</b>	<b>Mean Supporting Reads (Range)</b>	<b>PCR (%) n/N (95% CI)</b>	<b>NCR (%) n/N (95% CI)</b>
<i>RET</i>	25	1.8x*	24438 (740-175619)	100 46/46 (92.3, 100)	100 608/608 (99.4, 100)
<i>ROS1</i>	45	0.8x**	1273 (504-2328)	95 38/40 (90.8, 100)	100 614/614 (99.4, 100)

**f. Analysis of Source of Variance**

The Average Positive Agreement (APA) and Average Negative Agreement (ANA) was assessed to analyze the imprecision caused by different sources of variance. Data analysis is presented stratified by variant type and present for 1) overall, 2) instrument to instrument, 3) operator to operator, and 4) day to day (Table 11).

**Table 11. Reproducibility of the CellDx – Tissue Assay**

<b>Alteration Type</b>	<b>Metric</b>	<b>Overall (%) (95% CI)</b>	<b>Inter-Instrument (%) (95% CI)</b>	<b>Inter-Operator (%) (95% CI)</b>	<b>Inter-Day (%) (95% CI)</b>
Variants with Evidence of Clinical Significance	APA	99.56 (98.42, 99.88)	99.56 (95.19, 100)	99.57 (98.34, 100)	99.57 (96.74, 100)
	ANA	100 (100, 100)	100 (100, 100)	100 (100, 100)	100 (100, 100)
Hotspot Variants	APA	97.99 (96.98, 98.67)	97.94 (95.27, 99.44)	97.95 (96.21, 98.74)	97.99 (95.31, 98.99)
	ANA	100 (100, 100)	100 (100, 100)	100 (100, 100)	100 (100, 100)
Non-hotspot Variants	APA	99.12 (98.57, 99.46)	99.11 (97.62, 99.82)	99.11 (98.27, 99.55)	99.12 (97.76, 99.66)
	ANA	100 (100, 100)	100 (100, 100)	100 (100, 100)	100 (100, 100)
SNVs	APA	99.13 (98.65, 99.44)	99.11 (97.61, 99.72)	99.12 (98.33, 99.50)	99.13 (97.88, 99.61)
	ANA	100 (100, 100)	100 (100, 100)	100 (100, 100)	100 (100, 100)

<b>Alteration Type</b>	<b>Metric</b>	<b>Overall (%) (95% CI)</b>	<b>Inter-Instrument (%) (95% CI)</b>	<b>Inter-Operator (%) (95% CI)</b>	<b>Inter-Day (%) (95% CI)</b>
Insertions (All)	APA	100 (97.40, 100)	100 (86.20, 100)	100 (94.93, 100)	100 (90.36, 100)
	ANA	100 (100, 100)	100 (100, 100)	100 (100, 100)	100 (100, 100)
Insertions (1 – 5 bp)	APA	100 (97.40, 100)	100 (86.20, 100)	100 (94.93, 100)	100 (90.36, 100)
	ANA	100 (100, 100)	100 (100, 100)	100 (100, 100)	100 (100, 100)
Insertions (5 – 10 bp)	APA	N/A	N/A	N/A	N/A
	ANA	100 (100, 100)	100 (100, 100)	100 (100, 100)	100 (100, 100)
Insertions (11 – 30 bp)	APA	N/A	N/A	N/A	N/A
	ANA	100 (100, 100)	100 (100, 100)	100 (100, 100)	100 (100, 100)
Deletions (All)	APA	96.69 (94.89, 97.87)	96.69 (89.77, 98.37)	96.69 (93.71, 98.10)	96.69 (92.13, 98.51)
	ANA	100 (100, 100)	100 (100, 100)	100 (100, 100)	100 (100, 100)
Deletions (1 – 5 bp)	APA	96.39 (94.43, 97.68)	96.38 (88.89, 98.22)	96.38 (93.14, 97.92)	96.39 (91.44, 98.37)
	ANA	100 (100, 100)	100 (100, 100)	100 (100, 100)	100 (100, 100)
Deletions (5 – 10 bp)	APA	N/A	N/A	N/A	N/A
	ANA	100 (100, 100)	100 (100, 100)	100 (100, 100)	100 (100, 100)
Deletions (11 – 30 bp)	APA	100 (92.59, 100)	100 (67.56, 100)	100 (86.20, 100)	100 (75.75, 100)
	ANA	100 (100, 100)	100 (100, 100)	100 (100, 100)	100 (100, 100)

N/A: Not Available / Applicable

Alteration Type	Metric	Overall (%) (95% CI)	Inter-Instrument (%) (95% CI)	Inter-Operator (%) (95% CI)	Inter-Day (%) (95% CI)
<i>ERBB2</i> Amplification	APA	100 (92.59, 100)	100 (67.56, 100)	100 (86.20, 100)	100 (75.75, 100)
	ANA	100 (99.25, 100)	100 (95.68, 100)	100 (98.51, 100)	100 (97.06, 100)
<i>ALK</i> Fusion	APA	100 (92.13, 100)	100 (60.97, 100)	100 (79.61, 100)	100 (74.12, 100)
	ANA	100 (99.20, 100)	100 (94.65, 100)	100 (97.64, 100)	100 (96.87, 100)
<i>RET</i> Fusion	APA	100 (92.29, 100)	100 (64.57, 100)	100 (79.61, 100)	100 (75.75, 100)
	ANA	100 (99.20, 100)	100 (94.65, 100)	100 (97.64, 100)	100 (96.87, 100)
<i>ROS1</i> Fusion	APA	95.00 (83.50, 98.62)	91.67 (43.65, 96.99)	95.83 (77.19, 100)	95.00 (72.25, 100)
	ANA	100 (99.20, 100)	100 (94.65, 100)	100 (97.64, 100)	100 (96.87, 100)

N/A: Not Available / Applicable

**g. Lot-to-Lot Precision**

The performance of the CellDx – Tissue assay was assessed across two unique kit lots by determining the concordance of variant calls in 17 FFPE tissue samples. The two unique lots were used to process 372 test cases in duplicate for a total of 744 observations. All batches were sequenced on the same instrument. Table 12. lists the APA and ANA used to assess lot to lot performance. The overall panel-wide variants (SNVs, Insertions and Deletions) APA is 98.68% and the ANA is 100%. *ERBB2* amplification, *ALK*, and *RET* fusions had a concordance of 100%; and *ROS1* fusions demonstrated a concordance of 93.75%.

**Table 12. Lot-to-Lot Precision of CellDx – Tissue Assay**

Variant Type	Performance	Between Lot 1& Lot 2 (%) (95% CI)
Variants with Evidence of Clinical Significance	APA	99.57 (98.34, 100)
	ANA	100 (100, 100)
Panel-Wide (SNVs + Insertions + Deletions)	APA	98.68 (97.97, 99.16)
	ANA	100 (100, 100)
Hotspot SNVs (Including Clinically Significant Variants)	APA	97.95 (96.21, 98.74)
	ANA	100 (100, 100)
Non-hotspot SNVs (Including Clinically Significant Variants)	APA	99.11 (98.27, 99.55)
	ANA	100 (100, 100)
SNVs (Hotspot + Non-hotspot)	APA	99.12 (98.33, 99.50)
	ANA	100 (100, 100)
Insertions	APA	100 (94.93, 100)
	ANA	100 (100, 100)
Deletions	APA	96.69 (93.71, 98.10)
	ANA	100 (100, 100)
<i>ERBB2</i> Amplification	APA	100 (86.20, 100)
	ANA	100 (98.51, 100)

Variant Type	Performance	Between Lot 1& Lot 2 (%) (95% CI)
<i>ALK</i> Fusion	APA	100 (85.13, 100)
	ANA	100 (98.41, 100)
<i>RET</i> Fusion	APA	100 (85.69, 100)
	ANA	100 (98.41, 100)
<i>ROS1</i> Fusion	APA	93.75 (76.39, 99.11)
	ANA	100 (98.41, 100)

### 3. **Analytical Sensitivity – Limit of Detection (LoD)**

The LoD of the CellDx – Tissue assay for SNVs and indels is defined as the lowest concentration with  $\geq 95\%$  of replicates for a variant are reliably detected. The LoD of the CellDx – Tissue assay for *ERBB2* amplification is determined as the minimum tumor purity required for robust reporting of amplification status and incorporated the clinical reporting framework where  $CN \geq 8.5$  is Positive,  $CN \geq 4$  to  $< 8.5$  is Indeterminate / Low-Level Amplification (below LoD), and  $CN < 4$  is Negative. The LoD of the CellDx – tissue assay for *ALK*, *RET*, and *ROS1* fusions is determined as the minimum tumor purity required for robust reporting of fusion status. The recommended nucleic acid input for the CellDx – Tissue assay is 20 ng of DNA and 20 ng of RNA recovered from tissue with a minimum of 20% viable tumor nuclei. Details of the data are discussed and shown below.

#### ***a. LoD – SNVs, Insertions, and Deletions***

The analytical sensitivity of the CellDx – Tissue assay for SNVs, insertions, and deletions was evaluated by assessing 10 clinical FFPE specimens from seven different cancer types with 29 SNVs, 1 insertion, and 2 deletions. Each sample was diluted to at least five dilution levels and tested in 10 replicates per level using two reagent lots. The call rate was determined for each variant and the LoD was approximated between the call rate that was below 95% and the highest call rate (100%). A summary of the estimation of the LoD range for a set of representative variants is shown in Table 13. The established VAF range for each variant type (Hotspot SNVs, Non-hotspot SNVs, insertions, and deletions) is shown in Table 14.

**Table 13. Estimation of LoD Range for Representative Variants**

Tumor Type	Variant Type	Level	Gene	AA Change	Avg DP	Avg AD	VAF Range	Mean VAF (%)	PCR (%)
Lung	SNV	3	<i>TP53</i>	p.R175H	1869	68	0.026-0.048	3.7	100
Lung	SNV	5	<i>TET2</i>	c.3954+3T>G	1426	63	0.024-0.098	4.3	100
Lung	DEL	5	<i>EGFR</i>	p.E746_A750del	1939	69	0.026-0.058	3.5	100
Breast	SNV	2	<i>GPS2</i>	c.318-6C>T	1886	123	0.056-0.078	6.5	100
Breast	SNV	4	<i>PIK3CA</i>	p.H1047R	1678	58	0.022-0.049	3.5	100
Uterus	SNV	2	<i>TP53</i>	p.R273C	1955	75	0.024-0.053	3.8	100
Uterus	SNV	2	<i>CTLA4</i>	p.A54T	2000	132	0.054-0.078	6.6	100
Uterus	SNV	3	<i>CTNNB1</i>	p.G34E	1959	74	0.026-0.045	3.7	100
Head and Neck	SNV	2	<i>NRAS</i>	p.A59T	2000	107	0.044-0.071	5.4	100
Head and Neck	SNV	3	<i>DICER1</i>	c.3269+6C>T	1921	128	0.045-0.13	6.8	100
Head and Neck	SNV	4	<i>TP53</i>	p.G262V	1885	75	0.026-0.047	3.9	100
Head and Neck	SNV	3	<i>NF2</i>	p.R196*	659	48	0.059-0.104	6.7	100
Head and Neck	SNV	4	<i>EP300</i>	c.3728+5G>A	1271	56	0.035-0.052	4.0	100
Head and Neck	SNV	4	<i>FAT1</i>	p.R1205*	1919	118	0.05-0.078	6.1	100

Tumor Type	Variant Type	Level	Gene	AA Change	Avg DP	Avg AD	VAF Range	Mean VAF (%)	PCR (%)
Head and Neck	SNV	4	<i>LATS1</i>	p.Q953*	1603	56	0.028-0.047	3.4	100
Head and Neck	SNV	3	<i>PHF6</i>	p.Y301*	938	52	0.041-0.077	5.7	100
Bowel	SNV	2	<i>TP53</i>	p.R273H	1925	127	0.052-0.077	6.6	100
Esophagus / Stomach	INS	4	<i>ARID1A</i>	p.A1136Gfs*57	1503	92	0.042-0.088	6.2	100
Esophagus / Stomach	SNV	3	<i>TP53</i>	p.R209*	1924	165	0.064-0.16	8.5	100
Esophagus / Stomach	SNV	3	<i>POT1</i>	p.E456*	1281	71	0.032-0.107	5.5	100
Lung	SNV	1	<i>TP53</i>	p.R248W	1797	124	0.043-0.094	6.9	100
Uterus	SNV	2	<i>FANCA</i>	c.2015-8C>T	1765	101	0.04-0.066	5.7	100
Uterus	SNV	2	<i>TP53</i>	p.Y220C	1287	85	0.055-0.083	6.6	100
Uterus	SNV	2	<i>PPP2R1A</i>	p.P179R	1708	57	0.022-0.05	3.4	100
Prostate	SNV	5	<i>MLH1</i>	p.Q346H	1781	82	0.033-0.056	4.1	100
Prostate	SNV	5	<i>AR</i>	p.L702H	1579	83	0.044-0.072	5.4	100
Lung	INS	2	<i>ARID1A</i>	p.P1468Lfs*13	938	88	0.076-0.122	9.5	100
Lung	SNV	5	<i>TCF7L2</i>	c.1162-12C>T	1720	134	0.039-0.144	8.0	100
Lung	SNV	3	<i>TP53</i>	p.Y126*	1911	151	0.029-0.154	7.2	100
Lung	SNV	4	<i>FANCE</i>	c.855+6T>A	1029	64	0.026-0.095	5.2	100
Lung	SNV	2	<i>KMT2C</i>	p.E1333*	891	105	0.045-0.214	11.5	100
Lung	SNV	3	<i>KMT2D</i>	c.840-10G>A	1861	147	0.025-0.153	7.3	100

**Table 14. Analysis Sensitivity (LoD VAF) for Representative SNVs and Indels**

Variant	Established VAF Range (%)	Tumor Purity (%)	Average Allele Depth	Number of Variants in Clinical Cases in the Established Range
Hotspot SNVs	2.2 – 9.4	7	57 – 124	10
Non-hotspot SNVs	2.4 – 21.4	4 – 24	63 – 105	19
Insertions	4.2 – 8.8	6	92	1
Deletions	2.6 – 12.2	4 – 24	69 – 88	2
Homopolymer Context	4.2 – 8.8	6	92	1
Non-Homopolymer Context	2.6 – 12.2	4 – 24	69 – 88	2

**b. LoD – Amplifications and Fusions (Tumor Purity)**

**i. ERBB2 Amplification**

The analytical sensitivity of *ERBB2* was assessed by testing one clinical breast cancer sample diluted with normal FFPE DNA to five detection levels. Each dilution level was run in 10 replicates per level per reagent lot (two lots), resulting in a total of 50 observations. The cutoff for *ERBB2* detection is an observed copy number  $\geq 8.5$ . The data summarized in Table 15 showed 100% call rates for samples with a tumor purity of 15%.

**Table 15. Analysis Sensitivity for *ERBB2* Amplification by the CellDx-Tissue Assay**

Dilution Level	Tumor Purity (%)	Mean Observed Copy Number (Mean $\pm$ SD)	Observed Copy Number Range	Positive Call Rate (%) (n/N)
1	60	30.2 $\pm$ 0.8	29–32	100 (10/10)
2	30	15.0 $\pm$ 0.7	14–16	100 (10/10)
3	15	8.5 $\pm$ 0.2	8–11	100 (10/10)

Dilution Level	Tumor Purity (%)	Mean Observed Copy Number (Mean ± SD)	Observed Copy Number Range	Positive Call Rate (%) (n/N)
4	8	5.1 ± 0.9	4-7	20 (2/10)
5	4	3.15 ± 0.1	3-3	0 (0/10)

ii. *ALK*, *RET*, and *ROS1* RNA Fusions

The analytical sensitivity of *ALK*, *RET*, and *ROS1* was assessed by testing four clinical (three lung and one bowel) cancer samples. *ALK* and *RET* clinical samples were diluted with normal FFPE DNA to five dilution levels. *ROS1* clinical samples were diluted to nine dilution levels. Each dilution level was run in 10 replicates per level (across two lots), resulting in a total of 50 observations for *ALK* and *RET* and 90 observations for *ROS1*. The data is summarized in Table 16 – Table 19.

**Table 16. Analysis Sensitivity for *ALK* Amplification by the CellDx-Tissue Assay**

Dilution Level	Tumor Purity (%)	Fusion Read Counts	Mean Fusion Read Range	Positive Call Rate (%) (n/N)
1	20	4868	1640-8989	100 (10/10)
2	15	1048	603-1544	80 (8/10)
3	11	1104	734-1506	30 (3/10)
4	8	1014	642-1409	30 (3/10)
5	6	931	740-1049	30 (3/10)

**Table 17. Analysis Sensitivity for *RET* Amplification by the CellDx-Tissue Assay**

Dilution Level	Tumor Purity (%)	Fusion Read Counts	Mean Fusion Read Range	Positive Call Rate (%) (n/N)
1	25	11534	5788-19474	100 (10/10)
2	19	7321	3996-10730	100 (7/7)
3	14	4282	1366-8416	100 (10/10)
4	11	7604	2120-16136	100 (7/7)

Dilution Level	Tumor Purity (%)	Fusion Read Counts	Mean Fusion Read Range	Positive Call Rate (%) (n/N)
5	8	6095	600-18193	86 (6/7)

**Table 18. Analysis Sensitivity for *ROS1* Amplification by the CellDx-Tissue Assay**

Dilution Level	Tumor Purity (%)	Fusion Read Counts	Mean Fusion Read Range	Positive Call Rate (%) (n/N)
1	70	80775	28692-130814	100 (10/10)
2	53	41636	28054-58784	100 (10/10)
3	39	24685	5883-39701	100 (10/10)
4	30	46067	9422-123494	100 (10/10)
5	22	9000	5861-14743	100 (10/10)
6	12	11498	4925-33280	100 (9/9)
7	9	5692	1785-13236	100 (10/10)
8	7	4393	1015-13653	100 (10/10)
9	5	4346	1605-7058	100 (10/10)
10*	20	39,875	34,176 - 49,210	100 (5/5)
11*	10	30,333	22,285 - 42,599	100 (5/10)
12*	7	1,853	1,158 - 2,641	100 (5/5)
13*	5	1,599	892 - 3,587	100 (5/5)
14*	4	942	715 - 1,431	70 (7/10)
15*	3	883	554 - 1,141	30 (3/10)
16*	2	-	-	0 (0/0)

\*Samples and dilution levels are from a Supplemental Study to determine the tumor purity LoD for *ROS1*.

**Table 19. Analysis Sensitivity Summary of RNA Fusion**

<b>Gene</b>	<b>Tumor Purity Range (%)</b>	<b>Fusion Read Count Range</b>	<b>LoD Tumor Purity (%)</b>	<b>LoD Supporting Reads</b>
<i>ALK</i>	20 – 6	4868-931	20	>1640
<i>RET</i>	25 – 8	11533–6095	14	>1366
<i>ROSI*</i>	20 – 2	39875-883	5	>892

\*Supplemental LoD Study

**4. Linearity/Assay Reportable Range**

Not applicable

**5. Traceability (Controls, Calibrators, or Methods)**

**a. *Traceability***

The CellDx-Tissue assay is not traceable to any known standard. Controls and quality metrics are described in the device description section.

**b. *Stability***

Reagent stability is based on the manufacturer expiration dating, verified by Datar Cancer Genetics. Stability of the reagents is monitored through the use of consistent controls.

**c. *Expected Values (Controls, Calibrators, or Methods)***

The external positive control (commercial vendor qualified by Datar Cancer Genetics) contains different confirmed mutations, representing a range of mutation allele frequencies. A negative external control (commercial vendor qualified by Datar Cancer Genetics) contains a non-cancerous cell line derived DNA /RNA with no variants of interest. Both external controls are processed from library preparation through sequencing to serve as an end-to-end control to demonstrate assay performance. Failure of either external control to meet the quality control metrics will result in all test samples on the run be flagged and re-sequenced.

**6. Analytical Specificity**

**a. *Cut-Off/ False Positive Range (Limit of Blank)***

Non-cancerous FFPE tissues were evaluated for analytical specificity to assess the risk of false positives in normal tissue when detecting SNVs, indels, amplifications, and fusions using the CellDx – Tissue assay. A total of 37 normal or benign-adjacent tissue were processed across two reagent lots, multiple operators, and days. Ninety (90) false positive events were detected for SNVs and indels for a false positive rate (FPR) of 0.000054 (95 CI: 0.000043 - 0.000066). None of the detected variants were classified as Variants with Evidence of Clinical Significance (Level 2), while four were classified as

Variants with Potential Evidence of Clinical Significance (Level3). The FPR was determined to be 0 (95 CI: 0 - 3.09) for *ERBB2* amplifications, and *ALK*, *RET*, and *ROSI* fusion genes.

**b. Necrotic Tissue and Tumor Block Age**

The impact of necrosis on the performance of the CellDx – Tissue assay was evaluated by assessing the invalid rates for both DNA and RNA libraries, defined as the proportion of libraries that failed to meet the predefined quality control thresholds. All 378 FFPE DNA and 377 FFPE RNA samples were evaluated for necrosis over a range of 0 - 80 % and invalid rates were examined. The assay performance was assessed and compared to the predicate device. The DNA and RNA invalid rates ranged from 0 to 10.52% and 0 to 1.78%, respectively (Table 20). The PPA for overall variant (SNVs, insertions, and deletions) detected ranged from 97.8% to 100%, while an NPA of 99.5% was observed (Table 21).

**Table 20. Necrotic Tissue Invalid Rate for the CellDx-Tissue Assay**

Cancer Types	Sample Count	Necrotic Tissue Content (%)	Invalid Rate (%) (n/N)
DNA			
23	281	0 – 5	1.78 (5/281)
15	72	5 – 20	0 (0/72)
12	19	21 – 40	10.52 (2/19)
4	6	41 – 80	0 (0/6)
RNA			
23	281	0 – 5	1.42 (4/281)
15	71	5 – 20	1.40 (1/71)
12	19	21 – 40	0 (0/19)
4	6	41 – 80	0 (0/6)

**Table 21. Concordance of Overall Variant (SNV + Insertions + Deletions) Detection**

Necrotic Tissue Content (%)	Sample Count	PPA (%) (95 CI) (n/N)	NPA (%) (95 CI) (n/N)
0 – 5	86	97.8 (93.8, 99.3) (135/138)	99.6 (99.4, 99.7) (13140/13197)

<b>Necrotic Tissue Content (%)</b>	<b>Sample Count</b>	<b>PPA (%) (95 CI) (n/N)</b>	<b>NPA (%) (95 CI) (n/N)</b>
>5 – 20	28	100 (91.0, 100) (39/39)	99.2 (98.9, 99.4) (4268/4303)
21 – 40	8	100 (77.2, 100) (13/13)	99.8 (99.3, 99.9) (1224/1227)
41 – 53	1	100 (34.2, 100) (2/2)	99.3 (96.4, 99.9) (153/153)
Overall	123	98.4 (95.5, 99.5) (189/192)	99.5 (99.4, 99.6) (18784,18880)

**c. Interfering Substances**

The impact of interfering substances on the performance of the CellDx – Tissue assay was assessed by processing DNA and RNA from FFPE samples tested in the presence of each interfering substance at varying amounts (Table 22). The samples were evaluated for concordance of variant call when compared to samples processed without the interfering substances. Replicates for six test cases were analyzed for six experimental and one baseline condition. Performance was evaluated across nine samples x seven conditions x five replicates. Analysis of all variant types tested (SNVs, indels, amplifications, and fusions) showed no effect of exogenous and endogenous interferent for all conditions (Table 23 – Table 25).

**Table 22. Endogenous and Exogenous Interfering Substances Tested**

<b>Type</b>	<b>Substance</b>	<b>Concentration Spiked</b>	
		<b>At Extraction</b>	<b>At Library Preparation</b>
Exogenous	Ethanol	3x	5%
Exogenous	Index Adaptor	-	3x
Exogenous	Proteinase K	3x	0.04 mg/mL
Exogenous	Wash Buffer	3x	5%
Endogenous	Melanin	0.2 mg/mL	-
Endogenous	Hemoglobin	2 mg/mL	-

**Table 23. Interfering Substances Concordance by Test Condition for SNVs and Indels**

<b>Test Condition</b>	<b>PCR (%) (n/N) (95% CI)</b>	<b>NCR (%) (n/N) (95% CI)</b>
Ethanol	99.2 (119/120) (95.4, 99.9)	100 (106080/106080) (99.99, 100)
Adaptor	98.8 (115/120) (90.6, 98.2)	100 (106080/106080) (99.99, 100)
Proteinase K	94.9 (112/118) (89.3, 97.6)	100 (104312/104312) (99.99, 100)
Wash Buffer	98.3 (118/120) (94.1, 99.5)	100 (106080/106080) (99.99, 100)
Melanin	100 (110/110) (96.6, 100.0)	100 (97240/97240) (99.99, 100)
Hemoglobin	98.3 (118/120) (94.1, 99.5)	100 (106080/106080) (99.99, 100)

**Table 24. Interfering Substances Concordance by Test Condition for *ERBB2***

<b>Test Condition</b>	<b>PCR (%) (n/N) (95% CI)</b>	<b>NCR (%) (n/N) (95% CI)</b>
Ethanol	100 (10/10) (69.2, 100.0)	100 (35/35) (90.1, 100.0)
Adaptor	100 (10/10) (69.2, 100.0)	100 (35/35) (90.1, 100.0)
Proteinase K	100 (10/10) (69.2, 100.0)	100 (35/35) (90.1, 100.0)
Wash Buffer	100 (10/10) (69.2, 100.0)	100 (35/35) (90.1, 100.0)
Melanin	100 (10/10) (69.2, 100.0)	100 (30/30) (88.6, 100.0)

Test Condition	PCR (%) (n/N) (95% CI)	NCR (%) (n/N) (95% CI)
Hemoglobin	100 (10/10) (69.2, 100.0)	100 (35/35) (90.1, 100.0)

**Table 25. Interfering Substances Concordance by Test Condition for RNA Fusions (*ALK*, *RET*, and *ROS1*)**

Test Condition	PCR (%) (n/N) (95% CI)	NCR (%) (n/N) (95% CI)
Ethanol	100 (15/15) (78.2, 100)	100 (30/30) (88.4, 100)
Adaptor	100 (15/15) (78.2, 100)	100 (30/30) (88.4, 100)
Proteinase K	100 (15/15) (78.2, 100)	100 (29/29) (88.1, 100)
Wash Buffer	100 (15/15) (78.2, 100)	100 (30/30) (88.4, 100)
Melanin	100 (15/15) (78.2, 100)	100 (30/30) (88.4, 100)
Hemoglobin	100 (15/15) (78.2, 100)	100 (30/30) (88.4, 100)

**d. Sample Carryover and Cross-Contamination**

Cross-contamination (contamination from one sample to another within the same batch) and sample carryover (contamination from a previous sequencing run when using the same instrument) were assessed by evaluating the false positive rate in 24 FFPE samples. Twelve of the 24 FFPE samples had known positive variants, and the remaining samples were known negative samples. All FFPE samples were assessed across two batches to test for contamination within and between runs. The samples were processed in a checkerboard layout during library synthesis. In the first batch, 24 libraries and 12 libraries underwent a single DNA and RNA sequencing, respectively. The second batch contained 12 DNA libraries and six RNA libraries with known negative samples were processed after the first batch. No positive variant results were observed in the known negative samples tested.

## 7. Robustness Study

### a. *Sample Stability*

The stability of tissue curls sectioned from FFPE blocks were established using six FFPE tumor specimens stored under appropriate conditions (17.8°C – 26.8°C), for 35 days by comparing the detection of the intended variants within these samples to the baseline timepoint (T0). Three sets of curls from each block were prepared on Day 0. The first set of curls was immediately processed for DNA/RNA extraction and subsequent analysis to establish the T0 baseline. On Day 35 (T35), the remaining stored curls were used for DNA/RNA extraction and analysis for each claimed variant. FFPE curls stored for up to 35 days show concordant variant calls for all variants (Table 26) demonstrating a stability of 28 days for FFPE blocks stored under appropriate conditions.

**Table 26. FFPE Tissue Curl Stability Concordance by Variant**

Variant Type	PPA (%) (n/N) (95% CI)		NPA (%) (n/N) (95% CI)	
	T0	T35	T0	T35
SNVs + Indels	99.32 (147/148) (96.27, 99.88)	96.96 (278/296) (94.14, 98.47)	100 (884/884) (99.58, 100)	99.55 (880/884) (98.78, 99.85)
<i>ERBB2</i> Amplification	100 (4/4) (51.01, 100)	100 (8/8) (67.56, 100)	100 (20/20) (83.89, 100)	100 (40/40) (91.24, 100)
<i>ALK</i> and <i>ROS1</i> Fusions	100 (8/8) (67.56, 100)	85.7 (12/14) (60.1, 96.0)	100 (64/64) (94.3, 100)	98.4 (126/128) (94.5, 99.6)
RET Fusion	100 (8/8) (67.56, 100)	100 (16/16) (80.6, 100)	100 (8/8) (67.56, 100)	100 (16/16) (80.6, 100)

### b. *DNA Input*

The recommended amount of input DNA for the CellDx-Tissue assay is 20 ng. A DNA input study was performed to demonstrate the analytical performance of the CellDx-Tissue assay across a range of DNA inputs (2.5, 5, 10, 20, 30,40, and 50 ng) using three FFPE tumor samples panning multiple variant types including SNVs, insertions, deletions, and *ERBB2* amplification. Each sample was tested across 12 replicates, two operators, two lot kits, and seven instruments across multiple non-consecutive days. The QC metrics and variant reporting for each sample was evaluated per DNA input level. Table 27 and Table 28 list the PCR and NCR for each input level were analyzed, including SNVs, indels, and *ERBB2* amplifications.

**Table 27. Panel-wide (SNVs + Indels) Call Rates at each DNA Input Level**

DNA Input (ng)	PCR (%) (n/N) (95% CI)	NCR (%) (n/N) (95% CI)
2.5	85.8 (151/176) (77.8, 88.5)	99.9 (45812769/45813184) (99.9, 100)
5	80.1 (137/171) (69.4, 81.8)	99.9 (38654612/38654874) (99.9, 100)
10	97.8 (175/179) (93.7, 98.8)	99.9 (50108021/51539832) (99.9, 100)
20	97.2 (147/179) (92.9, 98.5)	99.9 (50108137/50108170) (99.9, 100)
30	97.8 (175/179) (93.7, 98.8)	99.9 (50108155/50108170) (99.9, 100)
40	98.3 (176/179) (94.4, 99.1)	99.9 (50108167/50108170) (99.9, 100)
50	100 (180/180) (97.9, 100)	99.9 (51539823/51539832) (99.9, 100)

**Table 28. ERBB2 Call Rates at each DNA Input Level**

DNA Input (ng)	Tumor Purity (%)	Mean CN (CN Range)	SD (±)	CV (%)	PCR (%) (n/N) (95% CI)
2.5	21	9.1 (8–11)	1.0	11.5	92 (11/12) (61.5, 99.8)
5	21	9.0 (8–10)	0.7	7.4	83 (10/12) (51.6, 97.9)
10	21	7.4 (6–9)	0.8	10.7	100 (12/12) (73.5, 100)
20	21	7.5 (6–9)	0.9	12.1	100 (12/12) (73.5, 100)
30	21	8.3 (7–10)	1.1	13.7	100 (12/12) (73.5, 100)

DNA Input (ng)	Tumor Purity (%)	Mean CN (CN Range)	SD (±)	CV (%)	PCR (%) (n/N) (95% CI)
40	21	7.1 (6–8)	0.5	7.6	92 (11/12) (61.5, 99.8)
50	21	7.6 (7–9)	0.9	11.8	100 (12/12) (73.5, 100)

**c. RNA Input**

The recommended amount of input RNA for the CellDx-Tissue assay is 20 ng. A RNA input study was performed to demonstrate the analytical performance of the CellDx-Tissue assay across a range of RNA inputs (2.5, 5, 10, 20, 30, 40, and 50 ng) using three FFPE tumor samples including *ALK*, *RET*, and *ROS1* gene fusions. Each sample was tested across 12 replicates, two operators, two lot kits, and seven instruments across multiple non-consecutive days. The QC metrics and variant reporting for each sample was evaluated per RNA input level. Table 29 - Table 31 list the PCR for each input level were analyzed.

**Table 29. *ALK* Call Rates at each RNA Input Level**

RNA Input (ng)	Tumor Purity (%)	Mean Read Counts	PCR (%) (n/N) (95% CI)
2.5	29	7358	33 (4/12) (13.8-60.9%)
5	29	9711	100 (12/12) (75.7-100%)
10	29	12428	100 (12/12) (75.7-100%)
20	29	14842	100 (12/12) (75.7-100%)
30	29	16876	100 (12/12) (75.7-100%)
40	29	19954	100 (12/12) (75.7-100%)
50	29	29524	100 (12/12) (75.7-100%)

**Table 30. *RET* Call Rates at each RNA Input Level**

<b>RNA Input (ng)</b>	<b>Tumor Purity (%)</b>	<b>Mean Read Counts</b>	<b>PCR (%) (n/N) (95% CI)</b>
2.5	73	0	0 (0/12) (0-24.3%)
5	73	527	33.3 (4/12) (13.8-60.9%)
10	73	636	100 (12/12) (75.7-100%)
20	73	943	100 (12/12) (75.7-100%)
30	73	1421	100 (12/12) (75.7-100%)
40	73	2115	100 (12/12) (75.7-100%)
50	73	4771	100 (12/12) (75.7-100%)

**Table 31. *ROS1* Call Rates at each RNA Input Level**

<b>RNA Input (ng)</b>	<b>Tumor Purity (%)</b>	<b>Mean Read Counts</b>	<b>PCR (%) (n/N) (95% CI)</b>
2.5	29	0	0 (0/12) (0-24.3%)
5	29	549	75 (9/12) (46.8-91.1%)
10	29	776	100 (12/12) (75.7-100%)
20	29	1491	100 (12/12) (75.7-100%)
30	29	2665	100 (12/12) (75.7-100%)

RNA Input (ng)	Tumor Purity (%)	Mean Read Counts	PCR (%) (n/N) (95% CI)
40	29	4569	100 (12/12) (75.7-100%)
50	29	137465	100 (12/12) (75.7-100%)

## 8. Comparison Studies

### a. *Method Comparison (Accuracy)*

The analytical accuracy of the CellDx-Tissue assay as a tumor profiling assay was evaluated using 252 FFPE tumor tissue samples from patients diagnosed with various solid tumor types. The clinical FFPE samples were selected using an orthogonal method based on the clinical history (known biomarker status) and the availability of tissue for testing. Data was aggregated at the variant level for SNVs and indels and at the gene level for amplification and gene fusions. The accuracy is summarized for the entire cohort of 252 samples for each of the assessed types (SNVs, indels, *ERBB2* amplification, *ALK*, *RET*, and *ROSI* gene fusions). Orthogonal methods used consisted of validated Next Generation Sequencing (NGS), and Fluorescence In Situ Hybridization (FISH). For all analyses, the PPA and NPA were calculated by comparing the concordance between the CellDx-Tissue assay to the appropriate comparator to evaluate the degree of concordance between the assays. (Table 32).

**Table 32. Accuracy of the CellDx-Tissue Assay**

Variant Category	Orthogonal Method	Analysis Category	PPA (95% CI)	NPA (95% CI)
Overall	NGS	All	98.4 (95.5, 99.5)	99.5 (99.4, 99.6)
		SNV	98.8 (95.6, 99.7)	99.7 (99.6, 99.8)
		Insertion	85.7 (48.7, 97.4)	99.9 (99.8, 99.9)
		Deletion	100 (86.2, 100)	99.9 (99.9, 99.9)

Variant Category	Orthogonal Method	Analysis Category	PPA (95% CI)	NPA (95% CI)
Variants with Evidence of Clinical Significance	NGS	All	98.2 (90.6, 99.7)	99.8 (99.6, 99.9)
		SNV	100 (92.6, 100)	99.8 (99.7, 99.9)
		Insertion	66.7 (20.8, 93.9)	100 (99.9, 100)
		Deletion	100 (56.6, 100)	100 (99.9, 100)
Variants with Potential Clinical Significance	NGS	All	98.5 (94.8, 99.6)	99.4 (99.3, 99.5)
		SNV	98.2 (93.8, 99.5)	99.6 (99.5, 99.7)
		Insertion	100 (51.0, 100)	99.9 (99.8, 99.9)
		Deletion	100 (83.2, 100)	99.9 (99.8, 99.9)
Hotspot	NGS	All	97.8 (92.3, 99.4)	99.9 (99.8, 99.9)
		SNV	97.8 (92.3, 99.4)	99.9 (99.8, 99.9)
		Insertion	-	100 (99.8, 100)
		Deletion	-	100 (99.8, 100)
Non-Hotspot	NGS	All	99 (94.6, 99.8)	99.6 (99.5, 99.7)
		SNV	100 (94.8, 100)	99.8 (99.7, 99.9)
		Insertion	85.7 (48.7, 97.4)	99.9 (99.8, 99.9)
		Deletion	100 (86.2, 100)	99.9 (99.9, 99.9)
<i>ERBB2</i>	FISH	Amplification	91.18 (77.04, 96.95)	100 (92.13, 100)
<i>ALK</i>	FISH	Fusion	94.10 (73.0, 99.0)	100 (89.8, 100)
<i>RET</i>	FISH	Fusion	100 (43.9, 100)	100 (92.7, 100)
<i>ROS1</i>	FISH	Fusion	91.70 (64.6, 98.5)	100 (91.0, 100)

**b. Accuracy – SNVs and Indels**

The CellDx-tissue accuracy study included a total of 189 unique variants in 88 exons over 41 genes. Variants included 159 SNVs, six insertions, and 24 deletions. Performance was stratified by mutation type and gene for positive percent agreement (PPA) and negative percent agreement (NPA) with two-sided 95% confidence interval (CI). Overall, the CellDx-Tissue assay yielded concordant analytical performance for variant calls across the SNVs and Indels with a PPA of 98.4% and an NPA  $\geq$  99.5%. The concordance between detected mutations for each category of variants (Variants with Evidence of Clinical Significance, Variants with Potential Clinical Significance, Hotspot, and Non-Hotspot) are shown in Table 33. The complete listing of data can be found in Appendix E.

**Table 33. SNV and Indel Concordance between CellDx-Tissue Assay and Orthogonal NGS Comparator**

<b>Variant Category</b>	<b>Analysis Category</b>	<b>PPA<sup>1</sup> (n/N) (95% CI)</b>	<b>NPA<sup>1</sup> (n/N) (95% CI)</b>
Overall	All	98.40% (189/192) (95.5%, 99.5%)	99.50% (18784/18880) (99.4%, 99.6%)
	SNV	98.80% (159/161) (95.6%, 99.7%)	99.70% (18851/18911) (99.6%, 99.8%)
	Insertions	85.70% (6/7) (48.7%, 97.4%)	99.90% (19046/19065) (99.8%, 99.9%)
	Deletions	100% (24/24) (86.2%, 100%)	99.90% (19031/19048) (99.9%, 99.9%)
Variants with Evidence of Clinical Significance	All	98.20% (55/56) (90.6%, 99.7%)	99.80% (4405/4414) (99.6%, 99.9%)
	SNV	100% (48/48) (92.6%, 100%)	99.80% (4415/4422) (99.7%, 99.9%)
	Insertions	66.70% (2/3) (20.8%, 93.9%)	100% (4466/4467) (99.9%, 100%)
	Deletions	100% (5/5) (56.6%, 100%)	100% (4464/4465) (99.9%, 100%)

<b>Variant Category</b>	<b>Analysis Category</b>	<b>PPA<sup>1</sup> (n/N) (95% CI)</b>	<b>NPA<sup>1</sup> (n/N) (95% CI)</b>
Variants with Potential Clinical Significance	All	98.50% (134/136) (94.8%, 99.6%)	99.40% (14379/14466) (99.3%, 99.5%)
	SNV	98.20% (111/113) (93.8%, 99.5%)	99.60% (14436/14489) (99.5%, 99.7%)
	Insertions	100% (4/4) (51.0%, 100%)	99.90% (14580/14598) (99.8%, 99.9%)
	Deletions	100% (19/19) (83.2%, 100%)	99.90% (14567/14583) (99.8%, 99.9%)
Hotspot	All	97.80% (89/91) (92.3%, 99.4%)	99.90% (18959/18981) (99.8%, 99.9%)
	SNV	97.80% (89/91) (92.3%, 99.4%)	99.90% (18959/18981) (99.8%, 99.9%)
	Insertions	-	100% (19072/19072) (100%, 100%)
	Deletions	-	100% (19072/19072) (100%, 100%)
Non-Hotspot	All	99.00% (100/101) (94.6%, 99.8%)	99.60% (18897/18971) (99.5%, 99.7%)
	SNV	100% (70/70) (94.8%, 100%)	99.80% (18964/19002) (99.7%, 99.9%)
	Insertions	85.70% (6/7) (48.7%, 97.4%)	99.90% (19046/19065) (99.8%, 99.9%)
	Deletions	100% (24/24) (86.2%, 100%)	99.90% (19031/19048) (99.9%, 99.9%)

Variant Category	Analysis Category	PPA <sup>1</sup> (n/N) (95% CI)	NPA <sup>1</sup> (n/N) (95% CI)
Insertions	1 – 5 bp	80.00% (4/5) (37.6%, 96.4%)	99.90% (19049/19067) (99.9%, 99.9%)
	6 – 10 bp	100% (1/1) (20.7%, 100%)	100% (19070/19071) (100%, 100%)
	11 – 20 bp	100% (1/1) (20.7%, 100%)	100% (19071/19071) (100%, 100%)
Deletions	1 – 5 bp	100% (19/19) (83.2%, 100%)	99.90% (19038/19053) (99.9%, 100%)
	6 – 10 bp	-	100% (19071/19072) (100%, 100%)
	11 – 20 bp	100% (3/3) (43.9%, 100%)	100% (19069/19069) (100%, 100%)
	21 – 30 bp	100% (2/2) (34.2%, 100%)	100% (19069/19070) (100%, 100%)

<sup>1</sup>Two-sided CIs were estimated using the Wilson method.

**c. Accuracy – Amplification and Translocations**

**i. ERBB2 Amplification Concordance**

In total, 83 different FFPE samples representing 11 different tumor types, including Breast, Bowel, Esophagus/Stomach, Lung, Head and Neck, Ovary/Fallopian Tube, Bladder/Urinary Tract, Uterus, Biliary Tract, Bone and Prostate were analyzed for concordance between FISH status and CellDx-Tissue *ERBB2* status. Four samples that had a positive *ERBB2* FISH result were classified as Indeterminate by the CellDx-Tissue (below assay’s limit of detection for *ERBB2*) and were excluded from analysis (see Table 34). Three samples had a positive *ERBB2* FISH result but were negative for *ERBB2* in the CellDx-Tissue assay. Of these three samples, one had a FISH score of 2.02 (borderline) and history of *HER2*-negative IHC findings. The remaining two discordant samples had a FISH score of 2.84 and 2.94, with a history of *HER2*-negative IHC reported. The PPA and NPA values for *ERBB2* amplification reflected the across borderline and non-borderline samples (Table 34). In non-borderline cases (excluding all cases of a FISH ratio 1.5 – 2.5), a PPA of 93.55% (95% CI: 79.28%, 98.21%) and an NPA of 100% (95% CI: 89.85%, 100%) (Table 35) was observed.

**Table 34. *ERBB2* Amplification Concordance between CellDx-Tissue Assay and FISH**

CellDx-Tissue	FISH		
	<i>ERBB2</i> (+)	<i>ERBB2</i> (-)	Total
<i>ERBB2</i> (+) (CN ≥ 8.5)	31	0	31
<i>ERBB2</i> Indeterminate* (CN ≥ 4 to < 8.5)	4 <sup>#</sup>	0	4
<i>ERBB2</i> (-) (CN < 4)	3	45	48
<b>Total</b>	38	45	83
<b>PPA (2-sided 95% CI)</b>	91.18% (77.04, 96.95)		
<b>NPA (2-sided 95% CI)</b>	100% (92.13, 100)		

\*Samples with CN ≥ 4 to <8.5 are classified as "Indeterminate / Low-Level Amplification (below LoD)" and require reflex testing with another FDA-cleared test. Per the tiered CellDx – Tissue device *ERBB2* reporting framework, these specimens are excluded from the primary PPA calculation as they are not definitive positive or negative calls by the assay.

<sup>#</sup>Four FISH-positive samples yielded Indeterminate / Low-Level Amplification (Below LOD, CN ≥4 to <8.5) by CellDx-Tissue.

**Table 35. *ERBB2* Amplification Concordance between CellDx-Tissue Assay and FISH Borderline Samples**

Category	Total Cases	TP	FN	TN	FN	I/LL*	PPA (%) (2-sided 95% CI)	NPA (%) (2-sided 95% CI)
All	83	31	0	45	3	4	91.18 (77.04, 96.95)	100 (92.13, 100)
Excluding FISH 1.8-2.2	78	31	0	42	2	3	93.94 (80.39, 98.32)	100 (91.62, 100)
Only FISH 1.8-2.2	5	0	0	3	1	1	0 (0, 79.35)	100 (43.85, 100)
Excluding FISH 1.5-2.5	68	29	0	34	2	3	93.55 (79.28, 98.21)	100 (89.85, 100)
Only FISH 1.5-2.5	15	2	0	11	1	1	66.67 (20.77, 93.85)	100 (74.12, 100)

\*Indeterminate / Low-Level Amplification (below LoD).

**ii. *ALK* Fusion Concordance**

A total of 51 FFPE samples from six tumor types, including Lung, Thyroid, Soft Tissue, Breast, Bowel and Prostate, were included in the analysis. Out of the 17 *ALK* positive samples, 16 were detected for *ALK* gene fusions using both FISH and CellDx-Tissue assay. One *ALK* positive sample by FISH was not reported as positive

by the CellDx-Tissue assay due to borderline tumor purity. The concordance between the CellDx-Tissue assay and FISH is shown in Table 36.

**Table 36. *ALK* Fusion Concordance between CellDx-Tissue Assay and FISH**

CellDx-Tissue	FISH		
	<i>ALK</i> (+)	<i>ALK</i> (-)	Total
<i>ALK</i> (+)	16	0	16
<i>ALK</i> (-)	1	34	35
<b>Total</b>	17	34	51
<b>PPA (2-sided 95% CI)</b>	94.10% (73.0, 99.0)		
<b>NPA (2-sided 95% CI)</b>	100% (89.8, 100)		

**iii. *RET* Fusion Concordance**

A total of 52 FFPE samples from six tumor types, including Lung, Thyroid, Bowel, Soft Tissue, Breast and Prostate, were included in the analysis. Out of the three *RET* positive samples, three were detected for *RET* gene fusions using both FISH and CellDx-Tissue assay. The concordance between the CellDx-Tissue assay and FISH is shown in Table 37.

**Table 37. *RET* Fusion Concordance between CellDx-Tissue Assay and FISH**

CellDx-Tissue	FISH		
	<i>RET</i> (+)	<i>RET</i> (-)	Total
<i>RET</i> (+)	3	0	3
<i>RET</i> (-)	0	49	49
<b>Total</b>	3	49	52
<b>PPA (2-sided 95% CI)</b>	100% (43.9, 100)		
<b>NPA (2-sided 95% CI)</b>	100% (92.7, 100)		

**iv. *ROS1* Fusion Concordance**

A total of 51 FFPE samples from six tumor types, including Lung, Thyroid, Bowel, Soft Tissue, Breast and Prostate, were included in the analysis. Out of the 12 *ROS1* positive samples, 11 were detected for *ROS1* gene fusions using both FISH and CellDx-Tissue assay. One *ROS1* positive sample by FISH was not reported as positive by the CellDx-Tissue assay; this discordant result was confirmed by an independent NGS test and classified as a false negative. The concordance between the CellDx-Tissue assay and FISH is shown in Table 38.

**Table 38. ROS1 Fusion Concordance between CellDx-Tissue Assay and FISH**

CellDx-Tissue	FISH		
	ROS1 (+)	ROS1 (-)	Total
ROS1 (+)	11	0	11
ROS1 (-)	1	39	40
Total	12	39	51
PPA (2-sided 95% CI)	91.70% (64.6, 98.5)		
NPA (2-sided 95% CI)	100% (91.0, 100)		

**d. Method Comparison – Wild-Type Calls**

The accuracy of wild-type calls was evaluated by observing the NPA of relevant positions of Variants with Evidence of Clinical Significance (Table 39).

**Table 39. CellDx-Tissue Wild-Type Accuracy for Variants with Evidence of Clinical Significance.**

Variant Category	NPA (%) (n/N) (95% CI Wilson Score)
SNVs	99.80 (4415/4422) (99.7, 99.9)
Insertions	99.98 (4466/4467) (99.9, 100)
Deletions	99.98 (4464/4465) (99.9, 100)

**9. Tissue Comparability**

A tissue comparability study was performed to evaluate the CellDx-Tissue assay’s analytical performance across a diverse range of solid tumor tissue types. The study assessed the CellDx-Tissue assay’s performance across more than 20 tumor types for all reportable biomarkers (SNVs, Indels, *ERBB2* gene amplification, and *ALK*, *RET*, and *ROS1* gene fusions) from FFPE samples. The analytical performance of the CellDx-Tissue assay is consistent across the full spectrum of tumor types within the intended use population. The number of false negative (FN) variant calls, false positive (FP) variant calls, true positive (TP) variant calls, and the positive percent agreement (PPA) and false negative rate (FNR) for each tumor type are shown for all reportable biomarkers in Table 40 – Table 44.

**Table 40. CellDx-Tissue Assay Performance in SNVs, Insertions, and Deletions**

Tissue	Sample Number	Variant Calls				PPA (%) (95% CI)	NPA (%) (95% CI)
		FN	FP	TN	TP		
Adrenal Gland	1			1410148	1	100 (20.65, 100)	100 (99.99, 100)
Biliary Tract	2		3	2820294	1	100 (20.65, 100)	99.99 (99.99, 100)
Bladder/Urinary Tract	2			2820295	3	100 (43.85,100)	100 (99.99, 100)
Bowel	25	1	29	35253643	52	98.11 (90.06, 99.67)	99.99 (99.99, 100)
Breast	26		10	36663825	39	100 (91.03, 100)	100 (99.99, 100)
Cervix	3		3	4230442	2	100 (34.24, 100)	100 (99.99, 100)
CNS/Brain	5			7050735	10	100 (72.25, 100)	100 (99.99, 100)
Esophagus/Stomach	3		2	4230442	3	100 (43.85, 100)	100 (99.99, 100)
Head and neck	6		3	8460884	7	100 (64.57, 100)	100 (99.99, 100)
Lung	15		9	21152197	29	100 (88.30, 100)	100 (99.99, 100)
Ovary/Fallopian Tube	8	1	9	11281177	5	83.33 (43.65, 96.99)	100 (99.99, 100)
Pancreas	5		9	7050731	5	100 (56.55, 100)	100 (99.99, 100)
Peritoneum	1		1	1410148		NA	100 (99.99, 100)
Prostate	3		3	4230440	4	100 (51.01, 100)	100 (99.99, 100)
Skin	3	1	1	4230443	2	66.67 (20.77, 93.85)	100 (99.99, 100)
Soft Tissue	3			4230440	7	100 (64.57, 100)	100 (99.99, 100)
Testis	1			1410147	2	100 (34.24, 100)	100 (99.99, 100)
Thymus	2		4	2820292	2	100 (34.24, 100)	100 (99.99, 100)
Thyroid	3			4230444	3	100 (43.85, 100)	100 (99.99, 100)
Unknown Primary	1			1410148	1	100 (20.65, 100)	100 (99.99, 100)
Uterus	5		10	7050724	11	100 (74.12, 100)	100 (99.99, 100)
<b>Overall</b>	<b>123</b>	<b>3</b>	<b>96</b>	<b>173448039</b>	<b>189</b>	<b>98.44</b> <b>(95.51, 99.47)</b>	<b>100</b> <b>(99.99, 100)</b>

**Table 41. CellDx-Tissue Assay Performance in *ERBB2* Amplification**

Tissue (N)	Variant Calls					PPA (%) (95% CI)	NPA (%) (95% CI)
	FN	FP	TN	TP	I/LL*		
Breast (4)	1		23	18		94.74 (75.36, 99.06)	100 (85.69, 100)
Bowl (32)	1		21	8	2	88.89 (56.50, 98.01)	100 (84.54, 100)
Lung (11)			4	6	1	100 (60.97, 100)	100 (51.01, 100)
Head and Neck (7)	1		4	1	1	100 (34.24, 100)	100 (51.01, 100)
Esophagus/ Stomach (6)			2	4		100 (34.24, 100)	100 (51.01, 100)
Ovary/ Fallopian Tube (6)			4	2		100 (51.01, 100)	100 (34.24, 100)
Pancreas (4)			3	1		100 (34.24, 100)	100 (34.24, 100)
Uterus (4)			2	2		100 (20.65, 100)	100 (43.85, 100)
Bladder/ Urinary Track (3)			2	1		NA	100 (43.85, 100)
Prostate (3)			3			100 (20.65, 100)	100 (34.24, 100)
Biliary Tract (2)				2		NA	100 (34.24, 100)
Cervix (2)			2			100 (34.24, 100)	NA
Bone (1)			1			NA	100 (20.65, 100)
CNS/ Brain (1)			1			NA	100 (20.65, 100)
Liver (1)			1			NA	100 (20.65, 100)

Tissue (N)	Variant Calls					PPA (%) (95% CI)	NPA (%) (95% CI)
	FN	FP	TN	TP	I/LL*		
Thymus (1)			1			NA	100 (20.65, 100)
<b>Overall (126)</b>	<b>3</b>	<b>0</b>	<b>74</b>	<b>45</b>	<b>4</b>	<b>93.75 (83.16, 97.85)</b>	<b>100 (95.06, 100)</b>

\*Indeterminate / Low-Level Amplification (below LoD)

**Table 42. CellDx-Tissue Assay Performance in *ALK* Fusion**

Tissue (N)	Variant Calls				PPA (%) (95% CI)	NPA (%) (95% CI)
	FN	FP	TN	TP		
Lung (64)			36	28	100 (84.94, 100)	100 (90.36, 100)
Thyroid (13)			13		NA	100 (77.19, 100)
Head and Neck (6)			6		NA	100 (60.97, 100)
Breast (5)			4	1	100 (20.65, 100)	100 (51.01, 100)
Bowel (4)			4		NA	100 (51.01, 100)
Cervix (3)			3		NA	100 (43.85, 100)
Ovary/ Fallopian Tube (3)			3		NA	100 (43.85, 100)
Pancreas (3)			3		NA	100 (43.85, 100)
Prostate (3)			3		NA	100 (43.85, 100)
Soft Tissue (3)	1		2		0 (0, 79.35)	100 (34.24, 100)
Uterus (3)			3		NA	100 (43.85, 100)
Esophagus/ Stomach (2)			2		NA	100 (34.24, 100)
Bladder/ Urinary Tract (1)			1		NA	100 (20.65, 100)

Tissue (N)	Variant Calls				PPA (%) (95% CI)	NPA (%) (95% CI)
	FN	FP	TN	TP		
Liver (1)			1		NA	100 (20.65, 100)
Thymus (1)			1		NA	100 (20.65, 100)
<b>Overall (115)</b>	<b>1</b>		<b>85</b>	<b>29</b>	<b>96.67 (83.33, 99.41)</b>	<b>100 (95.68, 100)</b>

**Table 43. CellDx-Tissue Assay Performance in *RET* Fusion**

Tissue (N)	Variant Calls				PPA (%) (95% CI)	NPA (%) (95% CI)
	FN	FP	TN	TP		
Lung (51)			45	6	100 (54.07, 100)	100 (92.13, 100)
Thyroid (15)			14	1	100 (20.65, 100)	100 (78.47, 100)
Head and Neck (6)			6		NA	100 (60.97, 100)
Breast (5)			5		NA	100 (56.55, 100)
Bowel (4)			4		NA	100 (51.01, 100)
Cervix (3)			3		NA	100 (43.85, 100)
Ovary/ Fallopian Tube (3)			3		NA	100 (43.85, 100)
Pancreas (3)			3		NA	100 (43.85, 100)
Prostate (3)			3		NA	100 (43.85, 100)
Soft Tissue (3)			3		NA	100 (43.85, 100)
Uterus (3)			3		NA	100 (43.85, 100)
Esophagus/ Stomach (2)			2		NA	100 (32.24, 100)
Bladder/ Urinary Tract (1)			1		NA	100 (20.65, 100)

Tissue (N)	Variant Calls				PPA (%) (95% CI)	NPA (%) (95% CI)
	FN	FP	TN	TP		
Liver (1)			1		NA	100 (20.65, 100)
Thymus (1)			1		NA	100 (20.65, 100)
<b>Overall (104)</b>			<b>97</b>	<b>7</b>	<b>100 (59.04, 100)</b>	<b>100 (96.26, 100)</b>

**Table 44. CellDx-Tissue Assay Performance in *ROS1* Fusion**

Tissue (N)	Variant Calls				PPA (%) (95% CI)	NPA (%) (95% CI)
	FN	FP	TN	TP		
Lung (59)	1		38	20	95.24 (77.33, 99.15)	100 (90.77, 100)
Thyroid (13)			13		NA	100 (77.19, 100)
Head and Neck (6)			6		NA	100 (60.97, 100)
Breast (5)			5		NA	100 (56.55, 100)
Bowel (4)			4		NA	100 (51.01, 100)
Cervix (3)			3		NA	100 (43.85, 100)
Ovary/ Fallopian Tube (3)			3		NA	100 (43.85, 100)
Pancreas (3)			3		NA	100 (43.85, 100)
Prostate			2	1	NA	66.67 (20.77, 93.85)
Soft Tissue (3)			3		NA	100 (43.85, 100)
Uterus (3)			3		NA	100 (43.85, 100)
Esophagus/ Stomach (2)			2		NA	100 (34.24, 100)

Tissue (N)	Variant Calls				PPA (%) (95% CI)	NPA (%) (95% CI)
	FN	FP	TN	TP		
Bladder/Urinary Tract (1)			1		NA	100 (20.65, 100)
Liver (1)			1		NA	100 (20.65, 100)
Thymus (1)			1		NA	100 (20.65, 100)
<b>Overall (110)</b>	<b>1</b>		<b>88</b>	<b>21</b>	<b>95.45 (72.20, 99.19)</b>	<b>100 (95.82, 100)</b>

## V. Instrument Name

Thermo Fisher Ion GeneStudio S5 Prime NGS Systems (qualified by Datar Cancer Genetics)

## VI. System Description

### 1. Modes of Operations:

Does the applicant's device contain the ability to transmit data to a computer, webserver, or mobile device?

Yes:  Or No:

Does the applicant's device transmit data to a computer, webserver, or mobile device using wireless transmission?

Yes:  Or No:

### 2. Software:

FDA has reviewed applicant's Hazard Analysis and software development processes for this line of product types:

Yes:  Or No:

### 3. Level of Concern:

Moderate

### 4. Specimen Handling:

Refer to Device Description section above

### 5. Calibration and Quality Controls

Refer to Device Description section above

## VII. Other Supportive Instrument Performance Characteristics

Not applicable

### **VIII. Proposed Labeling**

The labeling is sufficient, and it satisfies the requirements of 21 CFR Parts 801 and 809, as applicable.

### **IX. Conclusion**

The submitted information in this premarket notification is complete and supports a substantial equivalent conclusion.

## X. Appendix

### Appendix A. List of Genes on the CellDx-Tissue Assay Panel

Gene	Transcript ID	Full-Length Coverage
<i>AICF</i>	NM_138932.2	-
<i>ABCB1</i>	NM_000927.4	-
<i>ABL1</i>	NM_005157.6	-
<i>ABL2</i>	NM_005158.5	-
<i>ABRAXAS1</i>	NM_139076.3	(Full-length)
<i>ACSM2B</i>	NM_001105069.2	-
<i>ACVR1</i>	NM_001111067.4	-
<i>ACVR1B</i>	NM_020328.4	(Full-length)
<i>ACVR2A</i>	NM_001616.5	(Full-length)
<i>ADAM18</i>	NM_014237.3	-
<i>ADAMTS12</i>	NM_030955.4	(Full-length)
<i>ADAMTS2</i>	NM_014244.5	(Full-length)
<i>AKT1</i>	NM_001014431.2	-
<i>AKT2</i>	NM_001626.6	-
<i>AKT3</i>	NM_005465.7	-
<i>ALK</i>	NM_004304.5	-
<i>AMER1</i>	NM_152424.4	(Full-length)
<i>ANO4</i>	NM_178826.4	-
<i>APC</i>	NM_000038.6	(Full-length)
<i>AR</i>	NM_000044.6	-
<i>ARAF</i>	NM_001654.5	-
<i>ARHGAP35</i>	NM_004491.5	(Full-length)
<i>ARID1A</i>	NM_006015.6	(Full-length)
<i>ARID1B</i>	NM_001371656.1	(Full-length)
<i>ARID2</i>	NM_152641.4	(Full-length)
<i>ARID5B</i>	NM_032199.3	(Full-length)
<i>ARMC4</i>	NM_018076.5	-
<i>ASXL1</i>	NM_015338.6	(Full-length)
<i>ASXL2</i>	NM_018263.6	(Full-length)
<i>ATM</i>	NM_000051.3	(Full-length)
<i>ATP1A1</i>	NM_000701.8	-
<i>ATR</i>	NM_001184.4	(Full-length)
<i>ATRX</i>	NM_000489.5	(Full-length)
<i>AURKA</i>	NM_003600.4	-
<i>AURKB</i>	NM_004217.4	-
<i>AURKC</i>	NM_001015878.2	-
<i>AXIN1</i>	NM_003502.4	(Full-length)
<i>AXIN2</i>	NM_004655.4	(Full-length)
<i>AXL</i>	NM_021913.5	-
<i>B2M</i>	NM_004048.3	(Full-length)

<b>Gene</b>	<b>Transcript ID</b>	<b>Full-Length Coverage</b>
<i>BAP1</i>	NM_004656.4	(Full-length)
<i>BARD1</i>	NM_000465.4	(Full-length)
<i>BCL2</i>	NM_000633.2	-
<i>BCL2L12</i>	NM_138639.1	-
<i>BCL6</i>	NM_001706.5	-
<i>BCOR</i>	NM_001123385.2	(Full-length)
<i>BCR</i>	NM_004327.4	-
<i>BLM</i>	NM_000057.4	(Full-length)
<i>BMP5</i>	NM_021073.4	-
<i>BMPR2</i>	NM_001204.7	(Full-length)
<i>BRAF</i>	NM_004333.6	-
<i>BRCA1</i>	NM_007294.4	(Full-length)
<i>BRCA2</i>	NM_000059.3	(Full-length)
<i>BRINP3</i>	NM_199051.3	-
<i>BRIP1</i>	NM_032043.3	(Full-length)
<i>BTK</i>	NM_000061.3	-
<i>C6</i>	NM_000065.4	-
<i>C8A</i>	NM_000562.3	-
<i>C8B</i>	NM_000066.4	-
<i>CACNA1D</i>	NM_000720.4	-
<i>CALR</i>	NM_004343.4	(Full-length)
<i>CANX</i>	NM_001024649.2	-
<i>CARD11</i>	NM_032415.6	-
<i>CASP8</i>	NM_001080125.2	(Full-length)
<i>CASR</i>	NM_001178065.2	-
<i>CBFB</i>	NM_022845.3	(Full-length)
<i>CBL</i>	NM_005188.4	-
<i>CCND1</i>	NM_053056.3	-
<i>CCND2</i>	NM_001759.4	-
<i>CCND3</i>	NM_001760.5	-
<i>CCNE1</i>	NM_001238.4	-
<i>CD163</i>	NM_004244.5	-
<i>CD274</i>	NM_014143.4	(Full-length)
<i>CD276</i>	NM_001024736.2	(Full-length)
<i>CD79B</i>	NM_001039933.3	-
<i>CDC73</i>	NM_024529.5	(Full-length)
<i>CDH1</i>	NM_004360.5	(Full-length)
<i>CDH10</i>	NM_006727.5	(Full-length)
<i>CDK12</i>	NM_016507.4	(Full-length)
<i>CDK4</i>	NM_000075.4	-
<i>CDK6</i>	NM_001145306.2	-
<i>CDKN1A</i>	NM_078467.3	(Full-length)
<i>CDKN1B</i>	NM_004064.4	(Full-length)

<b>Gene</b>	<b>Transcript ID</b>	<b>Full-Length Coverage</b>
<i>CDKN2A</i>	NM_001195132.1	(Full-length)
<i>CDKN2B</i>	NM_004936.4	(Full-length)
<i>CDKN2C</i>	NM_078626.3	(Full-length)
<i>CHD4</i>	NM_001273.5	-
<i>CHEK1</i>	NM_001274.5	(Full-length)
<i>CHEK2</i>	NM_007194.4	(Full-length)
<i>CIC</i>	NM_015125.4	(Full-length)
<i>CIITA</i>	NM_000246.3	(Full-length)
<i>CNTN6</i>	NM_014461.4	-
<i>CNTNAP4</i>	NM_138994.5	-
<i>CNTNAP5</i>	NM_130773.4	-
<i>COL11A1</i>	NM_001854.4	-
<i>CREBBP</i>	NM_004380.3	(Full-length)
<i>CSF1R</i>	NM_005211.3	-
<i>CSMD3</i>	NM_198123.2	(Full-length)
<i>CTCF</i>	NM_006565.4	(Full-length)
<i>CTLA4</i>	NM_005214.5	(Full-length)
<i>CTNNB1</i>	NM_001904.4	-
<i>CTNND2</i>	NM_001332.4	-
<i>CUL1</i>	NM_003592.3	-
<i>CUL3</i>	NM_003590.5	(Full-length)
<i>CUL4A</i>	NM_001008895.4	(Full-length)
<i>CUL4B</i>	NM_003588.3	(Full-length)
<i>CYLD</i>	NM_001042355.2	(Full-length)
<i>CYP2C9</i>	NM_000771.4	(Full-length)
<i>CYP2D6</i>	NM_000106.6	(Full-length)
<i>CYSLTR2</i>	NM_020377.5	-
<i>DAXX</i>	NM_001141970.2	(Full-length)
<i>DCAF4L2</i>	NM_152418.4	-
<i>DCDC1</i>	NM_001367979.1	-
<i>DDR1</i>	NM_001954.4	-
<i>DDR2</i>	NM_006182.4	-
<i>DDX3X</i>	NM_001356.5	(Full-length)
<i>DGCR8</i>	NM_022720.7	-
<i>DICER1</i>	NM_030621.4	(Full-length)
<i>DNMT3A</i>	NM_022552.4	(Full-length)
<i>DOCK3</i>	NM_004947.5	(Full-length)
<i>DPYD</i>	NM_000110.4	(Full-length)
<i>DROSHA</i>	NM_013235.5	-
<i>DSC1</i>	NM_024421.2	(Full-length)
<i>DSC3</i>	NM_001941.5	(Full-length)
<i>E2F1</i>	NM_005225.3	-
<i>EGFR</i>	NM_005228.5	-

<b>Gene</b>	<b>Transcript ID</b>	<b>Full-Length Coverage</b>
<i>EIF1AX</i>	NM_001412.4	-
<i>ELF3</i>	NM_004433.5	(Full-length)
<i>EMSY</i>	NM_020193.4	-
<i>ENO1</i>	NM_001428.5	(Full-length)
<i>EP300</i>	NM_001429.4	(Full-length)
<i>EPAS1</i>	NM_001430.5	-
<i>EPCAM</i>	NM_002354.3	(Full-length)
<i>EPHA2</i>	NM_004431.5	(Full-length)
<i>ERAP1</i>	NM_016442.4	(Full-length)
<i>ERAP2</i>	NM_001130140.2	(Full-length)
<i>ERBB2</i>	NM_004448.3	-
<i>ERBB3</i>	NM_001982.4	-
<i>ERBB4</i>	NM_005235.3	-
<i>ERCC2</i>	NM_000400.4	(Full-length)
<i>ERCC4</i>	NM_005236.3	(Full-length)
<i>ERCC5</i>	NM_000123.4	(Full-length)
<i>ERG</i>	NM_182918.4	-
<i>ERRF11</i>	NM_018948.4	(Full-length)
<i>ESR1</i>	NM_001122740.1	-
<i>ETV1</i>	NM_001163147.1	-
<i>ETV4</i>	NM_001986.4	-
<i>ETV5</i>	NM_004454.3	-
<i>ETV6</i>	NM_001987.5	(Full-length)
<i>EZH2</i>	NM_004456.5	-
<i>FAM135B</i>	NM_015912.4	-
<i>FANCA</i>	NM_000135.4	(Full-length)
<i>FANCC</i>	NM_000136.3	(Full-length)
<i>FANCD2</i>	NM_033084.6	(Full-length)
<i>FANCE</i>	NM_021922.3	(Full-length)
<i>FANCF</i>	NM_022725.4	(Full-length)
<i>FANCG</i>	NM_004629.2	(Full-length)
<i>FANCI</i>	NM_001113378.2	(Full-length)
<i>FANCL</i>	NM_001114636.1	(Full-length)
<i>FANCM</i>	NM_020937.4	(Full-length)
<i>FAS</i>	NM_000043.6	(Full-length)
<i>FAT1</i>	NM_005245.4	(Full-length)
<i>FBXW7</i>	NM_033632.3	(Full-length)
<i>FGF19</i>	NM_005117.3	-
<i>FGF23</i>	NM_020638.3	-
<i>FGF3</i>	NM_005247.4	-
<i>FGF4</i>	NM_002007.4	-
<i>FGF7</i>	NM_002009.4	-
<i>FGF9</i>	NM_002010.3	-

<b>Gene</b>	<b>Transcript ID</b>	<b>Full-Length Coverage</b>
<i>FGFR1</i>	NM_001174067.1	-
<i>FGFR2</i>	NM_000141.5	-
<i>FGFR3</i>	NM_000142.4	-
<i>FGFR4</i>	NM_213647.3	-
<i>FLT3</i>	NM_004119.3	-
<i>FLT4</i>	NM_182925.5	-
<i>FOXA1</i>	NM_004496.5	-
<i>FOXL2</i>	NM_023067.4	-
<i>FOXO1</i>	NM_002015.4	-
<i>FUBP1</i>	NM_003902.4	(Full-length)
<i>FYN</i>	NM_153047.4	-
<i>GALNT17</i>	NM_022479.3	-
<i>GATA2</i>	NM_032638.5	-
<i>GATA3</i>	NM_001002295.2	(Full-length)
<i>GLI1</i>	NM_005269.3	-
<i>GLI3</i>	NM_000168.6	-
<i>GNAI1</i>	NM_002067.5	-
<i>GNAI3</i>	NM_006572.6	(Full-length)
<i>GNAQ</i>	NM_002072.5	-
<i>GNAS</i>	NM_000516.6	-
<i>GPR158</i>	NM_020752.3	-
<i>GPS2</i>	NM_004489.5	(Full-length)
<i>GRID2</i>	NM_001510.4	-
<i>H1-4</i>	NM_005321.3	-
<i>H2BC5</i>	NM_021063.4	-
<i>H3-3A</i>	NM_002107.7	-
<i>H3-3B</i>	NM_005324.5	-
<i>H3C2</i>	NM_003537.4	-
<i>HCN1</i>	NM_021072.4	-
<i>HDAC2</i>	NM_001527.4	(Full-length)
<i>HDAC9</i>	NM_178425.3	(Full-length)
<i>HIF1A</i>	NM_001530.4	-
<i>HLA-A</i>	NM_001242758.1	(Full-length)
<i>HLA-B</i>	NM_005514.8	(Full-length)
<i>HLA-C</i>	NM_001243042.1	-
<i>HNF1A</i>	NM_000545.8	(Full-length)
<i>HRAS</i>	NM_001130442.2	-
<i>ID3</i>	NM_002167.5	(Full-length)
<i>IDH1</i>	NM_005896.3	-
<i>IDH2</i>	NM_002168.4	-
<i>IGF1R</i>	NM_000875.5	-
<i>IKBKB</i>	NM_001556.3	-
<i>IL6ST</i>	NM_002184.4	-

<b>Gene</b>	<b>Transcript ID</b>	<b>Full-Length Coverage</b>
<i>IL7R</i>	NM_002185.5	-
<i>INPP4B</i>	NM_001101669.3	(Full-length)
<i>IRF4</i>	NM_002460.4	-
<i>IRS4</i>	NM_003604.2	-
<i>JAK1</i>	NM_002227.4	(Full-length)
<i>JAK2</i>	NM_004972.4	(Full-length)
<i>JAK3</i>	NM_000215.4	(Full-length)
<i>KCND2</i>	NM_012281.3	-
<i>KCNH7</i>	NM_033272.4	-
<i>KCNJ5</i>	NM_000890.5	-
<i>KDM5C</i>	NM_004187.5	(Full-length)
<i>KDM6A</i>	NM_021140.3	(Full-length)
<i>KDR</i>	NM_002253.3	-
<i>KEAP1</i>	NM_203500.2	(Full-length)
<i>KEL</i>	NM_000420.3	-
<i>KIR3DL1</i>	NM_001322168.1	-
<i>KIT</i>	NM_000222.3	-
<i>KLF4</i>	NM_004235.6	-
<i>KLF5</i>	NM_001730.5	-
<i>KLHL13</i>	NM_001168299.1	(Full-length)
<i>KMT2A</i>	NM_001197104.2	(Full-length)
<i>KMT2B</i>	NM_014727.3	(Full-length)
<i>KMT2C</i>	NM_170606.3	(Full-length)
<i>KMT2D</i>	NM_003482.4	(Full-length)
<i>KNSTRN</i>	NM_033286.3	-
<i>KRAS</i>	NM_033360.4	-
<i>KRTAP21-1</i>	NM_181619.2	-
<i>KRTAP6-2</i>	NM_181604.1	-
<i>LARP4B</i>	NM_015155.3	(Full-length)
<i>LATS1</i>	NM_004690.4	(Full-length)
<i>LATS2</i>	NM_014572.3	(Full-length)
<i>LRRC7</i>	NM_001370785.2	-
<i>MAGOH</i>	NM_002370.4	-
<i>MAP2K1</i>	NM_002755.4	-
<i>MAP2K2</i>	NM_030662.4	-
<i>MAP2K4</i>	NM_003010.4	(Full-length)
<i>MAP2K7</i>	NM_145185.4	(Full-length)
<i>MAP3K1</i>	NM_005921.2	(Full-length)
<i>MAP3K4</i>	NM_005922.4	(Full-length)
<i>MAP3K8</i>	NM_001244134.1	-
<i>MAPK1</i>	NM_002745.5	-
<i>MAPK8</i>	NM_139049.4	(Full-length)
<i>MARCO</i>	NM_006770.3	-

<b>Gene</b>	<b>Transcript ID</b>	<b>Full-Length Coverage</b>
<i>MAX</i>	NM_002382.5	-
<i>MCL1</i>	NM_021960.5	-
<i>MDM2</i>	NM_002392.5	-
<i>MDM4</i>	NM_002393.5	-
<i>MECOM</i>	NM_004991.4	-
<i>MED12</i>	NM_005120.3	-
<i>MEF2B</i>	NM_001145785.2	-
<i>MEN1</i>	NM_000244.3	(Full-length)
<i>MET</i>	NM_001127500.3	-
<i>MGA</i>	NM_001164273.1	(Full-length)
<i>MITF</i>	NM_198159.3	-
<i>MLH1</i>	NM_000249.4	(Full-length)
<i>MLH3</i>	NM_001040108.2	(Full-length)
<i>MPL</i>	NM_005373.3	-
<i>MRE11</i>	NM_005591.4	(Full-length)
<i>MSH2</i>	NM_000251.3	(Full-length)
<i>MSH3</i>	NM_002439.5	(Full-length)
<i>MSH6</i>	NM_000179.3	(Full-length)
<i>MTAP</i>	NM_002451.4	(Full-length)
<i>MTOR</i>	NM_004958.4	-
<i>MTUS2</i>	NM_001033602.3	(Full-length)
<i>MUTYH</i>	NM_001128425.2	(Full-length)
<i>MYB</i>	NM_001130173.2	-
<i>MYBL1</i>	NM_001080416.4	-
<i>MYC</i>	NM_002467.6	-
<i>MYCL</i>	NM_001033082.3	-
<i>MYCN</i>	NM_005378.6	-
<i>MYD88</i>	NM_001172567.2	-
<i>MYOD1</i>	NM_002478.5	-
<i>NBN</i>	NM_002485.5	(Full-length)
<i>NCOR1</i>	NM_006311.4	(Full-length)
<i>NF1</i>	NM_001042492.3	(Full-length)
<i>NF2</i>	NM_000268.4	(Full-length)
<i>NFE2L2</i>	NM_006164.5	-
<i>NLRC5</i>	NM_032206.4	-
<i>NOL4</i>	NM_003787.4	-
<i>NOTCH1</i>	NM_017617.5	(Full-length)
<i>NOTCH2</i>	NM_024408.4	(Full-length)
<i>NOTCH3</i>	NM_000435.3	(Full-length)
<i>NOTCH4</i>	NM_004557.4	(Full-length)
<i>NRAS</i>	NM_002524.5	-
<i>NRG1</i>	NM_013956.5	-
<i>NRXN1</i>	NM_004801.5	-

<b>Gene</b>	<b>Transcript ID</b>	<b>Full-Length Coverage</b>
<i>NSD2</i>	NM_001042424.3	-
<i>NT5C2</i>	NM_001134373.3	-
<i>NTRK1</i>	NM_002529.3	-
<i>NTRK2</i>	NM_006180.5	-
<i>NTRK3</i>	NM_001012338.2	-
<i>NUP93</i>	NM_014669.5	-
<i>NUTM1</i>	NM_175741.2	-
<i>NYAP2</i>	NM_020864.1	-
<i>OR10G8</i>	NM_001004464.1	-
<i>OR2G6</i>	NM_001013355.1	-
<i>OR2L13</i>	NM_175911.3	-
<i>OR2L2</i>	NM_001004686.2	-
<i>OR2L8</i>	NM_001001963.1	-
<i>OR2M3</i>	NM_001004689.1	-
<i>OR2T3</i>	NM_001005495.1	-
<i>OR2T33</i>	NM_001004695.1	-
<i>OR2T4</i>	NM_001004696.1	-
<i>OR2W3</i>	NM_001001957.2	-
<i>OR4A15</i>	NM_001005275.1	-
<i>OR4C15</i>	NM_001001920.2	-
<i>OR4C6</i>	NM_001004704.1	-
<i>OR4M1</i>	NM_001005500.1	-
<i>OR4M2</i>	NM_001004719.2	-
<i>OR5D18</i>	NM_001001952.1	-
<i>OR5F1</i>	NM_003697.1	-
<i>OR5L1</i>	NM_001004738.2	-
<i>OR5L2</i>	NM_001004739.1	-
<i>OR6F1</i>	NM_001005286.1	-
<i>OR8H2</i>	NM_001005200.1	-
<i>OR8I2</i>	NM_001003750.1	-
<i>OR8UI</i>	NM_001005204.1	-
<i>ORC4</i>	NM_181742.3	-
<i>PAK5</i>	NM_177990.3	-
<i>PALB2</i>	NM_024675.4	(Full-length)
<i>PARP1</i>	NM_001618.4	(Full-length)
<i>PARP2</i>	NM_005484.3	(Full-length)
<i>PARP3</i>	NM_001003931.3	(Full-length)
<i>PARP4</i>	NM_006437.4	(Full-length)
<i>PAX5</i>	NM_016734.3	-
<i>PBRM1</i>	NM_018313.5	(Full-length)
<i>PCBP1</i>	NM_006196.4	-
<i>PCDH17</i>	NM_001040429.3	-
<i>PDCD1</i>	NM_005018.3	(Full-length)

<b>Gene</b>	<b>Transcript ID</b>	<b>Full-Length Coverage</b>
<i>PDCD1LG2</i>	NM_025239.4	(Full-length)
<i>PDE1A</i>	NM_001258312.1	-
<i>PDE1C</i>	NM_001191058.4	-
<i>PDGFRA</i>	NM_006206.6	-
<i>PDGFRB</i>	NM_002609.4	-
<i>PDIA3</i>	NM_005313.5	(Full-length)
<i>PGD</i>	NM_002631.4	(Full-length)
<i>PHF6</i>	NM_032458.3	(Full-length)
<i>PIK3C2B</i>	NM_002646.4	-
<i>PIK3CA</i>	NM_006218.4	-
<i>PIK3CB</i>	NM_006219.3	-
<i>PIK3CD</i>	NM_005026.5	-
<i>PIK3CG</i>	NM_002649.3	-
<i>PIK3R1</i>	NM_181523.3	(Full-length)
<i>PIK3R2</i>	NM_005027.4	-
<i>PIMI</i>	NM_002648.4	-
<i>PLCG1</i>	NM_002660.3	-
<i>PLXDC2</i>	NM_032812.9	-
<i>PMS1</i>	NM_000534.5	(Full-length)
<i>PMS2</i>	NM_000535.7	(Full-length)
<i>POLD1</i>	NM_001256849.1	(Full-length)
<i>POLE</i>	NM_006231.4	(Full-length)
<i>POM121L12</i>	NM_182595.4	-
<i>POT1</i>	NM_015450.3	(Full-length)
<i>PPARG</i>	NM_015869.5	-
<i>PPF1A2</i>	NM_003625.5	-
<i>PPM1D</i>	NM_003620.4	(Full-length)
<i>PPP2R1A</i>	NM_014225.6	-
<i>PPP2R2A</i>	NM_002717.4	(Full-length)
<i>PPP6C</i>	NM_002721.5	-
<i>PRDMI</i>	NM_001198.4	(Full-length)
<i>PRDM9</i>	NM_020227.4	(Full-length)
<i>PRKACA</i>	NM_002730.4	-
<i>PRKACB</i>	NM_182948.4	-
<i>PRKARIA</i>	NM_212471.2	(Full-length)
<i>PSMB10</i>	NM_002801.4	(Full-length)
<i>PSMB8</i>	NM_148919.4	(Full-length)
<i>PSMB9</i>	NM_002800.5	(Full-length)
<i>PTCH1</i>	NM_000264.5	(Full-length)
<i>PTEN</i>	NM_000314.8	(Full-length)
<i>PTPN11</i>	NM_002834.5	-
<i>PTPRD</i>	NM_002839.4	-
<i>PTPRT</i>	NM_133170.4	(Full-length)

<b>Gene</b>	<b>Transcript ID</b>	<b>Full-Length Coverage</b>
<i>PXDNL</i>	NM_144651.5	-
<i>RAC1</i>	NM_018890.4	-
<i>RAD50</i>	NM_005732.4	(Full-length)
<i>RAD51</i>	NM_133487.4	(Full-length)
<i>RAD51B</i>	NM_133509.4	(Full-length)
<i>RAD51C</i>	NM_058216.3	(Full-length)
<i>RAD51D</i>	NM_133629.3	(Full-length)
<i>RAD52</i>	NM_134424.4	(Full-length)
<i>RAD54L</i>	NM_001142548.1	(Full-length)
<i>RAF1</i>	NM_002880.3	-
<i>RARA</i>	NM_000964.4	-
<i>RASAI</i>	NM_002890.3	(Full-length)
<i>RASA2</i>	NM_006506.5	(Full-length)
<i>RBI</i>	NM_000321.2	(Full-length)
<i>RBM10</i>	NM_001204468.1	(Full-length)
<i>RBP3</i>	NM_002900.3	-
<i>RECQL4</i>	NM_004260.4	(Full-length)
<i>REG1A</i>	NM_002909.5	-
<i>REG1B</i>	NM_006507.4	-
<i>REG3A</i>	NM_138937.3	-
<i>REG3G</i>	NM_001008387.3	-
<i>RELA</i>	NM_021975.4	-
<i>RET</i>	NM_020975.6	-
<i>RGS7</i>	NM_002924.6	-
<i>RHEB</i>	NM_005614.4	-
<i>RHOA</i>	NM_001664.4	-
<i>RICTOR</i>	NM_152756.5	-
<i>RIT1</i>	NM_001256821.2	-
<i>RNASEH2A</i>	NM_006397.3	(Full-length)
<i>RNASEH2B</i>	NM_024570.4	(Full-length)
<i>RNASEH2C</i>	NM_032193.4	(Full-length)
<i>RNF43</i>	NM_017763.5	(Full-length)
<i>ROSI</i>	NM_002944.2	-
<i>RPA1</i>	NM_002945.5	(Full-length)
<i>RPL10</i>	NM_006013.5	-
<i>RPL22</i>	NM_000983.4	(Full-length)
<i>RPL5</i>	NM_000969.5	(Full-length)
<i>RPS6KB1</i>	NM_003161.4	-
<i>RPTN</i>	NM_001122965.1	-
<i>RPTOR</i>	NM_020761.3	-
<i>RSPO2</i>	NM_178565.5	-
<i>RSPO3</i>	NM_032784.5	-
<i>RUNDC3B</i>	NM_138290.3	-

<b>Gene</b>	<b>Transcript ID</b>	<b>Full-Length Coverage</b>
<i>RUNX1</i>	NM_001754.4	(Full-length)
<i>RUNX1T1</i>	NM_001198634.2	(Full-length)
<i>SDHA</i>	NM_004168.4	(Full-length)
<i>SDHB</i>	NM_003000.3	(Full-length)
<i>SDHC</i>	NM_003001.5	(Full-length)
<i>SDHD</i>	NM_003002.4	(Full-length)
<i>SETBP1</i>	NM_015559.3	-
<i>SETD2</i>	NM_014159.6	(Full-length)
<i>SF3B1</i>	NM_012433.4	-
<i>SH3RF2</i>	NM_152550.4	-
<i>SIX1</i>	NM_005982.4	-
<i>SIX2</i>	NM_016932.5	-
<i>SLC15A2</i>	NM_021082.4	-
<i>SLC8A1</i>	NM_021097.4	-
<i>SLCO1B3</i>	NM_019844.4	-
<i>SLX4</i>	NM_032444.4	(Full-length)
<i>SMAD2</i>	NM_001003652.4	(Full-length)
<i>SMAD4</i>	NM_005359.6	(Full-length)
<i>SMARCA4</i>	NM_001128849.3	(Full-length)
<i>SMARCB1</i>	NM_003073.5	(Full-length)
<i>SMCIA</i>	NM_006306.4	-
<i>SMO</i>	NM_005631.5	-
<i>SNCAIP</i>	NM_005460.4	-
<i>SOCS1</i>	NM_003745.1	(Full-length)
<i>SOS1</i>	NM_005633.3	-
<i>SOX2</i>	NM_003106.4	-
<i>SOX9</i>	NM_000346.4	(Full-length)
<i>SPEN</i>	NM_015001.3	(Full-length)
<i>SPOP</i>	NM_001007228.2	-
<i>SRC</i>	NM_198291.2	-
<i>SRSF2</i>	NM_003016.4	-
<i>STAG2</i>	NM_001042749.2	(Full-length)
<i>STAT1</i>	NM_007315.4	(Full-length)
<i>STAT3</i>	NM_139276.2	-
<i>STAT5B</i>	NM_012448.4	-
<i>STAT6</i>	NM_003153.5	-
<i>STK11</i>	NM_000455.5	(Full-length)
<i>SUFU</i>	NM_016169.4	(Full-length)
<i>SYT10</i>	NM_198992.4	-
<i>SYT16</i>	NM_001367661.1	-
<i>TAF1</i>	NM_004606.5	-
<i>TAP1</i>	NM_000593.5	(Full-length)
<i>TAP2</i>	NM_018833.2	(Full-length)

<b>Gene</b>	<b>Transcript ID</b>	<b>Full-Length Coverage</b>
<i>TAPBP</i>	NM_172208.2	-
<i>TBX3</i>	NM_016569.4	(Full-length)
<i>TCF7L2</i>	NM_001146274.2	(Full-length)
<i>TERT</i>	NM_198253.3	-
<i>TET2</i>	NM_001127208.2	(Full-length)
<i>TFE3</i>	NM_006521.6	-
<i>TFEB</i>	NM_001167827.3	-
<i>TGFBR1</i>	NM_004612.4	-
<i>TGFBR2</i>	NM_001024847.2	(Full-length)
<i>TMEM132D</i>	NM_133448.3	(Full-length)
<i>TNFAIP3</i>	NM_001270507.2	(Full-length)
<i>TNFRSF14</i>	NM_003820.3	(Full-length)
<i>TOP1</i>	NM_003286.4	-
<i>TOP2A</i>	NM_001067.4	-
<i>TP53</i>	NM_000546.5	(Full-length)
<i>TP63</i>	NM_003722.5	(Full-length)
<i>TPMT</i>	NM_000367.4	-
<i>TPP2</i>	NM_003291.4	(Full-length)
<i>TPTE</i>	NM_199261.3	-
<i>TRHDE</i>	NM_013381.2	-
<i>TRIM48</i>	NM_024114.5	-
<i>TRIM51</i>	NM_032681.4	-
<i>TRRAP</i>	NM_001244580.1	-
<i>TSC1</i>	NM_000368.5	(Full-length)
<i>TSC2</i>	NM_000548.5	(Full-length)
<i>TSHR</i>	NM_000369.4	-
<i>U2AF1</i>	NM_006758.2	-
<i>UGT1A1</i>	NM_000463.3	(Full-length)
<i>USP8</i>	NM_005154.5	-
<i>USP9X</i>	NM_001039590.3	(Full-length)
<i>VHL</i>	NM_000551.4	(Full-length)
<i>WAS</i>	NM_000377.3	-
<i>WT1</i>	NM_024426.6	(Full-length)
<i>XPO1</i>	NM_003400.4	-
<i>XRCC2</i>	NM_005431.2	(Full-length)
<i>XRCC3</i>	NM_001100118.2	(Full-length)
<i>YAP1</i>	NM_001130145.3	-
<i>YES1</i>	NM_005433.4	-
<i>ZBTB20</i>	NM_001164342.2	(Full-length)
<i>ZFH3</i>	NM_006885.4	(Full-length)
<i>ZIM3</i>	NM_052882.1	-
<i>ZMYM3</i>	NM_201599.3	(Full-length)
<i>ZNF217</i>	NM_006526.2	-

<b>Gene</b>	<b>Transcript ID</b>	<b>Full-Length Coverage</b>
<i>ZNF429</i>	NM_001001415.4	-
<i>ZNF479</i>	NM_033273.2	-
<i>ZNF536</i>	NM_014717.3	-
<i>ZRSR2</i>	NM_005089.3	(Full-length)

**Appendix B. List of Excluded Regions in the CellDx-Tissue Assay Panel**

<b>Gene</b>	<b>Transcript ID</b>	<b>Chromosome Coordinates</b>	<b>Exon</b>	<b>Amino Acid</b>
<i>AICF</i>	NM_138932.2	10:52570909-52571036	11	451_460
<i>ACVR1B</i>	NM_020328.4	12:52378938-52379062	7	368_397
<i>ADAMTS2</i>	NM_014244.5	5:178580490-178580630	9	461_505
<i>APC</i>	NM_000038.6	5:112111416-112111495	5	171_177
<i>APC</i>	NM_000038.6	5:112175027-112175161	16	1245_1290
<i>APC</i>	NM_000038.6	5:112176400-112176535	16	1703_1748
<i>ARID1A</i>	NM_006015.6	1:27023349-27023486	1	152_197
<i>ARID1A</i>	NM_006015.6	1:27023804-27023940	1	304_349
<i>ARID1A</i>	NM_006015.6	1:27023554-27023693	1	221_266
<i>ARID1A</i>	NM_006015.6	1:27023113-27023274	1	74_127
<i>ARID1A</i>	NM_006015.6	1:27023248-27023395	1	119_167
<i>ARID1A</i>	NM_006015.6	1:27023895-27024032	1	334_379
<i>ARID1A</i>	NM_006015.6	1:27022975-27023123	1	28_76
<i>ARID1B</i>	NM_001371656.1	6:157099473-157099631	2	220_272
<i>ARID1B</i>	NM_001371656.1	6:157099605-157099761	2	264_316
<i>ARID1B</i>	NM_001371656.1	6:157099056-157099193	2	81_126
<i>ARID1B</i>	NM_001371656.1	6:157100200-157100359	2	463_515
<i>ARID1B</i>	NM_001371656.1	6:157100339-157100475	2	509_554
<i>ARID1B</i>	NM_001371656.1	6:157100427-157100540	2	538_575
<i>ARID2</i>	NM_152641.4	12:46215139-46215251	6	213_229
<i>ARID2</i>	NM_152641.4	12:46123677-46123796	1	19_31
<i>ARMC4</i>	NM_018076.5	10:28142255-28142356	13	534_547
<i>ASXL1</i>	NM_015338.6	20:31017638-31017723	8	189_195
<i>ASXL2</i>	NM_018263.6	2:26068393-26068530	2	20_33
<i>ATM</i>	NM_000051.3	11:108188033-108188127	43	2067_2076
<i>ATM</i>	NM_000051.3	11:108098390-108098524	3	25_32
<i>ATR</i>	NM_001184.4	3:142215112-142215235	34	1955_1966
<i>ATRX</i>	NM_000489.5	X:76891365-76891489	16	1539_1567
<i>ATRX</i>	NM_000489.5	X:76944398-76944496	7	162_169
<i>AXIN1</i>	NM_003502.4	16:338020-338157	11	851_862
<i>BLM</i>	NM_000057.4	15:91295168-91295270	4	317_320
<i>BRIP1</i>	NM_032043.3	17:59763375-59763500	19	868_910
<i>C6</i>	NM_000065.4	5:41201590-41201712	3	84_101

<b>Gene</b>	<b>Transcript ID</b>	<b>Chromosome Coordinates</b>	<b>Exon</b>	<b>Amino Acid</b>
<i>C8B</i>	NM_000066.4	1:57395076-57395210	12	549_592
<i>CBFB</i>	NM_022845.3	16:67063653-67063773	2	34_55
<i>CCND1</i>	NM_053056.3	11:69457730-69457842	2	67_81
<i>CD276</i>	NM_001024736.2	15:73994649-73994786	3	45_90
<i>CDH10</i>	NM_006727.5	5:24491838-24491965	11	542_575
<i>CDH10</i>	NM_006727.5	5:24535170-24535295	5	247_272
<i>CDKN2A</i>	NM_001195132.1	9:21994162-21994298	intronic	intronic
<i>CHD4</i>	NM_001273.5	12:6697015-6697117	24	1156_1189
<i>CHEK2</i>	NM_007194.4	22:29126414-29126556	intronic	intronic
<i>CIC</i>	NM_015125.4	19:42797804-42797935	16	1286_1329
<i>CIC</i>	NM_015125.4	19:42799149-42799314	20	1545_1599
<i>CIC</i>	NM_015125.4	19:42798952-42799102	20	1487_1529
<i>CIC</i>	NM_015125.4	19:42778280-42778419	intronic	intronic
<i>CNTNAP4</i>	NM_138994.5	16:76523547-76523669	12	556_585
<i>CNTNAP4</i>	NM_138994.5	16:76513305-76513422	11	516_551
<i>CSMD3</i>	NM_198123.2	8:113516097-113516215	30	1632_1669
<i>CSMD3</i>	NM_198123.2	8:113668349-113668473	18	972_1002
<i>CSMD3</i>	NM_198123.2	8:113841909-113842029	12	585_620
<i>CUL3</i>	NM_003590.5	2:225449580-225449705	1	7_22
<i>CUL4A</i>	NM_001008895.4	13:113863944-113864054	1	2_38
<i>CUL4A</i>	NM_001008895.4	13:113863976-113864079	1	12_46
<i>CYLD</i>	NM_001042355.2	16:50788298-50788400	5	292_305
<i>CYLD</i>	NM_001042355.2	16:50785802-50785947	4	264_269
<i>DDX3X</i>	NM_001356.5	X:41206876-41206991	17	637_662
<i>DNMT3A</i>	NM_022552.4	2:25459815-25459907	21	803_823
<i>DOCK3</i>	NM_004947.5	3:51394323-51394443	44	1502_1519
<i>DSCI</i>	NM_024421.2	18:28710584-28710709	16	830_860
<i>EPCAM</i>	NM_002354.3	2:47596574-47596708	1	1_22
<i>EPHA2</i>	NM_004431.5	1:16482351-16482467	1	1_27
<i>EPHA2</i>	NM_004431.5	1:16464547-16464685	5	327_372
<i>FANCA</i>	NM_000135.4	16:89882944-89883082	1	1_27
<i>FANCA</i>	NM_000135.4	16:89882874-89883013	1	4_27
<i>FANCM</i>	NM_020937.4	14:45624508-45624620	8	437_452
<i>FANCM</i>	NM_020937.4	14:45654539-45654654	18	1545_1558

<b>Gene</b>	<b>Transcript ID</b>	<b>Chromosome Coordinates</b>	<b>Exon</b>	<b>Amino Acid</b>
<i>FANCM</i>	NM_020937.4	14:45665689-45665788	21	1885_1906
<i>FGF19</i>	NM_005117.3	11:69517956-69518095	2	91_113
<i>FGF19</i>	NM_005117.3	11:69518462-69518558	1	30_62
<i>FGF3</i>	NM_005247.4	11:69633479-69633593	1	37_74
<i>FGF4</i>	NM_002007.4	11:69589456-69589591	1	88_114
<i>FGF4</i>	NM_002007.4	11:69589747-69589885	1	1_36
<i>FGF4</i>	NM_002007.4	11:69589661-69589802	1	18_65
<i>FGFR3</i>	NM_000142.4	4:1803463-1803603	7	248_262
<i>FGFR3</i>	NM_000142.4	4:1803463-1803603	7	247_261
<i>FLT4</i>	NM_182925.5	5:180039486-180039619	26	1144_1179
<i>FOXA1</i>	NM_004496.5	14:38060946-38061083	2	303_348
<i>GLI3</i>	NM_000168.6	7:42005740-42005858	15	939_978
<i>GNAI3</i>	NM_006572.6	17:63052330-63052447	1	89_95
<i>HCN1</i>	NM_021072.4	5:45262354-45262493	8	736_782
<i>HCN1</i>	NM_021072.4	5:45262187-45262327	8	791_837
<i>HCN1</i>	NM_021072.4	5:45695914-45696022	1	59_94
<i>HCN1</i>	NM_021072.4	5:45696048-45696189	1	3_50
<i>HDAC2</i>	NM_001527.4	6:114291926-114292056	1	6_18
<i>HDAC2</i>	NM_001527.4	6:114262826-114262939	13	460_479
<i>INPP4B</i>	NM_001101669.3	4:143129669-143129753	14	323_328
<i>JAK1</i>	NM_002227.4	1:65306894-65307021	19	852_883
<i>KDM6A</i>	NM_021140.3	X:44921920-44922026	15	485_509
<i>KDM6A</i>	NM_021140.3	X:44910947-44911016	9	219_240
<i>KEL</i>	NM_000420.3	7:142651264-142651401	8	266_309
<i>KLF5</i>	NM_001730.5	13:73633434-73633575	1	1_37
<i>KLF5</i>	NM_001730.5	13:73633567-73633707	1	35_81
<i>KMT2A</i>	NM_001197104.2	11:118342828-118342961	3	318_362
<i>KMT2A</i>	NM_001197104.2	11:118307485-118307648	1	87_140
<i>KMT2B</i>	NM_014727.3	19:36209127-36209285	1	69_121
<i>KMT2B</i>	NM_014727.3	19:36211486-36211594	3	413_448
<i>KMT2B</i>	NM_014727.3	19:36212068-36212218	3	607_656
<i>KMT2C</i>	NM_170606.3	7:152132733-152132891	1	1_47
<i>KMT2C</i>	NM_170606.3	7:151899993-151900085	26	1342_1364
<i>KMT2D</i>	NM_003482.4	12:49439932-49440058	18	1528_1537

<b>Gene</b>	<b>Transcript ID</b>	<b>Chromosome Coordinates</b>	<b>Exon</b>	<b>Amino Acid</b>
<i>LATS2</i>	NM_014572.3	13:21562584-21562722	4	399_445
<i>LATS2</i>	NM_014572.3	13:21562457-21562596	4	441_487
<i>MAP2K4</i>	NM_003010.4	17:11924214-11924350	1	4_39
<i>MAP2K4</i>	NM_003010.4	17:11924163-11924301	1	1_33
<i>MAP3K1</i>	NM_005921.2	5:56111532-56111687	1	45_96
<i>MAP3K1</i>	NM_005921.2	5:56111787-56111917	1	129_161
<i>MAP3K1</i>	NM_005921.2	5:56111472-56111587	1	25_62
<i>MAP3K1</i>	NM_005921.2	5:56111682-56111823	1	95_141
<i>MAP3K4</i>	NM_005922.4	6:161412930-161413070	1	1_36
<i>MAP3K4</i>	NM_005922.4	6:161412974-161413104	1	4_47
<i>MAP3K4</i>	NM_005922.4	6:161523738-161523852	19	1261_1286
<i>MCL1</i>	NM_021960.5	1:150551662-150551785	1	75_116
<i>MCL1</i>	NM_021960.5	1:150551871-150552002	1	2_46
<i>MEN1</i>	NM_000244.3	11:64572116-64572254	10	467_513
<i>MGA</i>	NM_001164273.1	15:41991035-41991118	4	672_691
<i>MLH1</i>	NM_000249.4	3:37067067-37067175	12	347_363
<i>MLH3</i>	NM_001040108.2	14:75506521-75506638	5	1182_1190
<i>MSH6</i>	NM_000179.3	2:48023229-48023327	intronic	intronic
<i>MSH6</i>	NM_000179.3	2:48032173-48032292	intronic	intronic
<i>MSH6</i>	NM_000179.3	2:48010430-48010570	1	20_66
<i>MYCN</i>	NM_005378.6	2:16082871-16082980	2	228_264
<i>NOTCH1</i>	NM_017617.5	9:139417228-139417354	4	231_248
<i>NOTCH1</i>	NM_017617.5	9:139401067-139401203	24	1301_1309
<i>NOTCH1</i>	NM_017617.5	9:139440064-139440226	1	5_21
<i>NOTCH1</i>	NM_017617.5	9:139396907-139397042	28	1723_1734
<i>NOTCH3</i>	NM_000435.3	19:15281367-15281501	26	1625_1631
<i>NOTCH3</i>	NM_000435.3	19:15271813-15271950	33	2164_2209
<i>NOTCH3</i>	NM_000435.3	19:15288737-15288876	24	1289_1334
<i>NOTCH3</i>	NM_000435.3	19:15288502-15288662	24	1360_1413
<i>NOTCH3</i>	NM_000435.3	19:15288601-15288767	24	1325_1380
<i>NOTCH3</i>	NM_000435.3	19:15311504-15311669	1	17_40
<i>NOTCH3</i>	NM_000435.3	19:15281148-15281279	27	1660_1703
<i>NRXN1</i>	NM_004801.5	2:50318374-50318501	19	1228_1241
<i>OR2G6</i>	NM_001013355.1	1:248685252-248685391	2	102_148

<b>Gene</b>	<b>Transcript ID</b>	<b>Chromosome Coordinates</b>	<b>Exon</b>	<b>Amino Acid</b>
<i>OR2L13</i>	NM_175911.3	1:248262843-248262981	3	56_101
<i>OR2L8</i>	NM_001001963.1	1:248113049-248113177	1	297_312
<i>OR2T4</i>	NM_001004696.1	1:248525216-248525351	1	112_516
<i>OR4A15</i>	NM_001005275.1	11:55136184-55136276	1	276_306
<i>OR4C15</i>	NM_001001920.2	11:55322096-55322225	1	1_94
<i>OR8H2</i>	NM_001005200.1	11:55872770-55872877	2	85_120
<i>OR8H2</i>	NM_001005200.1	11:55873244-55873352	2	243_278
<i>PARP4</i>	NM_006437.4	13:25009001-25009088	31	1398_1426
<i>PARP4</i>	NM_006437.4	13:25008757-25008887	31	1465_1507
<i>PARP4</i>	NM_006437.4	13:25033254-25033346	20	782_784
<i>PARP4</i>	NM_006437.4	13:25016031-25016129	30	1182_1207
<i>PARP4</i>	NM_006437.4	13:25020802-25020919	27	1095_1122
<i>PGD</i>	NM_002631.4	1:10464224-10464351	5	112_150
<i>PLCG1</i>	NM_002660.3	20:39766235-39766366	1	1_29
<i>POLD1</i>	NM_001256849.1	19:50918899-50919029	22	906_922
<i>POLE</i>	NM_006231.4	12:133263801-133263936	1	1_21
<i>POT1</i>	NM_015450.3	7:124493139-124493201	10	235_253
<i>PPP2R2A</i>	NM_002717.4	8:26149297-26149379	1	1_4
<i>PPP2R2A</i>	NM_002717.4	8:26150699-26150821	intronic	intronic
<i>PTCH1</i>	NM_000264.5	9:98270633-98270780	1	1_4
<i>PTCH1</i>	NM_000264.5	9:98270481-98270635	1	4_54
<i>PTEN</i>	NM_000314.8	10:89623876-89624006	intronic	intronic
<i>PTEN</i>	NM_000314.8	10:89623658-89623799	intronic	intronic
<i>RAD52</i>	NM_134424.4	12:1040414-1040536	3	29_53
<i>RAD52</i>	NM_134424.4	12:1036026-1036137	intronic	intronic
<i>RASA1</i>	NM_002890.3	5:86637114-86637222	6	342_350
<i>RASA2</i>	NM_006506.5	3:141205871-141206008	1	1_28
<i>RASA2</i>	NM_006506.5	3:141292020-141292098	13	439_453
<i>RASA2</i>	NM_006506.5	3:141300012-141300124	17	582_584
<i>RASA2</i>	NM_006506.5	3:141299927-141300017	17	559_584
<i>RBI</i>	NM_000321.2	13:48947486-48947588	12	376_392
<i>RBI</i>	NM_000321.2	13:48878051-48878191	1	1_46
<i>RBM10</i>	NM_001204468.1	X:47004746-47004881	1	1_13
<i>RBM10</i>	NM_001204468.1	X:47004871-47005006	1	9_24

<b>Gene</b>	<b>Transcript ID</b>	<b>Chromosome Coordinates</b>	<b>Exon</b>	<b>Amino Acid</b>
<i>RBM10</i>	NM_001204468.1	X:47006725-47006861	2	24_59
<i>RNASEH2B</i>	NM_024570.4	13:51517410-51517498	6	146_160
<i>RNASEH2B</i>	NM_024570.4	13:51484144-51484276	1	1_22
<i>RNASEH2B</i>	NM_024570.4	13:51484233-51484361	1	7_22
<i>rs1052325</i>	NM_001206997.1	4:81884679-81884751	7	223_246
<i>rs11009107</i>	NM_001026383.3	10:33123687-33123810	intronic	intronic
<i>rs11551095</i>	NM_020170.4	19:3186076-3186170	1	17_48
<i>rs12876018</i>	NM_020121.4	13:96540151-96540223	26	989_1012
<i>rs16864976</i>	NM_003469.5	2:224463781-224463901	2	34_73
<i>rs2229482</i>	NM_005529.7	1:22199162-22199284	32	1298_1328
<i>rs2295997</i>	NM_003738.5	1:45293652-45293762	14	605_641
<i>rs7044405</i>	NM_001004487.1	9:35869933-35870039	1	121_155
<i>rs8024779</i>	NM_178232.4	15:89436216-89436352	intronic	intronic
<i>RUNX1</i>	NM_001754.4	21:36164356-36164476	9	466_480
<i>SDHA</i>	NM_004168.4	5:218411-218535	1	1_22
<i>SETD2</i>	NM_014159.6	3:47205247-47205408	1	3_24
<i>SMARCA4</i>	NM_001128849.3	19:11098403-11098543	6	307_354
<i>SMARCB1</i>	NM_003073.5	22:24129272-24129377	1	1_8
<i>SMARCB1</i>	NM_003073.5	22:24129375-24129509	1	6_31
<i>SOCS1</i>	NM_003745.1	16:11348772-11348913	2	142_188
<i>SOCS1</i>	NM_003745.1	16:11349150-11349288	2	17_62
<i>SOX9</i>	NM_000346.4	17:70120011-70120150	3	338_384
<i>SPEN</i>	NM_015001.3	1:16262462-16262599	11	3242_3288
<i>STAG2</i>	NM_001042749.2	X:123202371-123202482	24	756_779
<i>SUFU</i>	NM_016169.4	10:104263902-104264024	1	1_39
<i>TBX3</i>	NM_016569.4	12:115112290-115112427	7	438_483
<i>TBX3</i>	NM_016569.4	12:115112075-115112231	7	503_555
<i>TCF7L2</i>	NM_001146274.2	10:114711012-114711156	2	79_86
<i>TMEM132D</i>	NM_133448.3	12:129822116-129822256	4	407_433
<i>TP63</i>	NM_003722.5	3:189455542-189455660	2	25_64
<i>TPP2</i>	NM_003291.4	13:103328654-103328765	28	1183_1220
<i>TSC1</i>	NM_000368.5	9:135772955-135773050	21	876_890
<i>TSC2</i>	NM_000548.5	16:2106675-2106815	8	226_258
<i>USP9X</i>	NM_001039590.3	X:41055607-41055677	27	1360_1362

<b>Gene</b>	<b>Transcript ID</b>	<b>Chromosome Coordinates</b>	<b>Exon</b>	<b>Amino Acid</b>
<i>USP9X</i>	NM_001039590.3	X:41077608-41077731	37	2070_2105
<i>USP9X</i>	NM_001039590.3	X:41026738-41026852	17	777_808
<i>USP9X</i>	NM_001039590.3	X:41057979-41058081	30	1526_1535
<i>WT1</i>	NM_024426.6	11:32456482-32456627	1	94_142
<i>WT1</i>	NM_024426.6	11:32456860-32456965	1	1_16
<i>XRCC3</i>	NM_001100118.2	14:104173409-104173542	5	69_113
<i>YAP1</i>	NM_001130145.3	11:101981593-101981759	1	5_60
<i>ZFH3</i>	NM_006885.4	16:72821562-72821701	10	3492_3538
<i>ZMYM3</i>	NM_201599.3	X:70466231-70466357	15	807_849
<i>ZMYM3</i>	NM_201599.3	X:70469263-70469340	7	480_490
<i>ZMYM3</i>	NM_201599.3	X:70472622-70472733	2	126_162
<i>ZMYM3</i>	NM_201599.3	X:70472683-70472778	2	111_142
<i>ZNF217</i>	NM_006526.2	20:52193028-52193147	4	720_759
<i>ZRSR2</i>	NM_005089.3	X:15841002-15841112	11	362_399

**Appendix C. List of Excluded Genes in the CellDx-Tissue Assay Panel**

<b>Gene</b>	<b>Transcript ID</b>	<b>Chromosome Coordinates</b>	<b>Exon</b>	<b>Amino Acid</b>
<i>AICF</i>	NM_138932.2	10:52570909-52571036	11	451_460
<i>ACVR1B</i>	NM_020328.4	12:52378938-52379062	7	368_397
<i>ADAMTS2</i>	NM_014244.5	5:178580490-178580630	9	461_505
<i>APC</i>	NM_000038.6	5:112111416-112111495	5	171_177
<i>APC</i>	NM_000038.6	5:112175027-112175161	16	1245_1290
<i>APC</i>	NM_000038.6	5:112176400-112176535	16	1703_1748
<i>ARID1A</i>	NM_006015.6	1:27023349-27023486	1	152_197
<i>ARID1A</i>	NM_006015.6	1:27023804-27023940	1	304_349
<i>ARID1A</i>	NM_006015.6	1:27023554-27023693	1	221_266
<i>ARID1A</i>	NM_006015.6	1:27023113-27023274	1	74_127
<i>ARID1A</i>	NM_006015.6	1:27023248-27023395	1	119_167
<i>ARID1A</i>	NM_006015.6	1:27023895-27024032	1	334_379
<i>ARID1A</i>	NM_006015.6	1:27022975-27023123	1	28_76
<i>ARID1B</i>	NM_001371656.1	6:157099473-157099631	2	220_272
<i>ARID1B</i>	NM_001371656.1	6:157099605-157099761	2	264_316
<i>ARID1B</i>	NM_001371656.1	6:157099056-157099193	2	81_126
<i>ARID1B</i>	NM_001371656.1	6:157100200-157100359	2	463_515
<i>ARID1B</i>	NM_001371656.1	6:157100339-157100475	2	509_554
<i>ARID1B</i>	NM_001371656.1	6:157100427-157100540	2	538_575
<i>ARID2</i>	NM_152641.4	12:46215139-46215251	6	213_229
<i>ARID2</i>	NM_152641.4	12:46123677-46123796	1	19_31
<i>ARMC4</i>	NM_018076.5	10:28142255-28142356	13	534_547
<i>ASXL1</i>	NM_015338.6	20:31017638-31017723	8	189_195
<i>ASXL2</i>	NM_018263.6	2:26068393-26068530	2	20_33
<i>ATM</i>	NM_000051.3	11:108188033-108188127	43	2067_2076
<i>ATM</i>	NM_000051.3	11:108098390-108098524	3	25_32
<i>ATR</i>	NM_001184.4	3:142215112-142215235	34	1955_1966
<i>ATRX</i>	NM_000489.5	X:76891365-76891489	16	1539_1567
<i>ATRX</i>	NM_000489.5	X:76944398-76944496	7	162_169
<i>AXIN1</i>	NM_003502.4	16:338020-338157	11	851_862
<i>BLM</i>	NM_000057.4	15:91295168-91295270	4	317_320
<i>BRIP1</i>	NM_032043.3	17:59763375-59763500	19	868_910
<i>C6</i>	NM_000065.4	5:41201590-41201712	3	84_101
<i>C8B</i>	NM_000066.4	1:57395076-57395210	12	549_592
<i>CBFB</i>	NM_022845.3	16:67063653-67063773	2	34_55
<i>CCND1</i>	NM_053056.3	11:69457730-69457842	2	67_81
<i>CD276</i>	NM_001024736.2	15:73994649-73994786	3	45_90
<i>CDH10</i>	NM_006727.5	5:24491838-24491965	11	542_575
<i>CDH10</i>	NM_006727.5	5:24535170-24535295	5	247_272
<i>CDKN2A</i>	NM_001195132.1	9:21994162-21994298	intronic	intronic
<i>CHD4</i>	NM_001273.5	12:6697015-6697117	24	1156_1189
<i>CHEK2</i>	NM_007194.4	22:29126414-29126556	intronic	intronic

<b>Gene</b>	<b>Transcript ID</b>	<b>Chromosome Coordinates</b>	<b>Exon</b>	<b>Amino Acid</b>
<i>CIC</i>	NM_015125.4	19:42797804-42797935	16	1286_1329
<i>CIC</i>	NM_015125.4	19:42799149-42799314	20	1545_1599
<i>CIC</i>	NM_015125.4	19:42798952-42799102	20	1487_1529
<i>CIC</i>	NM_015125.4	19:42778280-42778419	intronic	intronic
<i>CNTNAP4</i>	NM_138994.5	16:76523547-76523669	12	556_585
<i>CNTNAP4</i>	NM_138994.5	16:76513305-76513422	11	516_551
<i>CSMD3</i>	NM_198123.2	8:113516097-113516215	30	1632_1669
<i>CSMD3</i>	NM_198123.2	8:113668349-113668473	18	972_1002
<i>CSMD3</i>	NM_198123.2	8:113841909-113842029	12	585_620
<i>CUL3</i>	NM_003590.5	2:225449580-225449705	1	7_22
<i>CUL4A</i>	NM_001008895.4	13:113863944-113864054	1	2_38
<i>CUL4A</i>	NM_001008895.4	13:113863976-113864079	1	12_46
<i>CYLD</i>	NM_001042355.2	16:50788298-50788400	5	292_305
<i>CYLD</i>	NM_001042355.2	16:50785802-50785947	4	264_269
<i>DDX3X</i>	NM_001356.5	X:41206876-41206991	17	637_662
<i>DNMT3A</i>	NM_022552.4	2:25459815-25459907	21	803_823
<i>DOCK3</i>	NM_004947.5	3:51394323-51394443	44	1502_1519
<i>DSCI</i>	NM_024421.2	18:28710584-28710709	16	830_860
<i>EPCAM</i>	NM_002354.3	2:47596574-47596708	1	1_22
<i>EPHA2</i>	NM_004431.5	1:16482351-16482467	1	1_27
<i>EPHA2</i>	NM_004431.5	1:16464547-16464685	5	327_372
<i>FANCA</i>	NM_000135.4	16:89882944-89883082	1	1_27
<i>FANCA</i>	NM_000135.4	16:89882874-89883013	1	4_27
<i>FANCM</i>	NM_020937.4	14:45624508-45624620	8	437_452
<i>FANCM</i>	NM_020937.4	14:45654539-45654654	18	1545_1558
<i>FANCM</i>	NM_020937.4	14:45665689-45665788	21	1885_1906
<i>FGF19</i>	NM_005117.3	11:69517956-69518095	2	91_113
<i>FGF19</i>	NM_005117.3	11:69518462-69518558	1	30_62
<i>FGF3</i>	NM_005247.4	11:69633479-69633593	1	37_74
<i>FGF4</i>	NM_002007.4	11:69589456-69589591	1	88_114
<i>FGF4</i>	NM_002007.4	11:69589747-69589885	1	1_36
<i>FGF4</i>	NM_002007.4	11:69589661-69589802	1	18_65
<i>FGFR3</i>	NM_000142.4	4:1803463-1803603	7	248_262
<i>FGFR3</i>	NM_000142.4	4:1803463-1803603	7	247_261
<i>FLT4</i>	NM_182925.5	5:180039486-180039619	26	1144_1179
<i>FOXA1</i>	NM_004496.5	14:38060946-38061083	2	303_348
<i>GLI3</i>	NM_000168.6	7:42005740-42005858	15	939_978
<i>GNAI3</i>	NM_006572.6	17:63052330-63052447	1	89_95
<i>HCN1</i>	NM_021072.4	5:45262354-45262493	8	736_782
<i>HCN1</i>	NM_021072.4	5:45262187-45262327	8	791_837
<i>HCN1</i>	NM_021072.4	5:45695914-45696022	1	59_94
<i>HCN1</i>	NM_021072.4	5:45696048-45696189	1	3_50
<i>HDAC2</i>	NM_001527.4	6:114291926-114292056	1	6_18

<b>Gene</b>	<b>Transcript ID</b>	<b>Chromosome Coordinates</b>	<b>Exon</b>	<b>Amino Acid</b>
<i>HDAC2</i>	NM_001527.4	6:114262826-114262939	13	460_479
<i>INPP4B</i>	NM_001101669.3	4:143129669-143129753	14	323_328
<i>JAK1</i>	NM_002227.4	1:65306894-65307021	19	852_883
<i>KDM6A</i>	NM_021140.3	X:44921920-44922026	15	485_509
<i>KDM6A</i>	NM_021140.3	X:44910947-44911016	9	219_240
<i>KEL</i>	NM_000420.3	7:142651264-142651401	8	266_309
<i>KLF5</i>	NM_001730.5	13:73633434-73633575	1	1_37
<i>KLF5</i>	NM_001730.5	13:73633567-73633707	1	35_81
<i>KMT2A</i>	NM_001197104.2	11:118342828-118342961	3	318_362
<i>KMT2A</i>	NM_001197104.2	11:118307485-118307648	1	87_140
<i>KMT2B</i>	NM_014727.3	19:36209127-36209285	1	69_121
<i>KMT2B</i>	NM_014727.3	19:36211486-36211594	3	413_448
<i>KMT2B</i>	NM_014727.3	19:36212068-36212218	3	607_656
<i>KMT2C</i>	NM_170606.3	7:152132733-152132891	1	1_47
<i>KMT2C</i>	NM_170606.3	7:151899993-151900085	26	1342_1364
<i>KMT2D</i>	NM_003482.4	12:49439932-49440058	18	1528_1537
<i>LATS2</i>	NM_014572.3	13:21562584-21562722	4	399_445
<i>LATS2</i>	NM_014572.3	13:21562457-21562596	4	441_487
<i>MAP2K4</i>	NM_003010.4	17:11924214-11924350	1	4_39
<i>MAP2K4</i>	NM_003010.4	17:11924163-11924301	1	1_33
<i>MAP3K1</i>	NM_005921.2	5:56111532-56111687	1	45_96
<i>MAP3K1</i>	NM_005921.2	5:56111787-56111917	1	129_161
<i>MAP3K1</i>	NM_005921.2	5:56111472-56111587	1	25_62
<i>MAP3K1</i>	NM_005921.2	5:56111682-56111823	1	95_141
<i>MAP3K4</i>	NM_005922.4	6:161412930-161413070	1	1_36
<i>MAP3K4</i>	NM_005922.4	6:161412974-161413104	1	4_47
<i>MAP3K4</i>	NM_005922.4	6:161523738-161523852	19	1261_1286
<i>MCL1</i>	NM_021960.5	1:150551662-150551785	1	75_116
<i>MCL1</i>	NM_021960.5	1:150551871-150552002	1	2_46
<i>MEN1</i>	NM_000244.3	11:64572116-64572254	10	467_513
<i>MGA</i>	NM_001164273.1	15:41991035-41991118	4	672_691
<i>MLH1</i>	NM_000249.4	3:37067067-37067175	12	347_363
<i>MLH3</i>	NM_001040108.2	14:75506521-75506638	5	1182_1190
<i>MSH6</i>	NM_000179.3	2:48023229-48023327	intronic	intronic
<i>MSH6</i>	NM_000179.3	2:48032173-48032292	intronic	intronic
<i>MSH6</i>	NM_000179.3	2:48010430-48010570	1	20_66
<i>MYCN</i>	NM_005378.6	2:16082871-16082980	2	228_264
<i>NOTCH1</i>	NM_017617.5	9:139417228-139417354	4	231_248
<i>NOTCH1</i>	NM_017617.5	9:139401067-139401203	24	1301_1309
<i>NOTCH1</i>	NM_017617.5	9:139440064-139440226	1	5_21
<i>NOTCH1</i>	NM_017617.5	9:139396907-139397042	28	1723_1734
<i>NOTCH3</i>	NM_000435.3	19:15281367-15281501	26	1625_1631
<i>NOTCH3</i>	NM_000435.3	19:15271813-15271950	33	2164_2209

Gene	Transcript ID	Chromosome Coordinates	Exon	Amino Acid
<i>NOTCH3</i>	NM_000435.3	19:15288737-15288876	24	1289_1334
<i>NOTCH3</i>	NM_000435.3	19:15288502-15288662	24	1360_1413
<i>NOTCH3</i>	NM_000435.3	19:15288601-15288767	24	1325_1380
<i>NOTCH3</i>	NM_000435.3	19:15311504-15311669	1	17_40
<i>NOTCH3</i>	NM_000435.3	19:15281148-15281279	27	1660_1703
<i>NRXN1</i>	NM_004801.5	2:50318374-50318501	19	1228_1241
<i>OR2G6</i>	NM_001013355.1	1:248685252-248685391	2	102_148
<i>OR2L13</i>	NM_175911.3	1:248262843-248262981	3	56_101
<i>OR2L8</i>	NM_001001963.1	1:248113049-248113177	1	297_312
<i>OR2T4</i>	NM_001004696.1	1:248525216-248525351	1	112_516
<i>OR4A15</i>	NM_001005275.1	11:55136184-55136276	1	276_306
<i>OR4C15</i>	NM_001001920.2	11:55322096-55322225	1	1_94
<i>OR8H2</i>	NM_001005200.1	11:55872770-55872877	2	85_120
<i>OR8H2</i>	NM_001005200.1	11:55873244-55873352	2	243_278
<i>PARP4</i>	NM_006437.4	13:25009001-25009088	31	1398_1426
<i>PARP4</i>	NM_006437.4	13:25008757-25008887	31	1465_1507
<i>PARP4</i>	NM_006437.4	13:25033254-25033346	20	782_784
<i>PARP4</i>	NM_006437.4	13:25016031-25016129	30	1182_1207
<i>PARP4</i>	NM_006437.4	13:25020802-25020919	27	1095_1122
<i>PGD</i>	NM_002631.4	1:10464224-10464351	5	112_150
<i>PLCG1</i>	NM_002660.3	20:39766235-39766366	1	1_29
<i>POLD1</i>	NM_001256849.1	19:50918899-50919029	22	906_922
<i>POLE</i>	NM_006231.4	12:133263801-133263936	1	1_21
<i>POT1</i>	NM_015450.3	7:124493139-124493201	10	235_253
<i>PPP2R2A</i>	NM_002717.4	8:26149297-26149379	1	1_4
<i>PPP2R2A</i>	NM_002717.4	8:26150699-26150821	intronic	intronic
<i>PTCHI</i>	NM_000264.5	9:98270633-98270780	1	1_4
<i>PTCHI</i>	NM_000264.5	9:98270481-98270635	1	4_54
<i>PTEN</i>	NM_000314.8	10:89623876-89624006	intronic	intronic
<i>PTEN</i>	NM_000314.8	10:89623658-89623799	intronic	intronic
<i>RAD52</i>	NM_134424.4	12:1040414-1040536	3	29_53
<i>RAD52</i>	NM_134424.4	12:1036026-1036137	intronic	intronic
<i>RASA1</i>	NM_002890.3	5:86637114-86637222	6	342_350
<i>RASA2</i>	NM_006506.5	3:141205871-141206008	1	1_28
<i>RASA2</i>	NM_006506.5	3:141292020-141292098	13	439_453
<i>RASA2</i>	NM_006506.5	3:141300012-141300124	17	582_584
<i>RASA2</i>	NM_006506.5	3:141299927-141300017	17	559_584
<i>RBI</i>	NM_000321.2	13:48947486-48947588	12	376_392
<i>RBI</i>	NM_000321.2	13:48878051-48878191	1	1_46
<i>RBM10</i>	NM_001204468.1	X:47004746-47004881	1	1_13
<i>RBM10</i>	NM_001204468.1	X:47004871-47005006	1	9_24
<i>RBM10</i>	NM_001204468.1	X:47006725-47006861	2	24_59
<i>RNASEH2B</i>	NM_024570.4	13:51517410-51517498	6	146_160

Gene	Transcript ID	Chromosome Coordinates	Exon	Amino Acid
<i>RNASEH2B</i>	NM_024570.4	13:51484144-51484276	1	1_22
<i>RNASEH2B</i>	NM_024570.4	13:51484233-51484361	1	7_22
<i>rs1052325</i>	NM_001206997.1	4:81884679-81884751	7	223_246
<i>rs11009107</i>	NM_001026383.3	10:33123687-33123810	intronic	intronic
<i>rs11551095</i>	NM_020170.4	19:3186076-3186170	1	17_48
<i>rs12876018</i>	NM_020121.4	13:96540151-96540223	26	989_1012
<i>rs16864976</i>	NM_003469.5	2:224463781-224463901	2	34_73
<i>rs2229482</i>	NM_005529.7	1:22199162-22199284	32	1298_1328
<i>rs2295997</i>	NM_003738.5	1:45293652-45293762	14	605_641
<i>rs7044405</i>	NM_001004487.1	9:35869933-35870039	1	121_155
<i>rs8024779</i>	NM_178232.4	15:89436216-89436352	intronic	intronic
<i>RUNX1</i>	NM_001754.4	21:36164356-36164476	9	466_480
<i>SDHA</i>	NM_004168.4	5:218411-218535	1	1_22
<i>SETD2</i>	NM_014159.6	3:47205247-47205408	1	3_24
<i>SMARCA4</i>	NM_001128849.3	19:11098403-11098543	6	307_354
<i>SMARCB1</i>	NM_003073.5	22:24129272-24129377	1	1_8
<i>SMARCB1</i>	NM_003073.5	22:24129375-24129509	1	6_31
<i>SOCS1</i>	NM_003745.1	16:11348772-11348913	2	142_188
<i>SOCS1</i>	NM_003745.1	16:11349150-11349288	2	17_62
<i>SOX9</i>	NM_000346.4	17:70120011-70120150	3	338_384
<i>SPEN</i>	NM_015001.3	1:16262462-16262599	11	3242_3288
<i>STAG2</i>	NM_001042749.2	X:123202371-123202482	24	756_779
<i>SUFU</i>	NM_016169.4	10:104263902-104264024	1	1_39
<i>TBX3</i>	NM_016569.4	12:115112290-115112427	7	438_483
<i>TBX3</i>	NM_016569.4	12:115112075-115112231	7	503_555
<i>TCF7L2</i>	NM_001146274.2	10:114711012-114711156	2	79_86
<i>TMEM132D</i>	NM_133448.3	12:129822116-129822256	4	407_433
<i>TP63</i>	NM_003722.5	3:189455542-189455660	2	25_64
<i>TPP2</i>	NM_003291.4	13:103328654-103328765	28	1183_1220
<i>TSC1</i>	NM_000368.5	9:135772955-135773050	21	876_890
<i>TSC2</i>	NM_000548.5	16:2106675-2106815	8	226_258
<i>USP9X</i>	NM_001039590.3	X:41055607-41055677	27	1360_1362
<i>USP9X</i>	NM_001039590.3	X:41077608-41077731	37	2070_2105
<i>USP9X</i>	NM_001039590.3	X:41026738-41026852	17	777_808
<i>USP9X</i>	NM_001039590.3	X:41057979-41058081	30	1526_1535
<i>WT1</i>	NM_024426.6	11:32456482-32456627	1	94_142
<i>WT1</i>	NM_024426.6	11:32456860-32456965	1	1_16
<i>XRCC3</i>	NM_001100118.2	14:104173409-104173542	5	69_113
<i>YAP1</i>	NM_001130145.3	11:101981593-101981759	1	5_60
<i>ZFH3</i>	NM_006885.4	16:72821562-72821701	10	3492_3538
<i>ZMYM3</i>	NM_201599.3	X:70466231-70466357	15	807_849
<i>ZMYM3</i>	NM_201599.3	X:70469263-70469340	7	480_490

<b>Gene</b>	<b>Transcript ID</b>	<b>Chromosome Coordinates</b>	<b>Exon</b>	<b>Amino Acid</b>
<i>ZMYM3</i>	NM_201599.3	X:70472622-70472733	2	126_162
<i>ZMYM3</i>	NM_201599.3	X:70472683-70472778	2	111_142
<i>ZNF217</i>	NM_006526.2	20:52193028-52193147	4	720_759
<i>ZRSR2</i>	NM_005089.3	X:15841002-15841112	11	362_399

**Appendix D. Precision Summary of CellDx – Tissue Assay for Selected SNVs, Insertions, and Deletions**

<b>Gene Exon</b>	<b>Mutation (cDNA/AA Change)</b>	<b>Normalized Coverage Rance</b>	<b>VAF Range</b>	<b>VAF Mean</b>	<b>SD</b>	<b>%CV</b>	<b>n/N</b>	<b>PCR (%) (95% CI)</b>
<i>AMER1</i> Exon 2	c.2185G>T p.E729*	0.48-0.77	0.554-0.773	0.641	0.037	5.8	48/48	100 (92.6, 100)
<i>APC</i> Exon 16	c.2138C>G p.S713*	0.66-1.09	0.109-0.196	0.148	0.020	13.8	48/48	100 (92.6, 100)
<i>APC</i> Exon 16	c.4348C>T p.R1450*	0.94-1.68	0.186-0.286	0.233	0.023	10.1	48/48	100 (92.6, 100)
<i>APC</i> Exon 9	c.847C>T p.R283*	0.66-1.41	0.154-0.277	0.234	0.023	9.9	48/48	100 (92.6, 100)
<i>ASXL2</i> Exon 3	c.143G>A p.S48N	0.82-1.51	0.274-0.351	0.313	0.016	5.2	48/48	100 (92.6, 100)
<i>ATM</i> Exon 47	c.6908delA p.K2303Rfs*7	0.84-2.28	0.287-0.381	0.324	0.020	6.3	42/48	88 (75.3, 94.1)
<i>ATM</i> Exon 8	c.1009C>T p.R337C	0.55-1.19	0.256-0.392	0.326	0.028	8.6	48/48	100 (92.6, 100)
<i>AXIN1</i> Exon 11	c.2463-14C>T	1.23-1.93	0.436-0.533	0.491	0.017	3.6	48/48	100 (92.6, 100)
<i>B2M</i> Exon 1	c.3G>A	0.79-1.56	0.034-0.086	0.057	0.012	21.3	47/48	98 (89.1, 99.6)
<i>B2M</i> Exon 1	c.43_44delCT p.L15Ffs*41	0.69-1.45	0.035-0.102	0.063	0.012	19.8	48/48	100 (92.6, 100)
<i>B2M</i> Exon 2	c.244_247delTTCT p.F82Ifs*20	0.69-1.45	0.052-0.107	0.068	0.011	16.0	45/48	94 (83.2, 97.9)
<i>B2M</i> Exon 2	c.142_145delTCTG p.S48Gfs*12	0.94-1.73	0.274-0.363	0.322	0.017	5.4	48/48	100 (92.6, 100)
<i>B2M</i> Exon 2	c.68-1G>A	0.94-1.73	0.281-0.361	0.323	0.020	6.1	48/48	100 (92.6, 100)
<i>BRAF</i> Exon 15	c.1799_1800delTGinsAA p.V600E	0.39-0.83	0.359-0.47	0.415	0.023	5.5	48/48	100 (92.6, 100)
<i>BRCA2</i> Exon 11	c.4588A>T p.K1530*	0.36-0.78	0.413-0.546	0.481	0.030	6.3	46/46	100 (92.3, 100)

Gene Exon	Mutation (cDNA/AA Change)	Normalized Coverage Rance	VAF Range	VAF Mean	SD	%CV	n/N	PCR (%) (95% CI)
<i>BRCA2</i> Exon 20	c.8581A>T p.R2861*	1.20-2.00	0.738-0.796	0.772	0.013	1.7	48/48	100 (92.6, 100)
<i>CIC</i> Exon 14	c.3216_3217delGT p.S1073Tfs*77	0.27-0.776	0.702-0.855	0.773	0.035	4.6	48/48	100 (92.6, 100)
<i>CUL3</i> Exon 8	c.1162_1163delCT p.S389Ifs*4	0.37-1.35	0.093-0.182	0.132	0.020	15.2	48/48	100 (92.6, 100)
<i>DNMT3A</i> Exon 23	c.2701delC p.L901Sfs*5	0.78-1.25	0.05-0.13	0.072	0.014,	19.4	46/48	96 (86.0, 98.8)
<i>DNMT3A</i> Exon 6	c.541C>T p.R181C	0.05-0.68	0.227-0.55	0.331	0.086	26.0	31/31	100 (89.0, 100)
<i>EGFR</i> Exon 19	c.2236_2250delGAATTAAGAGAAGCA p.E746_A750del	0.62-1.55	0.201-0.333	0.237	0.023	9.6	48/48	100 (92.6, 100)
<i>ERBB2</i> Exon 19	c.2305G>T p.D769Y	0.76-1.64	0.265-0.407	0.342	0.032	9.4	46/46	100 (92.3, 100)
<i>ERBB3</i> Exon 23	c.2783A>G p.E928G	1.07-2.55	0.024-0.042	0.03	0.004	12.9	44/46	96 (84.6, 99.0)
<i>ERCC2</i> Exon 21	c.1975C>A p.H659N	0.08-0.47	0.262-0.412	0.352	0.030	8.6	48/48	100 (92.6, 100)
<i>ERCC4</i> Exon 8	c.1811+1G>A	0.86-1.36	0.269-0.382	0.331	0.024	7.3	48/48	100 (92.6, 100)
<i>ETV6</i> Exon 4	c.427_428delCA p.Q143Afs*10	0.48-1.27	0.27-0.367	0.322	0.023	7.3	48/48	100 (92.6, 100)
<i>FANCA</i> Exon 24	c.2222+7G>A	0.32-0.83	0.387-0.585	0.49	0.041	8.4	48/48	100 (92.6, 100)
<i>FANCA</i> Exon 38	c.3828+1G>C	0.25-0.75	0.686-0.795	0.75	0.020	2.7	48/48	100 (92.6, 100)
<i>FAT1</i> Exon 26	c.13002A>C p.T4334T	1.16-1.92	0.175-0.26	0.208	0.018	8.7	48/48	100 (92.6, 100)
<i>FUBP1</i> Exon 14	c.1213C>T p.Q405*	0.49-1.25	0.265-0.351	0.31	0.023	7.3	48/48	100 (92.6, 100)
<i>IDH1</i> Exon 4	c.395G>A p.R132H	1.60-2.19	0.414-0.503	0.458	0.019,	4.2	48/48	100 (92.6, 100)

Gene Exon	Mutation (cDNA/AA Change)	Normalized Coverage Rance	VAF Range	VAF Mean	SD	%CV	n/N	PCR (%) (95% CI)
<i>KEAP1</i> Exon 4	c.1531+8G>C	0.88-1.52	0.454-0.565	0.514	0.026,	5.1	48/48	100 (92.6, 100)
<i>KMT2A</i> Exon 3	c.1142delA p.K381Rfs*19	1.34-5.37	0.296-0.384	0.34	0.018	5.4	42/48	88 (75.3, 94.1)
<i>KMT2D</i> Exon 32	c.7061delC p.P2354Lfs*30	0.78-1.55	0.156-0.243	0.208	0.019	9.0	45/48	94 (83.2, 97.9)
<i>KMT2D</i> Exon 49	c.14710C>T p.R4904*	0.85-1.71	0.3-0.373	0.328	0.016	4.9	48/48	100 (92.6, 100)
<i>KRAS</i> Exon 2	c.35G>C p.G12A	1.03-2.52	0.023-0.061	0.039	0.008	20.4	44/46	96 (85.5, 98.8)
<i>KRAS</i> Exon 2	c.34G>T p.G12C	1.32-2.88	0.028-0.069	0.044	0.011	25.8	25/31	81 (63.7, 90.8)
<i>KRAS</i> Exon 2	c.38G>A p.G13D	1.32-2.88	0.281-0.354	0.307	0.018	6.0	31/31	100 (89.0, 100)
<i>KRAS</i> Exon 2	c.35G>T p.G12V	0.96-1.86	0.327-0.404	0.364	0.017	4.6	48/48	100 (92.6, 100)
<i>MAX</i> Exon 4	c.179G>A p.R60Q	0.86-2.06	0.286-0.369	0.325	0.021	6.4	48/48	100 (92.6, 100)
<i>MLH3</i> Exon 3	c.3367C>T p.Q1123*	1.34-2.07	0.415-0.573	0.489	0.039	8.0	31/31	100 (89.0, 100)
<i>MRE11</i> Exon 14	c.1563+1G>A	0.64-1.76	0.427-0.546	0.488	0.026	5.4	48/48	100 (92.6, 100)
<i>NBN</i> Exon 2	c.171+4T>C	0.76-1.44	0.419-0.538	0.474	0.025	5.3	46/46	100 (92.3, 100)
<i>NOTCH1</i> Exon 26	c.4777C>A p.L1593M	0.10-2.37	0.295-0.484	0.361	0.037	10.4	48/48	100 (92.6, 100)
<i>NOTCH3</i> Exon 4	c.348_349delCT p.D116Efs*45	0.71-1.90	0.328-0.46	0.388	0.028	7.3	47/48	98 (89.1, 99.6)
<i>PIK3CA</i> Exon 21	c.3140A>G p.H1047R	0.92-1.80	0.038-0.082	0.062	0.012	19.5	26/31	84 (67.4, 92.9)

Gene Exon	Mutation (cDNA/AA Change)	Normalized Coverage Rance	VAF Range	VAF Mean	SD	%CV	n/N	PCR (%) (95% CI)
<i>PIK3CA</i> Exon 21	c.3073A>G p.T1025A	0.93-1.58	0.194-0.273	0.233	0.016	6.7	48/48	100 (92.6, 100)
<i>PIK3CA</i> Exon 10	c.1624G>A p.E542K	0.25-1.60	0.502-0.608	0.564	0.018	3.2	48/48	100 (92.6, 100)
<i>POLD1</i> Exon 17	c.2066G>A p.R689Q	0.10-1.01	0.253-0.356	0.297	0.028	9.5	48/48	100 (92.6, 100)
<i>PPM1D</i> Exon 6	c.1281G>A p.W427*	0.62-1.02	0.245-0.359	0.303	0.026	8.4	48/48	100 (92.6, 100)
<i>PTEN</i> Exon 5	c.389G>A p.R130Q	0.71-1.09	0.032-0.115	0.063	0.014	22.7	48/48	100 (92.6, 100)
<i>PTEN</i> Exon 7	c.741_742insA p.P248Tfs*5	0.48-0.77	0.074-0.149	0.109	0.017	15.9	48/48	100 (92.6, 100)
<i>RAD51B</i> Exon 6	c.453-7C>T	0.65-1.29	0.514-0.645	0.573	0.027	4.7	48/48	100 (92.6, 100)
<i>RBM10</i> Exon 13	c.1444-2A>C	0.69-1.26	0.193-0.351	0.297	0.036	12.0	31/31	100 (89.0, 100)
<i>RECQL4</i> Exon 10	c.1704+6G>T	2.00-3.67	0.274-0.364	0.318	0.020	6.4	48/48	100 ( 92.6, 100)
<i>SETD2</i> Exon 20	c.7432-10C>T	0.92-1.31	0.103-0.205	0.144	0.020	14.1	48/48	100 (92.6, 100)
<i>SOX9</i> Exon 3	c.1280_1281dup p.Y428Tfs*43	0.40-1.27	0.096-0.145	0.117	0.010	8.8	48/48	100 (92.6, 100)
<i>TET2</i> Exon 3	c.1771C>T p.Q591*	0.42-1.16	0.051-0.084	0.07	0.008	11.7	47/48	98 (89.1, 99.6)
<i>TGFBR2</i> Exon 6	c.1411G>A p.D471N	0.62-1.33	0.297-0.366	0.334	0.014	4.3	48/48	100 (92.6, 100)
<i>TP53</i> Exon 6	c.614A>G p.Y205C	0.40-1.08	0.209-0.3	0.253	0.019	7.5	48/48	100 (92.6, 100)
<i>TP53</i> Exon 7	c.742C>T p.R248W	0.59-1.14	0.437-0.537	0.492	0.022	4.5	48/48	100 (92.6, 100)
<i>TP53</i> Exon 5	c.536A>G p.H179R	1.10-2.45	0.301-0.563	0.507	0.037	7.3	45/48	94 (83.2, 97.9)

<b>Gene Exon</b>	<b>Mutation (cDNA/AA Change)</b>	<b>Normalized Coverage Rance</b>	<b>VAF Range</b>	<b>VAF Mean</b>	<b>SD</b>	<b>%CV</b>	<b>n/N</b>	<b>PCR (%) (95% CI)</b>
<i>TSC2</i> Exon 27	c.3131+2T>A	0.68-2.98	0.299-0.394	0.344	0.023	6.6	48/48	100 (92.6, 100)

## Appendix E. Accuracy

### Appendix E.1 Concordance for SNV by Gene

Gene	Total Bases	Total Number of Variant per Gene	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
<i>ABRAXAS1</i>	1831	0	-	100% (100%, 100%)
<i>ACVR2A</i>	2610	0	-	100% (100%, 100%)
<i>AKT1</i>	1023	4	100% (51.01%, 100%)	100% (100%, 100%)
<i>AMER1</i>	3440	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>APC</i>	10461	22	100% (81.57%, 100%)	100% (100%, 100%)
<i>AR</i>	1003	1	-	100% (100%, 100%)
<i>ARHGAP35</i>	4709	0	-	100% (100%, 100%)
<i>ARID1A</i>	8463	11	80% (37.55%, 96.38%)	100% (100%, 100%)
<i>ARID1B</i>	8590	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>ARID2</i>	7469	10	100% (64.57%, 100%)	100% (100%, 100%)
<i>ARID5B</i>	4391	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>ASXL1</i>	5582	0	-	100% (100%, 100%)
<i>ASXL2</i>	5476	1	-	100% (100%, 100%)
<i>ATM</i>	15017	19	85.71% (48.69%, 97.43%)	100% (100%, 100%)
<i>ATR</i>	12585	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>ATRX</i>	11246	1	-	100% (100%, 100%)
<i>AXIN1</i>	3615	1	-	100% (100%, 100%)
<i>AXL</i>	1177	1	-	100% (100%, 100%)
<i>B2M</i>	769	4	100% (43.85%, 100%)	100% (100%, 100%)
<i>BAP1</i>	3550	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>BARD1</i>	3430	0	-	100% (100%, 100%)
<i>BCL6</i>	1911	1	-	100% (100%, 100%)
<i>BCOR</i>	6319	1	-	100% (100%, 100%)
<i>BRAF</i>	953	4	100% (43.85%, 100%)	100% (100%, 100%)
<i>BRCA1</i>	9233	3	100% (34.24%, 100%)	100% (100%, 100%)
<i>BRCA2</i>	13279	6	100% (56.55%, 100%)	100% (100%, 100%)
<i>BRIP1</i>	5546	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>CARD11</i>	1920	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>CASP8</i>	2487	3	100% (34.24%, 100%)	100% (100%, 100%)
<i>CDC73</i>	3019	0	-	100% (100%, 100%)
<i>CDH1</i>	4202	4	100% (34.24%, 100%)	100% (100%, 100%)
<i>CDK12</i>	5966	1	-	100% (100%, 100%)
<i>CDK4</i>	1251	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>CDKN1A</i>	1101	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>CDKN1B</i>	770	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>CDKN2A</i>	1416	6	80% (37.55%, 96.38%)	100% (100%, 100%)
<i>CDKN2B</i>	932	1	-	100% (100%, 100%)
<i>CHEK1</i>	2457	1	-	100% (100%, 100%)
<i>CHEK2</i>	3053	7	100% (43.85%, 100%)	100% (100%, 100%)
<i>CIC</i>	9290	4	100% (20.65%, 100%)	100% (100%, 100%)

Gene	Total Bases	Total Number of Variant per Gene	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
<i>CREBBP</i>	9902	4	100% (34.24%, 100%)	100% (100%, 100%)
<i>CTCF</i>	3086	3	100% (34.24%, 100%)	100% (100%, 100%)
<i>CTNNB1</i>	500	3	100% (43.85%, 100%)	100% (99.99%, 100%)
<i>CUL3</i>	4404	5	100% (34.24%, 100%)	100% (100%, 100%)
<i>CYLD</i>	4307	0	-	100% (100%, 100%)
<i>DAXX</i>	2919	1	-	100% (100%, 100%)
<i>DDR2</i>	893	2	-	100% (99.99%, 100%)
<i>DDX3X</i>	3080	3	100% (34.24%, 100%)	100% (100%, 100%)
<i>DICER1</i>	7883	4	100% (34.24%, 100%)	100% (100%, 100%)
<i>DNMT3A</i>	5341	6	100% (20.65%, 100%)	100% (100%, 100%)
<i>EGFR</i>	2022	9	100% (60.97%, 100%)	100% (100%, 100%)
<i>ELF3</i>	1686	0	-	100% (100%, 100%)
<i>EP300</i>	9915	5	100% (56.55%, 100%)	100% (100%, 100%)
<i>ERBB2</i>	3025	6	100% (51.01%, 100%)	100% (100%, 100%)
<i>ERBB3</i>	965	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>ERCC2</i>	4456	7	100% (56.55%, 100%)	100% (100%, 100%)
<i>ERCC4</i>	3838	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>ETV6</i>	1999	0	-	100% (100%, 100%)
<i>FANCA</i>	8471	6	100% (51.01%, 100%)	100% (100%, 100%)
<i>FANCC</i>	2965	0	-	100% (100%, 100%)
<i>FANCD2</i>	7785	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>FANCE</i>	2595	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>FANCF</i>	1239	0	-	100% (100%, 100%)
<i>FANCG</i>	3089	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>FANCI</i>	7067	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>FANCL</i>	2170	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>FAT1</i>	16016	5	100% (56.55%, 100%)	100% (100%, 100%)
<i>FBXW7</i>	3767	7	100% (43.85%, 100%)	100% (100%, 100%)
<i>FGF23</i>	1080	1	-	100% (100%, 100%)
<i>FGF3</i>	955	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>FGFR3</i>	1656	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>FGFR4</i>	1463	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>FOXA1</i>	920	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>FUBP1</i>	3648	1	-	100% (100%, 100%)
<i>GATA2</i>	1018	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>GATA3</i>	1902	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>GNAS</i>	1121	4	100% (43.85%, 100%)	100% (100%, 100%)
<i>HLA-A</i>	1588	3	100% (43.85%, 100%)	100% (100%, 100%)
<i>HLA-B</i>	1567	3	-	100% (100%, 100%)
<i>HNF1A</i>	2744	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>HRAS</i>	681	4	100% (43.85%, 100%)	100% (99.99%, 100%)
<i>IDH1</i>	125	3	100% (43.85%, 100%)	100% (99.98%, 100%)
<i>INPP4B</i>	4917	1	0% (0.0%, 79.35%)	100% (100%, 100%)

Gene	Total Bases	Total Number of Variant per Gene	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
<i>JAK1</i>	5570	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>JAK3</i>	5374	1	-	100% (100%, 100%)
<i>KDM5C</i>	7047	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>KDM6A</i>	7324	3	100% (20.65%, 100%)	100% (100%, 100%)
<i>KDR</i>	933	1	-	100% (100%, 100%)
<i>KEAP1</i>	2396	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>KEL</i>	1963	1	0% (0.0%, 79.35%)	100% (100%, 100%)
<i>KIT</i>	2012	3	100% (43.85%, 100%)	100% (100%, 100%)
<i>KMT2A</i>	14698	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>KMT2B</i>	11058	6	100% (56.55%, 100%)	100% (100%, 100%)
<i>KMT2C</i>	19770	3	100% (34.24%, 100%)	100% (100%, 100%)
<i>KMT2D</i>	20922	9	100% (67.56%, 100%)	100% (100%, 100%)
<i>KRAS</i>	929	30	96.43% (82.29%, 99.37%)	100% (99.99%, 100%)
<i>LATS1</i>	4133	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>LATS2</i>	4047	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>MAP2K1</i>	632	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>MAP3K1</i>	6425	3	100% (34.24%, 100%)	100% (100%, 100%)
<i>MAPK1</i>	1000	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>MAX</i>	1227	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>MEN1</i>	2606	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>MET</i>	1504	7	75.00% (30.06%, 95.44%)	100% (100%, 100%)
<i>MLH1</i>	5020	3	100% (20.65%, 100%)	100% (100%, 100%)
<i>MLH3</i>	5599	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>MPL</i>	946	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>MRE11</i>	3963	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>MSH2</i>	4576	3	100% (20.65%, 100%)	100% (100%, 100%)
<i>MSH3</i>	5529	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>MSH6</i>	5188	0	-	100% (100%, 100%)
<i>MTAP</i>	1674	0	-	100% (100%, 100%)
<i>MUTYH</i>	3016	1	-	100% (100%, 100%)
<i>MYC</i>	1154	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>NBN</i>	3558	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>NCOR1</i>	11641	0	-	100% (100%, 100%)
<i>NF1</i>	12882	11	100% (56.55%, 100%)	100% (100%, 100%)
<i>NF2</i>	3582	3	100% (34.24%, 100%)	100% (100%, 100%)
<i>NFE2L2</i>	936	3	100% (43.85%, 100%)	100% (100%, 100%)
<i>NOTCH1</i>	11747	5	100% (34.24%, 100%)	100% (100%, 100%)
<i>NOTCH2</i>	10611	3	100% (43.85%, 100%)	100% (100%, 100%)
<i>NOTCH3</i>	9364	5	100% (20.65%, 100%)	100% (100%, 100%)
<i>NOTCH4</i>	8205	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>NRAS</i>	829	5	100% (56.55%, 100%)	100% (100%, 100%)
<i>NTRK1</i>	1442	3	100% (20.65%, 100%)	100% (100%, 100%)
<i>PALB2</i>	4460	3	100% (43.85%, 100%)	100% (100%, 100%)

Gene	Total Bases	Total Number of Variant per Gene	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
<i>PARP1</i>	5183	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>PBRM1</i>	7314	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>PHF6</i>	1804	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>PIK3CA</i>	1750	29	100% (86.68%, 100%)	100% (100%, 100%)
<i>PIK3R1</i>	3282	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>PIK3R2</i>	1441	0	-	100% (100%, 100%)
<i>PMS1</i>	2271	0	-	100% (100%, 100%)
<i>PMS2</i>	4100	10	100% (34.24%, 100%)	100% (100%, 100%)
<i>POLD1</i>	5875	1	-	100% (100%, 100%)
<i>POLE</i>	11804	0	-	100% (100%, 100%)
<i>POT1</i>	3363	0	-	100% (100%, 100%)
<i>PPM1D</i>	2382	1	-	100% (100%, 100%)
<i>PPP2R1A</i>	1277	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>PPP6C</i>	974	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>PRDM1</i>	3323	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>PTCH1</i>	7082	0	-	100% (100%, 100%)
<i>PTEN</i>	2594	10	100% (70.09%, 100%)	100% (100%, 100%)
<i>PTPN11</i>	940	1	-	100% (100%, 100%)
<i>RAD50</i>	6073	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>RAD51B</i>	2324	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>RAD51C</i>	2038	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>RAD51D</i>	1997	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>RAF1</i>	1120	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>RASA1</i>	5763	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>RB1</i>	5388	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>RBM10</i>	5150	5	100% (43.85%, 100%)	100% (100%, 100%)
<i>RECQL4</i>	2856	4	100% (43.85%, 100%)	100% (100%, 100%)
<i>RET</i>	1424	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>RHOA</i>	162	5	100% (43.85%, 100%)	99.99% (99.97%, 100%)
<i>RNF43</i>	3231	0	-	100% (100%, 100%)
<i>SDHA</i>	3393	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>SDHC</i>	981	0	-	100% (100%, 100%)
<i>SETD2</i>	9510	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>SF3B1</i>	1083	2	-	100% (100%, 100%)
<i>SLX4</i>	6660	4	100% (34.24%, 100%)	100% (100%, 100%)
<i>SMAD4</i>	2846	4	100% (43.85%, 100%)	100% (100%, 100%)
<i>SMARCA4</i>	8148	6	100% (51.01%, 100%)	100% (100%, 100%)
<i>SMARCB1</i>	2129	2	-	100% (100%, 100%)
<i>SOX9</i>	1859	3	100% (34.24%, 100%)	100% (100%, 100%)
<i>SPEN</i>	12166	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>SPOP</i>	916	1	-	100% (100%, 100%)
<i>STAG2</i>	6507	8	100% (20.65%, 100%)	100% (100%, 100%)
<i>STAT3</i>	1023	1	100% (20.65%, 100%)	100% (100%, 100%)

Gene	Total Bases	Total Number of Variant per Gene	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
<i>STK11</i>	2331	6	100% (43.85%, 100%)	100% (100%, 100%)
<i>TBX3</i>	2732	0	-	100% (100%, 100%)
<i>TCF7L2</i>	3618	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>TERT</i>	1206	17	80% (54.81%, 92.95%)	100% (100%, 100%)
<i>TET2</i>	7164	1	-	100% (100%, 100%)
<i>TGFBR2</i>	2501	3	100% (34.24%, 100%)	100% (100%, 100%)
<i>TNFAIP3</i>	2903	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>TNFRSF14</i>	1637	1	-	100% (100%, 100%)
<i>TP53</i>	2344	74	97.01% (89.75%, 99.18%)	100% (100%, 100%)
<i>TP63</i>	1397	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>TSC1</i>	5349	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>TSC2</i>	10092	13	100% (72.25%, 100%)	100% (100%, 100%)
<i>U2AF1</i>	1008	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>UGT1A1</i>	146	4	100% (51.01%, 100%)	100% (99.98%, 100%)
<i>VHL</i>	889	1	-	100% (100%, 100%)
<i>WT1</i>	2792	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>XRCC2</i>	1157	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>ZFHX3</i>	10946	4	100% (43.85%, 100%)	100% (100%, 100%)

## Appendix E.2. Concordance for Insertions by Gene

Gene	Total Bases	Total Number of Variant per Gene	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
<i>ABRAXAS1</i>	1831	0	-	100% (100%, 100%)
<i>ACVR2A</i>	2610	0	-	100% (100%, 100%)
<i>AKT1</i>	1023	0	-	100% (100%, 100%)
<i>AMER1</i>	3440	0	-	100% (100%, 100%)
<i>APC</i>	10461	6	80% (37.55%, 96.38%)	100% (100%, 100%)
<i>AR</i>	1003	0	-	100% (100%, 100%)
<i>ARHGAP35</i>	4709	1	-	100% (100%, 100%)
<i>ARID1A</i>	8463	1	-	100% (100%, 100%)
<i>ARID1B</i>	8590	0	-	100% (100%, 100%)
<i>ARID2</i>	7469	0	-	100% (100%, 100%)
<i>ARID5B</i>	4391	1	-	100% (100%, 100%)
<i>ASXL1</i>	5582	2	0% (0.0%, 65.76%)	100% (100%, 100%)
<i>ASXL2</i>	5476	0	-	100% (100%, 100%)
<i>ATM</i>	15017	10	-	100% (100%, 100%)
<i>ATR</i>	12585	0	-	100% (100%, 100%)
<i>ATRX</i>	11246	0	-	100% (100%, 100%)
<i>AXIN1</i>	3615	0	-	100% (100%, 100%)
<i>AXL</i>	1177	0	-	100% (100%, 100%)
<i>B2M</i>	769	0	-	100% (100%, 100%)
<i>BAP1</i>	3550	0	-	100% (100%, 100%)
<i>BARD1</i>	3430	1	0% (0.0%, 79.35%)	100% (100%, 100%)
<i>BCL6</i>	1911	0	-	100% (100%, 100%)
<i>BCOR</i>	6319	0	-	100% (100%, 100%)
<i>BRAF</i>	953	0	-	100% (100%, 100%)
<i>BRCA1</i>	9233	5	50% (9.45%, 90.55%)	100% (100%, 100%)
<i>BRCA2</i>	13279	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>BRIP1</i>	5546	0	-	100% (100%, 100%)
<i>CARD11</i>	1920	0	-	100% (100%, 100%)
<i>CASP8</i>	2487	0	-	100% (100%, 100%)
<i>CDC73</i>	3019	1	-	100% (100%, 100%)
<i>CDH1</i>	4202	1	0% (0.0%, 79.35%)	100% (100%, 100%)
<i>CDK12</i>	5966	0	-	100% (100%, 100%)
<i>CDK4</i>	1251	0	-	100% (100%, 100%)
<i>CDKN1A</i>	1101	0	-	100% (100%, 100%)
<i>CDKN1B</i>	770	0	-	100% (100%, 100%)
<i>CDKN2A</i>	1416	0	-	100% (100%, 100%)
<i>CDKN2B</i>	932	0	-	100% (100%, 100%)
<i>CHEK1</i>	2457	0	-	100% (100%, 100%)
<i>CHEK2</i>	3053	0	-	100% (100%, 100%)
<i>CIC</i>	9290	1	-	100% (100%, 100%)
<i>CREBBP</i>	9902	0	-	100% (100%, 100%)
<i>CTCF</i>	3086	1	0% (0.0%, 79.35%)	100% (100%, 100%)

Gene	Total Bases	Total Number of Variant per Gene	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
<i>CTNNB1</i>	500	0	-	100% (99.99%, 100%)
<i>CUL3</i>	4404	0	-	100% (100%, 100%)
<i>CYLD</i>	4307	1	-	100% (100%, 100%)
<i>DAXX</i>	2919	0	-	100% (100%, 100%)
<i>DDR2</i>	893	0	-	100% (100%, 100%)
<i>DDX3X</i>	3080	0	-	100% (100%, 100%)
<i>DICER1</i>	7883	1	-	100% (100%, 100%)
<i>DNMT3A</i>	5341	2	-	100% (100%, 100%)
<i>EGFR</i>	2022	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>ELF3</i>	1686	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>EP300</i>	9915	0	-	100% (100%, 100%)
<i>ERBB2</i>	3025	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>ERBB3</i>	965	0	-	100% (100%, 100%)
<i>ERCC2</i>	4456	0	-	100% (100%, 100%)
<i>ERCC4</i>	3838	0	-	100% (100%, 100%)
<i>ETV6</i>	1999	0	-	100% (100%, 100%)
<i>FANCA</i>	8471	0	-	100% (100%, 100%)
<i>FANCC</i>	2965	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>FANCD2</i>	7785	0	-	100% (100%, 100%)
<i>FANCE</i>	2595	0	-	100% (100%, 100%)
<i>FANCF</i>	1239	0	-	100% (100%, 100%)
<i>FANCG</i>	3089	0	-	100% (100%, 100%)
<i>FANCI</i>	7067	0	-	100% (100%, 100%)
<i>FANCL</i>	2170	0	-	100% (100%, 100%)
<i>FAT1</i>	16016	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>FBXW7</i>	3767	0	-	100% (100%, 100%)
<i>FGF23</i>	1080	0	-	100% (100%, 100%)
<i>FGF3</i>	955	0	-	100% (100%, 100%)
<i>FGFR3</i>	1656	0	-	100% (100%, 100%)
<i>FGFR4</i>	1463	0	-	100% (100%, 100%)
<i>FOXA1</i>	920	0	-	100% (100%, 100%)
<i>FUBP1</i>	3648	1	-	100% (100%, 100%)
<i>GATA2</i>	1018	0	-	100% (100%, 100%)
<i>GATA3</i>	1902	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>GNAS</i>	1121	0	-	100% (100%, 100%)
<i>HLA-A</i>	1588	6	100% (20.65%, 100%)	100% (100%, 100%)
<i>HLA-B</i>	1567	1	-	100% (100%, 100%)
<i>HNF1A</i>	2744	1	-	100% (100%, 100%)
<i>HRAS</i>	681	0	-	100% (100%, 100%)
<i>IDH1</i>	125	0	-	100% (99.98%, 100%)
<i>INPP4B</i>	4917	1	-	100% (100%, 100%)
<i>JAK1</i>	5570	0	-	100% (100%, 100%)
<i>JAK3</i>	5374	0	-	100% (100%, 100%)

Gene	Total Bases	Total Number of Variant per Gene	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
<i>KDM5C</i>	7047	3	-	100% (100%, 100%)
<i>KDM6A</i>	7324	1	-	100% (100%, 100%)
<i>KDR</i>	933	0	-	100% (100%, 100%)
<i>KEAP1</i>	2396	0	-	100% (100%, 100%)
<i>KEL</i>	1963	0	-	100% (100%, 100%)
<i>KIT</i>	2012	0	-	100% (100%, 100%)
<i>KMT2A</i>	14698	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>KMT2B</i>	11058	0	-	100% (100%, 100%)
<i>KMT2C</i>	19770	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>KMT2D</i>	20922	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>KRAS</i>	929	0	-	100% (100%, 100%)
<i>LATS1</i>	4133	2	-	100% (100%, 100%)
<i>LATS2</i>	4047	0	-	100% (100%, 100%)
<i>MAP2K1</i>	632	0	-	100% (100%, 100%)
<i>MAP3K1</i>	6425	5	100% (51.01%, 100%)	100% (100%, 100%)
<i>MAPK1</i>	1000	0	-	100% (100%, 100%)
<i>MAX</i>	1227	0	-	100% (100%, 100%)
<i>MEN1</i>	2606	0	-	100% (100%, 100%)
<i>MET</i>	1504	1	-	100% (100%, 100%)
<i>MLH1</i>	5020	0	-	100% (100%, 100%)
<i>MLH3</i>	5599	0	-	100% (100%, 100%)
<i>MPL</i>	946	0	-	100% (100%, 100%)
<i>MRE11</i>	3963	0	-	100% (100%, 100%)
<i>MSH2</i>	4576	1	-	100% (100%, 100%)
<i>MSH3</i>	5529	0	-	100% (100%, 100%)
<i>MSH6</i>	5188	1	-	100% (100%, 100%)
<i>MTAP</i>	1674	1	-	100% (100%, 100%)
<i>MUTYH</i>	3016	0	-	100% (100%, 100%)
<i>MYC</i>	1154	0	-	100% (100%, 100%)
<i>NBN</i>	3558	1	-	100% (100%, 100%)
<i>NCOR1</i>	11641	1	-	100% (100%, 100%)
<i>NF1</i>	12882	0	-	100% (100%, 100%)
<i>NF2</i>	3582	0	-	100% (100%, 100%)
<i>NFE2L2</i>	936	0	-	100% (100%, 100%)
<i>NOTCH1</i>	11747	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>NOTCH2</i>	10611	1	-	100% (100%, 100%)
<i>NOTCH3</i>	9364	0	-	100% (100%, 100%)
<i>NOTCH4</i>	8205	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>NRAS</i>	829	0	-	100% (100%, 100%)
<i>NTRK1</i>	1442	0	-	100% (100%, 100%)
<i>PALB2</i>	4460	3	0% (0.0%, 79.35%)	100% (100%, 100%)
<i>PARP1</i>	5183	1	-	100% (100%, 100%)
<i>PBRM1</i>	7314	1	100% (20.65%, 100%)	100% (100%, 100%)

Gene	Total Bases	Total Number of Variant per Gene	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
<i>PHF6</i>	1804	0	-	100% (100%, 100%)
<i>PIK3CA</i>	1750	0	-	100% (100%, 100%)
<i>PIK3R1</i>	3282	1	-	100% (100%, 100%)
<i>PIK3R2</i>	1441	0	-	100% (100%, 100%)
<i>PMS1</i>	2271	0	-	100% (100%, 100%)
<i>PMS2</i>	4100	1	0% (0.0%, 79.35%)	100% (100%, 100%)
<i>POLD1</i>	5875	0	-	100% (100%, 100%)
<i>POLE</i>	11804	0	-	100% (100%, 100%)
<i>POT1</i>	3363	2	-	100% (100%, 100%)
<i>PPM1D</i>	2382	0	-	100% (100%, 100%)
<i>PPP2R1A</i>	1277	0	-	100% (100%, 100%)
<i>PPP6C</i>	974	0	-	100% (100%, 100%)
<i>PRDM1</i>	3323	0	-	100% (100%, 100%)
<i>PTCHI</i>	7082	0	-	100% (100%, 100%)
<i>PTEN</i>	2594	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>PTPN11</i>	940	0	-	100% (100%, 100%)
<i>RAD50</i>	6073	0	-	100% (100%, 100%)
<i>RAD51B</i>	2324	0	-	100% (100%, 100%)
<i>RAD51C</i>	2038	0	-	100% (100%, 100%)
<i>RAD51D</i>	1997	0	-	100% (100%, 100%)
<i>RAF1</i>	1120	0	-	100% (100%, 100%)
<i>RASA1</i>	5763	0	-	100% (100%, 100%)
<i>RBI</i>	5388	1	-	100% (100%, 100%)
<i>RBM10</i>	5150	0	-	100% (100%, 100%)
<i>RECQL4</i>	2856	0	-	100% (100%, 100%)
<i>RET</i>	1424	0	-	100% (100%, 100%)
<i>RHOA</i>	162	0	-	100% (99.98%, 100%)
<i>RNF43</i>	3231	0	-	100% (100%, 100%)
<i>SDHA</i>	3393	1	-	100% (100%, 100%)
<i>SDHC</i>	981	1	-	100% (100%, 100%)
<i>SETD2</i>	9510	1	-	100% (100%, 100%)
<i>SF3B1</i>	1083	0	-	100% (100%, 100%)
<i>SLX4</i>	6660	2	-	100% (100%, 100%)
<i>SMAD4</i>	2846	0	-	100% (100%, 100%)
<i>SMARCA4</i>	8148	0	-	100% (100%, 100%)
<i>SMARCB1</i>	2129	0	-	100% (100%, 100%)
<i>SOX9</i>	1859	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>SPEN</i>	12166	0	-	100% (100%, 100%)
<i>SPOP</i>	916	0	-	100% (100%, 100%)
<i>STAG2</i>	6507	2	-	100% (100%, 100%)
<i>STAT3</i>	1023	0	-	100% (100%, 100%)
<i>STK11</i>	2331	0	-	100% (100%, 100%)
<i>TBX3</i>	2732	1	100% (20.65%, 100%)	100% (100%, 100%)

<b>Gene</b>	<b>Total Bases</b>	<b>Total Number of Variant per Gene</b>	<b>PPA (2-sided 95% CI)</b>	<b>NPA (2-sided 95% CI)</b>
<i>TCF7L2</i>	3618	0	-	100% (100%, 100%)
<i>TERT</i>	1206	0	-	100% (100%, 100%)
<i>TET2</i>	7164	1	-	100% (100%, 100%)
<i>TGFBR2</i>	2501	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>TNFAIP3</i>	2903	0	-	100% (100%, 100%)
<i>TNFRSF14</i>	1637	0	-	100% (100%, 100%)
<i>TP53</i>	2344	4	100% (43.85%, 100%)	100% (100%, 100%)
<i>TP63</i>	1397	0	-	100% (100%, 100%)
<i>TSC1</i>	5349	1	-	100% (100%, 100%)
<i>TSC2</i>	10092	0	-	100% (100%, 100%)
<i>U2AF1</i>	1008	0	-	100% (100%, 100%)
<i>UGT1A1</i>	146	0	-	100% (99.98%, 100%)
<i>VHL</i>	889	0	-	100% (100%, 100%)
<i>WT1</i>	2792	1	-	100% (100%, 100%)
<i>XRCC2</i>	1157	0	-	100% (100%, 100%)
<i>ZFHX3</i>	10946	1	-	100% (100%, 100%)

### Appendix E.3. Concordance for Deletions by Gene

Gene	Total Bases	Total Number of Variant per Gene	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
<i>ABRAXAS1</i>	1831	1	-	100% (100%, 100%)
<i>ACVR2A</i>	2610	2	-	100% (100%, 100%)
<i>AKT1</i>	1023	0	-	100% (100%, 100%)
<i>AMER1</i>	3440	0	-	100% (100%, 100%)
<i>APC</i>	10461	8	80% (37.55%, 96.38%)	100% (100%, 100%)
<i>AR</i>	1003	0	-	100% (100%, 100%)
<i>ARHGAP35</i>	4709	1	-	100% (100%, 100%)
<i>ARID1A</i>	8463	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>ARID1B</i>	8590	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>ARID2</i>	7469	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>ARID5B</i>	4391	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>ASXL1</i>	5582	0	-	100% (100%, 100%)
<i>ASXL2</i>	5476	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>ATM</i>	15017	4	100% (20.65%, 100%)	100% (100%, 100%)
<i>ATR</i>	12585	0	-	100% (100%, 100%)
<i>ATRX</i>	11246	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>AXIN1</i>	3615	2	-	100% (100%, 100%)
<i>AXL</i>	1177	0	-	100% (100%, 100%)
<i>B2M</i>	769	6	100% (51.01%, 100%)	100% (99.99%, 100%)
<i>BAP1</i>	3550	0	-	100% (100%, 100%)
<i>BARD1</i>	3430	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>BCL6</i>	1911	0	-	100% (100%, 100%)
<i>BCOR</i>	6319	3	0% (0.0%, 79.35%)	100% (100%, 100%)
<i>BRAF</i>	953	0	-	100% (100%, 100%)
<i>BRCA1</i>	9233	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>BRCA2</i>	13279	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>BRIP1</i>	5546	0	-	100% (100%, 100%)
<i>CARD11</i>	1920	0	-	100% (100%, 100%)
<i>CASP8</i>	2487	0	-	100% (100%, 100%)
<i>CDC73</i>	3019	0	-	100% (100%, 100%)
<i>CDH1</i>	4202	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>CDK12</i>	5966	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>CDK4</i>	1251	0	-	100% (100%, 100%)
<i>CDKN1A</i>	1101	0	-	100% (100%, 100%)
<i>CDKN1B</i>	770	0	-	100% (100%, 100%)
<i>CDKN2A</i>	1416	3	100% (34.24%, 100%)	100% (100%, 100%)
<i>CDKN2B</i>	932	0	-	100% (100%, 100%)
<i>CHEK1</i>	2457	0	-	100% (100%, 100%)
<i>CHEK2</i>	3053	0	-	100% (100%, 100%)
<i>CIC</i>	9290	3	100% (43.85%, 100%)	100% (100%, 100%)
<i>CREBBP</i>	9902	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>CTCF</i>	3086	0	-	100% (100%, 100%)

Gene	Total Bases	Total Number of Variant per Gene	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
<i>CTNNB1</i>	500	0	-	100% (99.99%, 100%)
<i>CUL3</i>	4404	1	-	100% (100%, 100%)
<i>CYLD</i>	4307	0	-	100% (100%, 100%)
<i>DAXX</i>	2919	0	-	100% (100%, 100%)
<i>DDR2</i>	893	0	-	100% (100%, 100%)
<i>DDX3X</i>	3080	2	-	100% (100%, 100%)
<i>DICER1</i>	7883	0	-	100% (100%, 100%)
<i>DNMT3A</i>	5341	0	-	100% (100%, 100%)
<i>EGFR</i>	2022	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>ELF3</i>	1686	0	-	100% (100%, 100%)
<i>EP300</i>	9915	1	-	100% (100%, 100%)
<i>ERBB2</i>	3025	0	-	100% (100%, 100%)
<i>ERBB3</i>	965	0	-	100% (100%, 100%)
<i>ERCC2</i>	4456	0	-	100% (100%, 100%)
<i>ERCC4</i>	3838	0	-	100% (100%, 100%)
<i>ETV6</i>	1999	1	-	100% (100%, 100%)
<i>FANCA</i>	8471	1	-	100% (100%, 100%)
<i>FANCC</i>	2965	0	-	100% (100%, 100%)
<i>FANCD2</i>	7785	0	-	100% (100%, 100%)
<i>FANCE</i>	2595	0	-	100% (100%, 100%)
<i>FANCF</i>	1239	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>FANCG</i>	3089	0	-	100% (100%, 100%)
<i>FANCI</i>	7067	0	-	100% (100%, 100%)
<i>FANCL</i>	2170	0	-	100% (100%, 100%)
<i>FAT1</i>	16016	0	-	100% (100%, 100%)
<i>FBXW7</i>	3767	0	-	100% (100%, 100%)
<i>FGF23</i>	1080	0	-	100% (100%, 100%)
<i>FGF3</i>	955	0	-	100% (100%, 100%)
<i>FGFR3</i>	1656	0	-	100% (100%, 100%)
<i>FGFR4</i>	1463	0	-	100% (100%, 100%)
<i>FOXA1</i>	920	0	-	100% (100%, 100%)
<i>FUBP1</i>	3648	1	-	100% (100%, 100%)
<i>GATA2</i>	1018	0	-	100% (100%, 100%)
<i>GATA3</i>	1902	1	-	100% (100%, 100%)
<i>GNAS</i>	1121	0	-	100% (100%, 100%)
<i>HLA-A</i>	1588	0	-	100% (100%, 100%)
<i>HLA-B</i>	1567	2	-	100% (100%, 100%)
<i>HNF1A</i>	2744	3	0% (0.0%, 79.35%)	100% (100%, 100%)
<i>HRAS</i>	681	0	-	100% (100%, 100%)
<i>IDH1</i>	125	0	-	100% (99.98%, 100%)
<i>INPP4B</i>	4917	1	-	100% (100%, 100%)
<i>JAK1</i>	5570	0	-	100% (100%, 100%)
<i>JAK3</i>	5374	0	-	100% (100%, 100%)

Gene	Total Bases	Total Number of Variant per Gene	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
<i>KDM5C</i>	7047	0	-	100% (100%, 100%)
<i>KDM6A</i>	7324	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>KDR</i>	933	0	-	100% (100%, 100%)
<i>KEAP1</i>	2396	0	-	100% (100%, 100%)
<i>KEL</i>	1963	0	-	100% (100%, 100%)
<i>KIT</i>	2012	3	100% (43.85%, 100%)	100% (100%, 100%)
<i>KMT2A</i>	14698	1	-	100% (100%, 100%)
<i>KMT2B</i>	11058	0	-	100% (100%, 100%)
<i>KMT2C</i>	19770	4	100% (43.85%, 100%)	100% (100%, 100%)
<i>KMT2D</i>	20922	3	-	100% (100%, 100%)
<i>KRAS</i>	929	0	-	100% (100%, 100%)
<i>LATS1</i>	4133	0	-	100% (100%, 100%)
<i>LATS2</i>	4047	0	-	100% (100%, 100%)
<i>MAP2K1</i>	632	0	-	100% (100%, 100%)
<i>MAP3K1</i>	6425	1	-	100% (100%, 100%)
<i>MAPK1</i>	1000	0	-	100% (100%, 100%)
<i>MAX</i>	1227	0	-	100% (100%, 100%)
<i>MEN1</i>	2606	0	-	100% (100%, 100%)
<i>MET</i>	1504	0	-	100% (100%, 100%)
<i>MLH1</i>	5020	0	-	100% (100%, 100%)
<i>MLH3</i>	5599	1	-	100% (100%, 100%)
<i>MPL</i>	946	0	-	100% (100%, 100%)
<i>MRE11</i>	3963	0	-	100% (100%, 100%)
<i>MSH2</i>	4576	0	-	100% (100%, 100%)
<i>MSH3</i>	5529	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>MSH6</i>	5188	3	0% (0.0%, 79.35%)	100% (100%, 100%)
<i>MTAP</i>	1674	0	-	100% (100%, 100%)
<i>MUTYH</i>	3016	0	-	100% (100%, 100%)
<i>MYC</i>	1154	0	-	100% (100%, 100%)
<i>NBN</i>	3558	0	-	100% (100%, 100%)
<i>NCOR1</i>	11641	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>NF1</i>	12882	12	100% (43.85%, 100%)	100% (100%, 100%)
<i>NF2</i>	3582	0	-	100% (100%, 100%)
<i>NFE2L2</i>	936	0	-	100% (100%, 100%)
<i>NOTCH1</i>	11747	8	100% (34.24%, 100%)	100% (100%, 100%)
<i>NOTCH2</i>	10611	1	-	100% (100%, 100%)
<i>NOTCH3</i>	9364	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>NOTCH4</i>	8205	2	-	100% (100%, 100%)
<i>NRAS</i>	829	0	-	100% (100%, 100%)
<i>NTRK1</i>	1442	0	-	100% (100%, 100%)
<i>PALB2</i>	4460	2	50% (9.45%, 90.55%)	100% (100%, 100%)
<i>PARP1</i>	5183	0	-	100% (100%, 100%)
<i>PBRM1</i>	7314	1	-	100% (100%, 100%)

Gene	Total Bases	Total Number of Variant per Gene	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
<i>PHF6</i>	1804	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>PIK3CA</i>	1750	3	100% (34.24%, 100%)	100% (100%, 100%)
<i>PIK3R1</i>	3282	4	100% (51.01%, 100%)	100% (100%, 100%)
<i>PIK3R2</i>	1441	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>PMS1</i>	2271	1	0% (0.0%, 79.35%)	100% (100%, 100%)
<i>PMS2</i>	4100	0	-	100% (100%, 100%)
<i>POLD1</i>	5875	0	-	100% (100%, 100%)
<i>POLE</i>	11804	1	-	100% (100%, 100%)
<i>POT1</i>	3363	2	-	100% (100%, 100%)
<i>PPM1D</i>	2382	0	-	100% (100%, 100%)
<i>PPP2R1A</i>	1277	0	-	100% (100%, 100%)
<i>PPP6C</i>	974	0	-	100% (100%, 100%)
<i>PRDM1</i>	3323	0	-	100% (100%, 100%)
<i>PTCH1</i>	7082	2	-	100% (100%, 100%)
<i>PTEN</i>	2594	3	100% (43.85%, 100%)	100% (100%, 100%)
<i>PTPN11</i>	940	0	-	100% (100%, 100%)
<i>RAD50</i>	6073	3	100% (20.65%, 100%)	100% (100%, 100%)
<i>RAD51B</i>	2324	0	-	100% (100%, 100%)
<i>RAD51C</i>	2038	0	-	100% (100%, 100%)
<i>RAD51D</i>	1997	0	-	100% (100%, 100%)
<i>RAF1</i>	1120	0	-	100% (100%, 100%)
<i>RASA1</i>	5763	0	-	100% (100%, 100%)
<i>RBI</i>	5388	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>RBM10</i>	5150	0	-	100% (100%, 100%)
<i>RECQL4</i>	2856	0	-	100% (100%, 100%)
<i>RET</i>	1424	0	-	100% (100%, 100%)
<i>RHOA</i>	162	0	-	100% (99.98%, 100%)
<i>RNF43</i>	3231	2	-	100% (100%, 100%)
<i>SDHA</i>	3393	0	-	100% (100%, 100%)
<i>SDHC</i>	981	0	-	100% (100%, 100%)
<i>SETD2</i>	9510	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>SF3B1</i>	1083	0	-	100% (100%, 100%)
<i>SLX4</i>	6660	2	-	100% (100%, 100%)
<i>SMAD4</i>	2846	0	-	100% (100%, 100%)
<i>SMARCA4</i>	8148	1	100% (20.65%, 100%)	100% (100%, 100%)
<i>SMARCB1</i>	2129	0	-	100% (100%, 100%)
<i>SOX9</i>	1859	4	100% (20.65%, 100%)	100% (100%, 100%)
<i>SPEN</i>	12166	2	-	100% (100%, 100%)
<i>SPOP</i>	916	0	-	100% (100%, 100%)
<i>STAG2</i>	6507	2	-	100% (100%, 100%)
<i>STAT3</i>	1023	0	-	100% (100%, 100%)
<i>STK11</i>	2331	2	100% (34.24%, 100%)	100% (100%, 100%)
<i>TBX3</i>	2732	0	-	100% (100%, 100%)

<b>Gene</b>	<b>Total Bases</b>	<b>Total Number of Variant per Gene</b>	<b>PPA (2-sided 95% CI)</b>	<b>NPA (2-sided 95% CI)</b>
<i>TCF7L2</i>	3618	0	-	100% (100%, 100%)
<i>TERT</i>	1206	1	-	100% (100%, 100%)
<i>TET2</i>	7164	2	100% (20.65%, 100%)	100% (100%, 100%)
<i>TGFBR2</i>	2501	0	-	100% (100%, 100%)
<i>TNFAIP3</i>	2903	0	-	100% (100%, 100%)
<i>TNFRSF14</i>	1637	0	-	100% (100%, 100%)
<i>TP53</i>	2344	11	100% (70.09%, 100%)	100% (100%, 100%)
<i>TP63</i>	1397	0	-	100% (100%, 100%)
<i>TSC1</i>	5349	0	-	100% (100%, 100%)
<i>TSC2</i>	10092	0	-	100% (100%, 100%)
<i>U2AF1</i>	1008	0	-	100% (100%, 100%)
<i>UGT1A1</i>	146	0	-	100% (99.98%, 100%)
<i>VHL</i>	889	1	-	100% (100%, 100%)
<i>WT1</i>	2792	0	-	100% (100%, 100%)
<i>XRCC2</i>	1157	0	-	100% (100%, 100%)
<i>ZFHX3</i>	10946	2	-	100% (100%, 100%)