



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

DEC 14 2001

Mr. Jacob Møllenbach
R & D Manager
K-Systems, Kivex Biotec Ltd.
Klintehøj Vænge 21
DK-3460 Birkerød
DENMARK

Re: K013401
Trade/Device Name: K-Systems Culture Bag Sealer
and Culture Bags
Regulation Number: 21 CFR 884.6120
Regulation Name: Assisted reproduction accessories
Regulatory Class: II
Product Code: 85 MQG
Dated: October 12, 2001
Received: October 15, 2001

Dear Mr. Møllenbach:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.


This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon
Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

K013401

Statement of Indications for Use


The K-Systems Culture bag system is intended for use during Assisted Reproduction Technology (ART).

The K-SYSTEMS Culture Bag System is a device for storing culture dishes (with embryos and gametes) and to provide a stable CO₂ and pH microenvironment.

The K-Systems Culture Bag Sealer is a bench top device which inflates the K-Systems Culture Bags with correct gas mixture and which seals the K-Systems Culture Bags. The K-Systems Culture Bags are specially designed Culture Bags which have low penetration value for Oxygen and CO₂ and which are flexible and transparent making it possible to evaluate the cells within the culture bag.

When the K-Systems Culture Bags have been filled with the appropriate gas mixture the K-Systems Culture Bags are placed in a conventional incubator for maintaining the temperature.

The K-Systems culture bag system is a closed system that reduces the risk of cross contamination within the incubator.



(Division Sign-Off)

Division of Reproductive, Abdominal,
and Radiological Devices

510(k) Number

K013401

Prescription Use 

Please note new address and phone/fax number

K-SYSTEMS

KIVEX Blotac Ltd. Klintehøj Vænge 21

DK-3480 Birkerød, Denmark

Tel: +45 4599 5600

Fax: +45 4599 5619