



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Abraham Lavi, Ph.D., MBA
President
Vilex, Inc.
345 Old Curry Hollow Road
Pittsburgh, PA 15236

MAR 18 2002

Re: K014154

Trade/Device Name: Vilex/Duval/Orthex Cannulated Bone Screw Double Thread
Regulation Number: 21 CFR 888.3040
Regulation Name: Smooth or threaded metallic bone fixation fastener
Regulatory Class: Class II
Product Code: HWC
Dated: December 1, 2001
Received: December 18, 2001

Dear Dr. Lavi:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Dr. Abraham Lavi

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

Miriam C. Provost
for Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure



Inc.

Distributors of Orthex™, DuVal™, and Vilex™
orthopedic implant products, UAM™
Thomas Rose™ surgical instruments.

Phone: (412) 384-7570
FAX: (412) 384-7560
Orders: (800) 872-4911

Corporate Offices
345 Old Curry Hollow Road
Pittsburgh, Pennsylvania 15236 US/
e-mail: vilex@usaor.net

K991151, K991197

page 1 of 1

510 (K) NUMBER (IF KNOWN): K014154

DEVICE NAME: MULTITHREAD VILEX BONE SCREW

INDICATIONS FOR USE:

The Multithread Vilex Cannulated Bone Screw, as designed, has the following Indications for Use: Bone Fractures, Osteotomies, Arthrodeses, and Tendon Attachment. It is intended for, but not limited to, Hand Surgery, Orthopedic and Podiatric Surgery. The materials used to manufacture this screw are 316L, implant-quality stainless steel and Ti6Al4V, implant quality titanium alloy.

(Revised 3/5/02)

Prescription Use X
(Per 21 CFR 801.109)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Miriam C. Provost
(Division Sign-Off)
Division of General, Restorative
and Neurological Devices

510(k) Number K014154