

K 955000

510(k) SUMMARY

23.02.96

APR 28 1997

Datex Division  
Instrumentarium Corp.  
P.O. Box 446  
FIN-00101 Helsinki  
Finland  
Tel: +358 0 39411  
Fax: +358 0 1463310

Contact person: Rauno Ruoho

### PRODUCT NAMES

Device name: AS/3™ M-NMT

Common name: Neuromuscular transmission Device

### DESCRIPTION OF THE DEVICE

The AS/3™ NMT Module is a part of the modular anesthesia monitoring system AS/3™. It provides relaxation measurement of patient.

The product is designed to be used for monitoring relaxation of patient during and after anesthesia and for regional block stimulation for nerve location. The intended environment of use is anesthesia department, including operation rooms (OR) and post-anesthesia care units (PACU) etc.

The NMT measurement hardware is in the NMT module. The NMT Module stimulates the patient controlled by the AS/3™ Anaesthesia or Compact Monitor central unit. Then the NMT Module measures the response of the stimulation and sends it to the AS/3™ Anaesthesia or Compact Monitor central unit. The AS/3™ monitor software receives the responses from the NMT module, calculates the NMT values, displays the NMT values, handles NMT trending and NMT monitor features.