

K955930

A SUMMARY OF SAFETY AND EFFECTIVENESS INFORMATION

NAME OF DEVICE: Beltone Composer 2000 BTE HP

TYPE OF DEVICE: The Beltone Composer 2000 BTE HP Air Conduction Hearing Aid is identical to the Maico Gamma RD BTE HP programmable hearing aid (K933291) with the exception of changes to the labeling.

INTENDED USE: To amplify and transmit sound to the ear

FEATURES: Dispenser programmable
User remote control

ASSEMBLY: Assembled from standard components that are widely used by other hearing aid manufacturers

TECHNICAL CHARACTERISTICS: Technical specifications comply with S3.22-1987 ANSI Standards

FIT: Audiological characteristics dictated by individual audiogram

CONTROLS: Volume control, switches, and optional dispenser operated trimmer controls similar to those of other devices

POWER: Standard hearing aid battery

A USER'S MANUAL AND OTHER INFORMATION
IS SUPPLIED WITH EACH HEARING AID.

