

K960353

510(k) SUMMARY

K960353

January 22, 1996

APR 18 1996

This 510(k) summary of safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990.

1. Submitter:  
CeraMed Corporation  
12860 West Cedar Drive  
Lakewood, CO 80228  
(303) 985-0800
2. Device Name:  
OsteoGraf/LD-300 Hydroxylapatite  
Classification Name: Endosseous implant for bone filling and/or augmentation
3. Predicate Device:  
OsteoGraf/D-700 (previously OsteoGraf/AR) and others
4. Device Description:  
OsteoGraf/LD-300 is a high purity, radiopaque, polycrystalline form of hydroxylapatite, the major mineral phase of bone and dental enamel. It is manufactured as rounded, irregular shaped synthetic hydroxylapatite particles, sized at 250-420 microns.
5. Intended Use:  
The intended use of OsteoGraf/LD-300 is for the filling of periodontal defects and augmentation of bony defects of the alveolar ridge, including tooth extraction sites.
6. Comparison of Product Characteristics:  
OsteoGraf/LD-300 consists of 100% anorganic hydroxylapatite,  $\text{Ca}_{10}(\text{PO}_4)_6\text{OH}_2$ .

X-ray diffraction and infrared analysis (FTIR) show OsteoGraf/LD-300 to be 100% hydroxylapatite. OsteoGraf/LD-300 conforms to the requirements of ASTM standard #F1185, "Composition of Ceramic Hydroxylapatite for Surgical Implants."