

MAY - 7 1996

K961143

510(K) SUMMARY

Trade Name: Telex Select 2-40 with Adaptive Compression®
Common Name: FM at-the-ear Auditory Trainer
Classification Name: Hearing Aid, Group and Auditory Trainer, 77EPF
Equivalent to: Phonic Ear PE700 Free Ear Auditory Trainer, or
AVR/Sonovation Extend Ear, or
Telex Model TDR-7AA K913495

9600 Aldrich Avenue South

Minneapolis, Minnesota 55420 USA

Telephone 612-884-4051

Fax 612-884-0043

Description:

Intended use: to amplify and transmit sound to the ear.
Features: Adaptive Compression; True FM Override (patent pending); two channel switch selectable receiver (patent pending); user replaceable channel crystals (patent pending); function switch allows hearing aid only, FM reception only, or both;
Assembly: assembled from standard components that are widely used by other hearing aid manufacturers.
Technical characteristics: technical specifications comply with S3.2-1987 ANSI standards.
Fit: frequency response per S3.2-1987 ANSI standard as shown on specification filed with 510(k).
Controls: Volume control, Mode Switch, Channel Selector Switch, Audiological Controls (output, gain, and tone).
Power: Standard hearing aid battery (675)

Comparison to predicate device: identical to the Phonic Ear PE700, except for two channel switch selectable FM reception, user replaceable channel crystals, and an additional audiological control (gain).

Submitted by:



Tom Scheller
Chief Engineer
Hearing Instruments Group
Telex Communications, Inc.
9600 Aldrich Ave S.
Minneapolis, MN 55420
(612) 884-4051 voice
(612) 884-0043 fax

Contact: Tom Scheller
Prepared: 19 March, 1996