

**Summary of Safety and Effectiveness Information**

- Name of device** : Phonak Piconet2 P2 AZ
- Type of device** : Programmable behind-the-ear hearing instrument. Substantially equivalent to other programmable behind-the-ear hearing instruments
- Intended Use** : To amplify and transmit sound to the ear
- Features** : 3 programmable memories, telecoil, 5 forms of output limiting, remote control
- Assembly** : Assembled from components available to hearing instrument manufacturer. Delivered completely assembled to the dispenser
- Tech. Characteristics** : Technical specifications comply with ASA/ANSI S3.2 - 1987
- Fit** : Frequency response, gain, output are fit to the individual audiogram
- Controls** : Both manually operated volume control with on/off and/or remote control is available
- Power Source** : Standard hearing instrument battery, size 13

A user's manual and other information is supplied with each hearing instrument.