

Exhibit 8**Summary of Safety and Effectiveness Information**

Name of device : Phonak Sona P2

Type of device : Programmable in-the-canal hearing instrument. Substantially equivalent to other programmable in-the-canal hearing instruments

Intended Use : To amplify and transmit sound to the ear

Features : 3 programmable memories, 5 forms of output limiting, remote control

Assembly : Assembled from components available to hearing instrument manufacturer. Delivered completely assembled to the dispenser

Tech. Characteristics : Technical specifications comply with ASA/ANSI S3.2 - 1987

Fit : Frequency response, gain, output are fit to the individual audiogram

Controls : Operated with remote control

Power Source : Standard hearing instrument battery, size 312 or 10

A user's manual and other information is supplied with each hearing instrument.