

MAY 29 1963

K961663

SUMMARY OF SAFETY AND EFFECTIVENESS

You will find enclosed literature on Five (5) Lift chair companies already marketing very similar products.

As you will see the motors and switches are all UL listed and CSA approved on each of the five different chairs. I will be using the same motor and switch made by Hubbell Corp.

The lift chair market has been in existence for 25 years, so the product has been fully tested through many cycles.

My product is constructed of the same materials, as that of all the lift chair manufacturers I have provided. The chairs consist of oak wood, foam and fabric. The steel frame is constructed of 1"x 1", 16 gauge steel and also 1"x 2" 14 gauge steel.

The weight capacity of all their units is 300 lbs. which will be the same on mine.

The frame bases on all the units are basically the same width for stability. My frame will be the same.

As you can see my products are actually the same as the products of the five manufacturers enclosed.

As a summary of safety, my products do not raise any safety requirements.

Thank you very much and if you have any questions, my telephone number is (717)-822-8629.

My address is as follows:

788 Keating Street
Wilkes Barre Twp, PA. 18702