

JUL 19 1996

K961880

(1) 510(k) Summary:

This summary is submitted in conformance with the format described in the 21 CFR Part 807.92 interim rule dated April 28, 1992.

(1) Submitted By:

MedCare Medical Group, Inc.
234 Old Homestead Highway
E. Swanzey, NH 03446

Contact Person:

Craig J. Bell
Director of Research and Development
(603)352-3230

(2) Name of Device:

Suction Connecting Tubing

Classification Name:

Class II device, reference 21 CFR, Part 880.6740

(3) Identification of Predicate Devices:

Suction Connecting Tubing

(a) Kendall Healthcare Products Company
Mansfield, MA 02048

(4) Description of Device Subject to Premarket Notification:

The MedCare Suction Connecting Tubing is a single use disposable device indicated for general suction procedures as a connecting conduit from a suction device (i.e. Yankauer suction device, suction catheter, etc.) to a suction collection canister (vacuum source). The Suction Connecting Tubing is constructed of polyvinyl chloride (PVC) plastic extruded tubing with injection molded tapered connectors solvent bonded to each end. The connectors and made of PVC and have tapered internal diameters with internal ribs

(5) **Intended Use of the Device:**

The MedCare Suction Connecting Tubing is a single use disposable device indicated for general suction procedures as a connecting conduit from a suction device (i.e. Yankauer suction device, suction catheter, etc.) to a suction collection canister (vacuum source).

(6) **Technological Characteristics**

The MedCare Suction Connecting Tubing is made of similar construction and similar design to the predicate device for, as the name suggests, providing a tubing connection from the suction canister to a suction device. The MedCare Suction Connecting Tubing is provided with tapered fittings to allow secure fitment to a variety of suction devices.