

Karl Storz
Endoscopy-America, Inc.

600 Corporate Pointe
Culver City, California 90230-7600
Phone 310 558 1500

Toll Free 800 421 0837
Fax 310 410 5527

K962090

SEP 30 1996

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of the Safe Medical Devices Act (SMDA) of 1990 and 21 CFR 807.92. All data included in this document are accurate and complete to the best of KSEA's knowledge.

Applicant: Karl Storz Endoscopy - America, Inc.
600 Corporate Pointe
Culver City, CA 90230
(310) 558-1500

Contact: Betty M. Johnson
Manager, Regulatory Affairs

Device Identification: Common Name
Suction/irrigation Pump

Trade Name
Karl Storz Unimat Plus

Indication: The Karl Storz Unimat Plus system is designed to provide irrigation and aspiration of the surgical site during general, endoscopic and laparoscopic surgical procedures.

Device Description: The Karl Storz Unimat Plus and accessories is a membrane type irrigation/aspiration pump with fixed irrigation/aspiration flow rate of 3.5 liters/min. and a maximum irrigation pressure of 400 mmHg.

Substantial Equivalence: The Karl Storz Unimat Plus system is substantially equivalent to the predicate devices since the basic features and intended uses are similar. The minor differences between the Karl Storz Unimat Plus and the predicate devices raise no new issues of safety and effectiveness, as these differences have no effect on the performance, function or intended use of these devices.

Signed: Betty M. Johnson
Betty M. Johnson
Manager, Regulatory Affairs

000060