

AUG 15 1996

**510(K) SUMMARY**  
**(as required by 807.92(c))**

K 962161

**Submitter of 510(k):**

Tamima Itani, Ph.D.  
Baxter Healthcare Corp.  
I.V. Systems Division  
Route 120 & Wilson Rd.  
Round Lake, IL 60073

Regulatory & Marketing Services, Inc. (RMS)  
40178 U.S. 19 North  
Tarpon Springs, FL 34689

Phone: 813-942-3908  
Fax: 813-942-3828

**Contact Person:**

Ed Ransom

**Date of Summary:**

04/26/96

**Trade Name:**

Baxter Nebulizer Cap

**Classification Name:**

Nebulizer Cap

**Predicate Device:**

Baxter Healthcare Corporation K780406

**Device Description/  
Comparison:**

All blueprints dimension and specifications are exactly the same in this device as found in the predicate device. The only change is the aluminum alloy and the coating. The rest of the device remains the same as currently marketed.

**Intended Use:**

The Baxter Healthcare Corporation's Nebulizer Cap is for use with Sterile Water for Inhalation, USP, in 500 ml, 1000 ml, and 1500 ml plastic bottles, for inhalation therapy only.

MAY 30 1996

12