



**Health o meter**

**510(k) Summary  
K962252**

**NOV 25 1996**

Submission Date: September 5, 1996

Contact: Frank Brichetto  
Sr. Product Manager  
(708) 298-9100 x 401 - tel.  
(708) 233-5472 - fax

Classification Name: Powered Heating Pad (21 CFR, 890.5740)  
Common/Usual Name: Electric Heating Pads  
Proprietary Name: Health o meter Electric Heating Pad with massage,  
Model EHP4  
Classification: Class II

Description: Electric Heating Pad, 11.5" x 14.5" conforming to ANSI / UL Standard #130, Electric Heating Pads, operating on DC 9V current supplied through a standard adapter. The units are supplied with removable covers and sponge inserts for moist heat applications, with instructions, in conventional 4-color printed boxes, and are marketed through consumer retail channels.

Intended Use: *Use wherever applications of heat, massage, or both are desirable for personal comfort, and wherever recommended by your physician, dentist or health care professional for the relief of minor muscular or joint pain.*

Substantial Equivalence Claim: *The Model EHP4 is substantially equivalent to the following models which are currently in distribution:*

- |    |  |         |
|----|--|---------|
| 1. | Conair Moist Heat Pad,                         | K924799 |
| 2. | Conair Moist Heat Pad with Automatic Shut-off, | K924800 |
| 3. | Sunbeam / Walgreens Heating Pad                | none    |

Substantial equivalence is claimed because: manufacturing technology, operating principles, intended uses and safety standards are the same for all models. Functional differences, i.e. massage action, do not adversely affect safety or efficacy. Other differences are cosmetic in nature.

Revised September 4, 1996