

K962378

AUG - 2 1996

510(k) SUMMARY
FOR
OLYMPUS FG-45/46-1 GRASPING FORCEPS

Device Name: Olympus FG-45/46-1 Grasping Forceps

Common/Usual Name: Grasping Forceps

Classification Name: Bronchoscope and Accessories
21 CFR 876.4680, Class II

Predicate Device: Olympus FG-2 Reusable Grasping Forceps

Contact Person: Barry Sands
Olympus America, Inc.
Endoscope Division
2 Corporate Center Drive
Melville, NY 11747-3157
(516) 844-5747

Summary Preparation Date: June 13, 1996

Statement of Intended Use: The Olympus FG-45/46-1 Grasping Forceps is intended to be used to retrieve foreign bodies and excised tissue specimens endoscopically.

Device Description: The Olympus FG-45/46-1 Grasping Forceps is specifically designed for foreign body and tissue specimen collection under endoscopic visualization. The FG-45/46-1 consists of a flexible shaft and a proximal control handle. Operation of the proximal control handle actuates the distal tip grasping prong mechanism.