

K962533

AUG - 9 1996

**510(k) SUMMARY  
FOR  
OLYMPUS FG GRASPING FORCEPS**

**Device Name:** Olympus FG-4/14/20/21/22/26/32/38/40/##-1 Grasping Forceps (Urology)

**Common/Usual Name:** Grasping Forceps for bronchoscopy

**Classification Name:** Bronchoscope and Accessories  
21 CFR 874.4680, Class II

**Predicate Device:** Olympus pre-Amendment and 510(k) cleared grasping forceps  
Olympus FG Reusable Grasping Forceps (K955066) for GI

**Contact Person:** Barry Sands  
Olympus America, Inc.  
Endoscope Division  
2 Corporate Center Drive  
Melville, NY 11747-3157  
(516) 844-5474

**Summary Preparation Date:** June 26, 1996

**Statement of Intended Use:** The Olympus FG series of grasping forceps are intended to be used to grasp tissue and/or retrieve foreign body and excised tissue endoscopically.

**Device Description:** The Olympus FG series of grasping forceps are designed for grasping tissue and/or retrieving foreign body and excised tissue under endoscopic visualization. These forceps consist of a flexible shaft and a proximal control handle. Operation of the proximal control handle actuates the distal tip grasping jaws.