

SEP 27 1986

K963519

SAFETY AND EFFECTIVENESS SUMMARY - Euroline A•675 Series

NAME OF DEVICES: A•675 TSP, A•675 Sequel

TYPE OF DEVICE: Air Conduction Behind-the-Ear hearing aids -Substantially equivalent to the Starkey A-312 Series BTE Hearing Aids (K960494)

INTENDED USE: To assist in compensating for hearing loss by amplifying sound as it comes near the ear .

FEATURES: Adjustable volume control, Microphone/Telephone/Off switch, direct audio input.

ASSEMBLY: Assembled from standard components that are widely used in hearing aid devices.

TECHNICAL CHARACTERISTICS: Technical specifications comply with ANSI S3.22-1987 and IEC 118-0 1983.

FIT: For conductive and/or sensori-neural hearing loss where BTE type hearing aids are preferred or required.

CONTROLS: Volume control, Mic/Telephone/Off switch similar to those used on other hearing aid devices. Potentiometers and are incorporated as circuit options require. Any pots or switches are similar to those used in other BTE hearing aid devices.

POWER: Standard 675 or equivalent hearing aid battery.

8/26/96

SUMMARY -- Starkey Euroline A 675 Series SUBSTANTIAL EQUIVALENCE

The Starkey Euroline A•675 Series BTE Hearing aids are behind-the-ear air conduction hearing aids. This form of hearing aid is intended to be received by a Hearing Health Professional in a completed form, not being intended for final assembly.

The Starkey Euroline A•675 Series BTE Hearing aids are substantially equivalent to the Starkey Laboratories Euroline A-13 Series BTE Hearing aids.

Euroline A 675 Series

- BTE Air conduction
- Adjustable Gain Control
- Telephone/Mic/Off Switch
- Adjustable potentiometers for low frequency, peak clipping, gain, high frequency attenuation, high or low freq compression, and crossover frequency
- 675 or equiv Battery
- High output, High gain
- Colors: Beige or Dark brown

Euroline A-13 Series

- BTE Air conduction
- Adjustable Gain Control
- Telephone/Mic/Off Switch
- Adjustable potentiometers for low frequency, peak clipping or AGC, and high frequency attenuation
- 13 or equiv Battery
- Miniature size
- Colors: Choice of 15 different skin ton

6/26/96