

K964451

NOV 27 1996

**A SUMMARY OF SAFETY AND EFFECTIVENESS INFORMATION**

**Name of Device:** Piccolo K-amp VC+, Piccolo P-A1, Piccolo ATC and Piccolo P-AO

**Type of Device:** BTE Aid – Substantially equivalent to other BTE aids.

**Intended Use:** To amplify and transmit sound to the ear.

**Features:** The Piccolo series of BTE's are a small compact hearing aid. It has incorporated in it's package a very slim housing that can be worn with glasses. The use of the new modern and small receiver enables the integrated amplifier to have better battery life with less distortion. The aid also has AGC-I, active low frequency tone control and a three position, rotating switch for Off-Telecoil-Microphone.

**Assembly:** Assembled from standard components which are widely utilized by other hearing aid manufacturers.

**Technical Characteristics:** Technical specifications comply with S3.22-1987 ANSI Specifications. Preliminary data sheets for the Piccolo K-amp VC +, Piccolo P-A1, Piccolo ATC and the Piccolo P-AO are enclosed.

**Fit:** The frequency response of this product is dictated by the individual audiogram from each client.

**Power:** Standard hearing aid battery: size 13

A user's manual and other information is supplied with each hearing aid (enclosed).