

NOV 27 1996

A Summary of Safety and Effectiveness Information

Name of Device: Unitron Industries, Sound F/XD ⁺4

Type of Device: Behind-The-Ear Hearing Aid. Substantially equivalent to other behind-the-ear hearing aids.

Intended Use: To amplify and transmit sound via air conduction to the ear.

Features: Three fitter controls - Cross over frequency, low channel gain, and threshold knee point.

Two channel compression circuitry.

Class D circuitry for excellent battery life.

Telecoil with pre-amplifier.

Directional Microphone.

Assembly: Assembled from standard components that are used by Unitron and other hearing aid manufacturers.

Technical Characteristics: Technical specifications comply with ANSI Standard S 3.22-1987.

Fit: Hearing aid adjustment dictated by individual audiogram.

Controls: Volume control, similar to other devices. Switch for selecting microphone, telecoil and off.

Power: Standard hearing aid battery - size 13.

A user's guide and other information is supplied with each hearing aid.