

K 973384

## SECTION 7

**SUMMARY OF SAFETY AND EFFECTIVENESS**

---

**510(k) Summary of  
Safety and Effectiveness**

Information supporting claims of substantial equivalence, as defined under the Federal Food, Drug and Cosmetic Act, respecting safety and effectiveness is summarized below. For the convenience of the Reviewer, this summary is formatted in accordance with the Agency's final rule "...510(k) Summaries and 510(k) Statements..." (21 CFR 807) and can be used to provide a substantial equivalence summary to anyone requesting it from the Agency.

NEW DEVICE NAME: ETHICON Non-Stick Bipolar Forceps

PREDICATE DEVICE NAME: Kirwan Surgical Products  
Bipolar Forceps

---

**510(K) SUMMARY**

---

**Device Description**

The ETHICON Non-Stick Bipolar forceps are available in various standard sizes and shapes similar to conventional forceps. These forceps can be connected to the bipolar output mode on electrosurgical generators to facilitate bipolar coagulation and bipolar point coagulation.

The ETHICON Non-Stick Bipolar Forceps are designed for use with the ETHICON Bipolar Cable (K#960476). The ETHICON Bipolar Cable is used to connect the ETHICON Non-Stick Bipolar Forceps to the generator (ESU).

---

**Intended Use**

The intended use of the ETHICON Non-Stick Bipolar Forceps is to facilitate tissue grasping/manipulation and bipolar coagulation.

---

Continued on next page

## SUMMARY OF SAFETY AND EFFECTIVENESS, Continued

---

**Indications  
Statement**

The ETHICON Non-Stick Bipolar Forceps are non-sterile, reusable devices intended to facilitate tissue grasping/manipulating bipolar point coagulation and bipolar coagulation while grasping soft tissue, including vessels up to 3mm in diameter in open surgical procedures.

---

**Technological  
Characteristics**

The new device is technologically the same as the predicate.

---

**Performance Data**

Preclinical laboratory evaluations (complies with IEC-602-2.2) were performed to ensure that the device functions as intended. Clinical data was deemed unnecessary to support the Premarket Notification. Sufficient data has been gathered from preclinical testing to assess the safety and effectiveness characteristics of the new device.

---

**Conclusions**

Based on the 510(k) summaries and 510(k) statements (21 CFR 807) and the information provided herein, we conclude that the new device is substantially equivalent to the Predicate Device under the Federal Food, Drug and Cosmetic Act.

---

**Contact**

Gregory Jones  
Director, Regulatory Affairs  
ETHICON, Inc.  
Rt. #22 West  
Somerville, NJ 08876-0151

---

**Date**

September 5, 1997

---



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

DEC - 5 1997

Mr. Gregory R. Jones  
Director, Regulatory Affairs  
Ethicon, Incorporated  
P.O. Box 151  
Somerville, New Jersey 08876-0151

Re: K973384  
Trade Name: Ethicon Non-Stick Bipolar Forceps  
Regulatory Class: II  
Product Code: GEI  
Dated: September 5, 1997  
Received: September 8, 1997

Dear Mr. Jones:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for

devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4595. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

  
fr Celia M. Witten, Ph.D., M.D.

Director  
Division of General and  
Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

INDICATION FOR USE

510(k) Number (if known): K973384

Device Name: ETHICON Non-Stick Bipolar Forceps

Indications for Use: The ETHICON Non-Stick Bipolar Forceps are non-sterile, reusable devices intended to facilitate tissue grasping/manipulation, bipolar point coagulation, and bipolar coagulation while grasping soft tissue, including vessels up to 3mm in diameter in open surgical procedures.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON  
ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒ OR Over-The Counter Use ☐  
(Per 21 CFR 801.109)

  
(Division Sign-Off)  
Division of General Restorative Devices K973384  
510(k) Number \_\_\_\_\_ (Optional Format 1-2-9G)