



**510(k) SUMMARY  
AS REQUIRED BY 21 CFR 807.92**

K980601

1. Submitter: Varian Oncology Systems  
3045 Hanover Street  
Palo Alto, CA 94304  
  
Contact: Linda S. Nash, Manager  
Regulatory Compliance & Radiation Safety  
Phone (650) 424-6990  
FAX (650) 424-4830  
linda.nash@os.varian.com  
  
Prepared: February 12, 1998  
Revised: July 23, 1998
2. Device Name: VariSource Wright Vaginal Cuff Applicator for Varian  
VariSource™ Remote High Dose Rate Afterloader.
3. Predicate Device: Mick Radio-Nuclear Instruments, Inc., Wang Front Loading  
Applicator, K890485.
4. Description: Applicators for the Varian VariSource Remote High Dose Rate  
Afterloader are a part of a remote controlled radionuclide  
applicator system, including an electromechanical device to  
enable an operator to apply, by remote control, a radionuclide  
source of high activity at various internal or surface body  
locations for radiation brachytherapy. The shape and materials  
of the applicator determine where it will be utilized for  
treatment.
5. Intended Use: The Varian VariSource Remote High Dose Rate  
Afterloader [system, including applicators and  
accessories] is a device intended to be used by  
properly trained and licensed medical personnel to  
provide radiation brachytherapy. The VariSource  
Wright Vaginal Cuff Applicator which is the subject of  
this 510(k) is a component of the VariSource system.
6. Technological Characteristics: See attached comparison chart.

## Comparison to Predicate Device

#	Feature	MRNI Wang Front Loading Applicator, K89045	VariSource Wright Vaginal Cuff Applicator K980601
1	Afterloading Method	Manual	Remote HDR
2	Coupling Catheter Fittings	No	Yes
3	Vaginal Cylinder		
	Diameter and Length	3.0 cm X 7.0cm	3.4 cm X 12.0 cm
	Material	Polysulfone	Solid Water inside Polysulfone Shell
4	Irradiation Tubes	2	4
	Material	Stainless steel	Stainless Steel
	Configuration	1 central and 1 lateral with source positioned perpendicular to central.	4 central at 0, 90, 180, and 270 degrees. Ends come together to form dome (eggbeater) shape.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

OCT 21 1998

Linda Nash  
Regulatory Compliance  
and Radiation Safety Manager  
Varian Associates, Inc.  
3045 Hanover Street  
Palo Alto, CA 94304

Re: K980601  
VariSource Wright Cuff Applicator  
for VariSource HDR Afterloader  
Dated: June 24, 1998  
Received: August 4, 1998  
Regulatory class: II  
21 CFR 892.5700/Procode: 90 JAQ

Dear Ms. Nash:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4613. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

Lillian Yin, Ph.D.  
Director, Division of Reproductive,  
Abdominal, Ear, Nose and Throat  
and Radiological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure



**STATEMENT of INDICATIONS for USE\***

I state in my capacity as Manager, Regulatory Compliance and Safety, of Varian Oncology Systems that the Product which is the subject of this premarket notification, is intended to be used for the following:

The Varian VariSource™ Remote High Dose Rate Afterloader [system, including applicators and accessories] is a device intended to be used by properly trained and licensed medical personnel to provide radiation brachytherapy. The VariSource Wright Vaginal Cuff applicator which is the subject of this 510(k) is a component of the VariSource system.

A handwritten signature in cursive script, appearing to read "Charles H. Will", written over a horizontal line.

Charles H. Will, Manager  
Regulatory Compliance & Safety

February 12, 1998

Date

\*Suggested language and format to meet the requirements of section 513(i) of the Federal Food, Drug, and Cosmetic Act, as amended, and 21 CFR sections 801.4 and 809.92(a)(5).

A handwritten number "K980601" written in a simple, blocky style.

510(k) Number

A handwritten signature in cursive script, appearing to read "William Yip", written over a horizontal line.

Division Sign-off  
Office of Device Evaluation

Prescription Use ☒  
(Per 21 CFR 801.109)

Over-the-Counter Use ☐