Lovytex Sdn. Bhd.

(Company No. 52031-D) Lot 8, Jalan Suasa, 42500 Telok Panglima Garang, Kuala Langat, Selangor Darul Ehsan, Malaysia. Tel: 603-3526711 Fax: 603-3527733



MAR 13 1998

510 (k) SUMMARY

1. Submission Applicant

| Name: Lovytex Sdn Bhd | |
|--|------------------------------------|
| Street Address: Lot 8, Jalan Suasa ,4250 | 0 Telok Panglima Garang, |
| Selangor Darul Ehssn | _ Country <u>Malaysia</u> |
| Phone No <u>+ 6-03-352 7402</u> | FAX No <u>+ 6-03-352 7733</u> |
| Contact Person: M.G. Desai | |
| | |
| 8040475 | A 444401 |
| [Registration Number, Form 2891(a)] | (Device Listing Number, Form 2892) |
| | |
| Activity: | |
| [X] Manufacturer | |
| | |
| Applicant Lovytex Sdn Bhd | |
| 510(k) Number (if known): K98 | 0618 * |

2. Device Particulars

Device Name: Powder-free Latex Examination Gloves

Trade/Proprietary Name: N/A

Common Name: Powder-free Examination Gloves

Classification Name: Patient Examination Gloves

3. Device Classification

Device Class: Class 1

Product Code: Latex (Powder-free). 80LYY

4. Device Description

Device Class: Class 1 Product Code: Latex (Powder-free). 80LYY

5. Summary of Intended Use

A medical glove is worn on the hand of healthcare and similar personnel to prevent contamination between heathcare personnel and the patient.

6. Summary of Technological Characteristics

| Characteristic I | Reference Document | Device Performance |
|---------------------------|--------------------|-------------------------|
| | | |
| Water Leak | ASTM D 3578/FDA | Meets or Exceeds |
| Residual Powder | ASTM DRAFT/FDA | Meets or Exceeds |
| Tensile | | |
| - Unaged | ASTM D 3578-95 | Meets or Exceeds |
| - Aged | ASTM D 3578-95 | Meets or Exceeds |
| Elongation @ Break | | |
| - Unaged | ASTM D 3578-95 | Meets or Exceeds |
| - Aged | ASTM D 3578-95 | Meets or Exceeds |
| Water Extractable Protein | ASTM D 5712-95 | 50 μ grams per gram |

7. Assessment of Non-Clinical Performance Data

The device :-

- meets or ecxeeds the ASTM standard or equivalent standard
- meets FDA pinhole requirements; and
- meets the labelling claim as shown by the data in Section 6.

8. Assessment of Biocompatibility Performance Data

SUMMARY OF RESULTS OF BIOCOMPATIBILITY TESTS

CYTOTOXICITY

-MEM test extract was NON-TOXIC at dilution 1:8 at 24 hours..

RABBIT SKIN IRRITATION

- PASSES

GUINEA SENSITIZATION STUDY

- PASSES

9. CONCLUSIONS OF NON-CLINICAL & BIOCOMPATIBILITY PERFORMANCE DATA

The device has been carefully compared to existing performance criteria of the ASTM and FDA to which all legally marketed devices are required to be compared. The data summaries indicate that the proposed device meets or exceeds all acceptable scores and satisfies all requirements for powder-free natural latex examination gloves.

Pursuant to 21 CFR 807.87 (j), I M.G. Desai, the Chief Executive Officer of Lovytex Sdn Bhd, certify that to the best of my knowledge and belief and based upon the data and information submitted to me in the course of my responsibilities as the Chief Executive Officer of Lovytex Sdn Bhd, and in reliance thereupon, the data and information submitted in this premarket notification are truthful and accurate and that no facts material to a review of the substantial equivalence of this device have been knowingly omitted from this submission.

(Signature)

M. G. Desai

(Typed Name)

2 - 16 - 98

(Dated)



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

MAR | 3 1998

Mr. M. G. Desai
Chief Executive Officer
Lovytex Sdn. Bhd.
Lot #8 Jalan Suasa
Telok Panglima Garang
Selangor Darul Ehsan
Malaysia

Re: K980618

Trade Name: Lovytex Sdn Bdh Powder-Free Latex

Examination Gloves with Protein Content Labeling Claim

(50 microgram per gram)

Regulatory Class: Product Code: LYY

Dated: February 16, 1998 Received: February 18, 1998

Dear Mr. Desai:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class_III_(Premarket_Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements

concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4618. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "http://www.fda.gov/cdrh/dsmamain.html".

Timo X A. Ulatowski

Director

Division of Dental, Infection Control and General Hospital Devices Office of Device Evaluation Center for Devices and

Radiological Health

Enclosure

3.0 Indication For Use

INDICATION FOR USE

| Applicant : Lovytex Sdn Bhd 510(k) Number (if known): | 980618 | * ; |
|--|--|----------------------------|
| Device Name: Powder-free Latex Exa | umination Gloves 1 (50 mcrogrem pu | glime) |
| Indication For Use: | | |
| A medical glove is worn on the hand o | of healthcare and si | milar personnel to prevent |
| contamination between heathcare pe | rsonnel and the pat | <u>ient.</u> |
| | | |
| • | | |
| | | |
| (PLEASE DO NOT WRITE BELOW THIS | LINE - CONTINUE | ON ANOTHER PAGE IF |
| NEEDED) | | |
| Concurrence of CDRH | Office of Device Ev | aluation (ODE) |
| Prescription Use | OR | Over-The-Counter X |
| Per 21 CFR 801.109 | | |
| | | (optional Format 1-2-96) |
| (Division Sign Officers) | Chin (- | Lin |
| and General Hospital De | (Division Sign-Off) Division of Dental, In | |
| 510(k) Number | and General Hospital 510(k) Number | Devices 98 0 6 1 8 |

CONFIDENTIAL- DO NOT COPY