

JUN 25 1998

K981275



Attachment VI:

**Summary of Safety and Effectiveness Information
[510(k) Summary]**

SUBMITTER

Synthes (USA)
1690 Russell Road
Paoli, PA 19301
(610) 647-9700

Contact: Sheri L. Musgnung

**COMMON OR USUAL
NAME**

Plate, Bone
Screw, Fixation, Intraosseous

**DEVICE
CLASSIFICATION:**

Class II, 21 CFR 872.4760 and 872.4880

PREDICATE DEVICE:

Synthes Mini Lengthening Apparatus (MLA)

DESCRIPTION:

The Synthes Ratcheting Distractor is an intra-oral distraction device. It features two telescoping components, with contoured legs that have threaded screw holes, that are fixed to bone via 1.3 mm – 3.0 mm screws. The distractor is available in right and left versions; has a ratcheting mechanism that is activated by an "activation instrument"; and is capable of distraction lengths of up to 30 mm. It is manufactured from a titanium alloy.

INTENDED USE:

Synthes Ratcheting Distractor is intended for craniofacial surgery, reconstructive procedures, and selective orthognathic surgery of the mandible and maxilla. Specifically, it is indicated for correcting deficiencies in craniofacial bones, especially orthognathic mandible lengthenings.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUN 25 1998

Ms. Sheri L. Musgnung
Regulatory Affairs Associate
Synthes® (USA)
1690 Russell Road
Post Office Box 1766
Paoli, Pennsylvania 19301

Re: K981275
Trade Name: Synthes® (USA) Ratcheting Distractor
Regulatory Class: II
Product Code: MQN
Dated: April 7, 1998
Received: April 8, 1998

Dear Ms. Musgnung:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of

the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control,
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure



K981275

2.0 Indications for Use Statement

Page 1 of 1

510(k) Number (if known): _____

Device Name: Synthes (USA) Ratcheting Distractor

Indications For Use:

Synthes Ratcheting Distractor is intended for craniofacial surgery, reconstructive procedures, and selective orthognathic surgery of the mandible and maxilla. Specifically, it is indicated for correcting deficiencies in craniofacial bones, especially orthognathic mandible lengthenings.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓ OR Over-The-Counter Use _____
(Per 21 CFR 801.109)

Gerald Shuppr
(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number K981275