

MAY 29 1998

K 981389

10. SMDA Summary of Safety and Effectiveness - "510(k) Summary"

A. Submitter Information

Microline, Inc.  
1100 Cummings Center  
Suite 350-G  
Beverly, MA 01915

Telephone: (978) 922-9810  
Contact Person: Mr. Hughes de Laforcade  
President

Date Prepared: April 13, 1997

B. Device Identification

Common/Usual Name: Manual Detachable Surgical Instruments  
Proprietary Name: Disposable Re-New Forceps Laparoscopic Surgical Tips

C. Identification of Predicate Device(s)

The Disposable Re-New Forceps Laparoscopic Surgical Tips are substantially equivalent to their predicate, Re-New Forceps (K974066), previously cleared and currently marketed.

D. Device Description

The Disposable Re-New Forceps Laparoscopic Surgical Tips are an extension of the Re-New Forceps Laparoscopic Instruments product line used to cut, grasp, and dissect various abdominal tissue during Endoscopic (inclusive of laparoscopic) surgical procedures.

The Disposable Re-New Forceps Laparoscopic Surgical Tips consist of a series disposable tip configurations. Tip configurations include: Grasper, Fenestreded Forceps, "Dolphin Nose" Dissector, "Babcock", and "Maryland" Dissector.

E. Substantial Equivalence

The technical characteristics of the Disposable Re-New Forceps Laparoscopic Surgical Tips are identical to those of the Microline Re-New Forceps (K974066). Differences that exist between these devices relating to technical specifications, materials, and physical appearance do not affect the relative safety or effectiveness of the Disposable Re-New Forceps relative to its predicate.

The Disposable Re-New Forceps Laparoscopic Surgical Tips are intended for use to cut, grasp, and dissect various abdominal tissue in endoscopic, including laparoscopic surgical procedures where instruments are introduced into the body through a cannula .

000009



MAY 29 1998

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Microline, Incorporated  
c/o Ms. Jacqueline E. Masse  
Senior Consultant  
Interactive Consulting  
70 Walnut Street  
Wellesley, Massachusetts 02181

Re: K981389  
Trade Name: Disposable Re-New Forceps  
Regulatory Class: II  
Product Code: GEI  
Dated: April 13, 1998  
Received: April 16, 1998

Dear Ms. Masse:


We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4595. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Celia M. Witten, Ph.D., M.D.  
Director  
Division of General and  
Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number (if known): K 981389


Device Name: **Disposable Re-New Forceps Laparoscopic Surgical Tips**

Indications For Use:

**Cutting, Grasping, and Dissecting Various Abdominal Tissue during Endoscopic (inclusive of laparoscopic) Surgical Procedures.**

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)  
Division of General Restorative Devices  
510(k) Number K981389

Prescription Use   
(Per 21 CFR 801.109)

OR

Over-The-Counter Use