

JUL 9 1998

K981706

## 510(k) Summary

**Submitter's Name/Address**

Abbott Laboratories  
1920 Hurd Drive  
Irving, Texas 75038

**Contact Person**

Mark Littlefield  
Section Manager MS 1-8  
ADD Regulatory Affairs  
(972) 518-6062  
Fax (972) 753-3367

**Date of Preparation of this Summary:**

May 13, 1998

**Device Trade or Proprietary Name:**

Amm Cal, Bil Cal, CO<sub>2</sub> Cal,  
HDL Cal, Iron/Mg Cal, LDL Cal,  
MC Cal, UPro Cal, and ISE  
Calibrator Levels 1 and 2.

**Device Common/Usual Name or Classification Name:** Abbott Clinical Chemistry  
Calibrators

**Classification Number/Class:**

Class II

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is: \_\_\_\_\_.

**Intended Use:**

The Ammonia Calibrator is intended for *in vitro* diagnostic use in clinical chemistry assays for Ammonia.

The Bilirubin Calibrator is intended for *in vitro* diagnostic use in clinical chemistry assays for Bilirubin.

The Carbon Dioxide Calibrator is intended for *in vitro* diagnostic use in clinical chemistry assays for Carbon Dioxide.

The HDL Calibrator is intended for *in vitro* diagnostic use in clinical chemistry assays for HDL.

0000010

The Iron/Magnesium Calibrator is intended for *in vitro* diagnostic use in clinical chemistry assays for Iron and Magnesium.

The LDL Calibrator is intended for *in vitro* diagnostic use in clinical chemistry assays for LDL.

The MC Calibrator is intended for *in vitro* diagnostic use in clinical chemistry assays for Albumin, Calcium, Cholesterol, Creatinine, Glucose, Phosphorus, Total Protein, Triglyceride, Urea Nitrogen, and Uric Acid.

The UPro Calibrator is intended for *in vitro* diagnostic use in clinical chemistry assays for protein in urine and CSF.

The ISE Calibrator Levels 1 and 2 are intended for *in vitro* diagnostic use in clinical chemistry assays for Sodium, Potassium, and Chloride.

**Conclusion:** Abbott Clinical Chemistry Calibrators were used for the calibration of each clinical chemistry reagent on the ALCYON 300/300i. The calibration curve generated was used for the quantitation of each analyte for the purpose of collecting performance data in support of the reagent 510(k) as outlined on the next page.

<b>Calibrator Name</b>	<b>Reagent Name</b>	<b>Reagent 510(k) Number</b>
Ammonia Calibrator	Ammonia	K981467
Liquid Bilirubin Calibrator	Total Bilirubin	K981336
	Direct Bilirubin	K981222
Carbon Dioxide Calibrator	Carbon Dioxide	K981231
HDL Calibrator	Direct HDL	K981224
Iron/Magnesium Calibrator	Iron	K981241
	Magnesium	K981192
LDL Calibrator	Direct LDL	K981303
ISE Calibrator Levels I and Level 2	Sodium Potassium Chloride	K974779
Multiconstituent Calibrator	Albumin BCG Albumin BCP Calcium Cholesterol Creatinine Glucose Phosphorus Total Protein Triglyceride Urea Nitrogen Uric Acid	K981468 K981457 K981232 K981476 K981240 K981185 K981118 Exempt Class II K981223 (Exempt Class I) K981123 K981189
Urine/CSF Protein Calibrator	Urine/CSF Protein	K981295



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

JUL 9 1998

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Mark Littlefield  
• Section Manager, Regulatory Affairs  
Abbott Laboratories  
1920 Hurd Drive  
Irving, Texas 75038

Re: K981706  
Amm Cal, Bil Cal, CO<sub>2</sub> Cal, HDL Cal, Iron/Mg Cal, LDL Cal,  
MC Cal, Upro Cal, and ISE Calibrator Levels 1 and 2  
Regulatory Class: II  
Product Code: JIX  
Dated: May 13, 1998  
Received: May 14, 1998


Dear Mr. Littlefield:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market. If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "dsma@fdadr.cdrh.fda.gov".

Sincerely yours,



Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical  
Laboratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health


Enclosure

510(k) Number (if known): \_\_\_\_\_

Device Name: Abbott Clinical Chemistry Calibrators

Indications For Use:

An Abbott Clinical Chemistry Calibrator is a device intended for medical purposes for use in a test system to establish points of reference that are used in the determination of values in the measurement of substances in human specimens.

  
(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number 12981706

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_

(Optional Format 1-2-96)

0000007