

**Ansell  
Edmont  
Industrial**

16982131

**Ansell Edmont Industrial Inc.**

1300 Walnut Street  
Coshocton, OH 43812

Technology Center  
Telephone: (740) 623-3593

Technology Center  
Fax: (740) 623-3515

OCT 29 1998

**Nitra-Touch™**  
**Ansell Edmont Industrial, Inc.**  
1300 Walnut Street  
Coshocton, Ohio 43812  
Telephone: 740-622-4311  
Fax: 740-623-3515

**Checklist  
Section 21.0**

- [1] 510 (k) Summary -
- [2] Ansell Edmont Industrial, Inc.  
1300 Walnut Street  
Coshocton, Ohio 43812
- Telephone: 740-622-4311  
Fax: 740-623-3515
- Contact: Mike W. Hagans  
Telephone: 740-623-3595  
Fax: 740-623-3515
- June 12, 1998
- [3] Trade Name: Nitra-Touch™  
Common Name: Exam Gloves  
Classification Name: Patient Examination Glove
- [4] Nitra-Touch™ examination gloves, meet all of the requirements of ASTM D 3578 with the following variation.
- [5] Nitra-Touch™ examination gloves meet all the current specifications for ASTM D 3578 Rubber Examination Gloves except for the ultimate elongation percentage before and after aging.
- [6] Nitra-Touch™ examination gloves are disposable device intended for medical purposes that is worn on the examiners hand to prevent contamination between patient and examiner and for use handling chemotherapy drugs.
- [7] Nitra-Touch™ examination gloves are summarized with the following technological characteristics compared to ASTM or equivalent standards.

Characteristics	Standard
Dimensions	Meets ASTM D 3578
Physical Properties	Meets ASTM D 3578

K982131

Nitra-Touch™  
Ansell Edmont Industrial, Inc.  
1300 Walnut Street  
Coshocton, Ohio 43812  
Telephone: 614-622-4311  
Fax: 614-623-3515

Except Ultimate Elongation meets  
specifications described under  
Attachment II

Original Ultimate Elongation 500%  
minimum  
Aged Ultimate Elongation 400% minimum

Freedom from holes  
Meets ASTM D 5151

Meets ASTM D 3578

Powder-Free  
Meets described test in Attachment V

Not more than 2 mg residue by mass.

Tensile Strength MPa  
Tensile Strength 60% higher than  
vinyl (synthetic) examination gloves  
Vinyl = 10.5  
Nitra-Touch™ = 17.1

ASTM D 5250  
ASTM D 3578

Puncture lbf/in  
Puncture resistance of Nitra-Touch™  
is eight (8) times greater than vinyl  
and nearly twice that of competitive  
nitrile gloves

ASTM D 120-87  
Nitra-Touch™ = 1625  
Blue Nitrile = 929  
Vinyl (Synthetic) = 184

Glutaraldehyde Resistance  
Nitrile film provides excellent chemical  
resistance to glutaraldehyde.

ASTM T F739-91  
Source: TRI/Environmental, Inc.  
Breakthrough Time = >480

Permeation Rate = <0.013 Fg/cm<sup>2</sup>\* min  
Biocompatibility

Primary Skin Irritation in Rabbits  
Guinea Pig Sensitization

Passes  
Passes

- [8] The performance test data of the non clinical tests are the same as mentioned immediately above.
- [9] Clinical data is not needed for medical gloves or for most devices cleared by the 510(k) process.
- [10] It is concluded that Nitra-Touch™ examination gloves are as safe, as effective, and perform as well as the glove performance standards referenced in Section 7 above and therefore meet:

ASTM listed standards,  
FDA hole requirements, and  
labeling claims for the product.

- [11] This summary will include any other information reasonably deemed necessary by The FDA.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

OCT 29 1998

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Mike W. Hagans  
Ansell Edmont Industrial, Incorporated  
1300 Walnut Street  
Coshocton, Ohio 43812

Re: K982131  
Trade Name: Nitra-Touch™ Powder-Free Nitrile Medical  
Examination Glove  
Regulatory Class: I  
Product Code: LZA  
Dated: September 9, 1998  
Received: September 21, 1998

Dear Mr. Hagans:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Mr. Hagans

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski  
Director

Division of Dental, Infection Control,  
and General Hospital Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**3.0 Indications for Use Statement:**

**INDICATIONS FOR USE**

**Applicant:** Ansell Edmont Industrial, Inc.

**510(K) Number (if known):** K 982131 \*  
**Device Name:** NITRA TOUCH POWDER-FREE NITRILE MEDICAL EXAMINATION GLOVE  
Patient Examination Glove

**Indications For Use:**

A disposable device intended for medical purposes that is worn on the examiners hand to prevent contamination between patient and examiner and for use handling chemotherapy drugs.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH Office of Device Evaluation (ODE)

Chin S. Lin  
(Division Sign-Off)  
Division of Dental, Infection Control,  
and General Hospital Devices  
510(k) Number K 982131

Prescription Use \_\_\_\_\_  
21 CFR 801.109

OR

Over-The-Counter X

(Optional Format 1-2-96)