

JUL 29 1999

PT. BERLIAN GLOVINDO

JLN. SUTOMO UJUNG 96
MEDAN 20235 - INDONESIA

Tel : 62-61-616944
Tel : 62-4-92777
Fax : 62-61-8459180

Page Numbers 1 of 2

"510 (K)" SUMMARY

K 992177

(1) Name of applicant : Mr. Tony Jap
Address : PT. Berlian Glovindo
Jl. Sutomo Ujung No. 96
Medan
Indonesia
Phone No. 62-61-616944 62-4-92777
Fax No. 62-61-8459180

The contact persons within the firm as well as in U.S.A are given below:

Contact person in firm : Mr. Tony Jap
Fax No.: 62-61-8459170

Contact person in U.S.A : Emmy Tjoeng
Fax No.: 626-913-1498

(2) Device details
Trade Name : Private lebel-Nirile Examination Gloves Powder Free
Classification Name : Patient Examination Gloves
Product Code : Nitrile 80 LZA

(3) Equivalent device legally marketed : Class I Nitrile Examination Gloves 80 LZA
Pre-powdered meeting ASTM D 3578-95

(4) Intended use : A patient examination glove is a disposable device intended for medical purpose that is worn on the examiner's hand or finger to prevent contamination between and patient examiners.

PT. BERLIAN GLOVINDO

JLN. SUTOMO UJUNG 96
MEDAN 20235 - INDONESIA

Tel : 62-61-616944
Tel : 62-4-92777
Fax : 62-61-8459180

Page Numbers 2 of 2

(5) Technological characteristic of the gloves.

a. Dimensions

Sizes	XS	S	M	L	XL
Length	240 mm	240 mm	240 mm	240 mm	240 mm
Width	80 < mm	80±10 mm	95±10 mm	111±10 mm	>111mm

Thickness

1. Cuff (min)	0.10 mm				
2. Palm (min)	0.10 mm				
3. Finger Tip (min)	0.10 mm				

b. Physical Properties

	Before aging	After aging at 100°C 22 hrs.
Tensile Strength :	14 Mpa (min.)	14 Mpa (min.)
Ultimate Elongation :	700 % (min.)	500 % (min.)

(6) Performance data is the same as mentioned immediately above.

(7) Clinical data is not needed for gloves or for most devices cleared by the 510 (K) process.

(8) Non-clinical data

We certify that these gloves meet ASTM D 3578 Standard.

Meets FDA pinhole requirement.

Meets labeling claim.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUL 29 1999

Mr. Tony Jap
PT. Berlian Glovindo
Jln. Raya Medan Rantau Prapat
Km +230 Desa Damuli
Kec. Labuhan Batu - North Sumatra
INDONESIA

Re: K992177
Trade Name: Nitrile Examination Gloves - Pre-Powdered,
Small, Medium, and Large
Regulatory Class: I
Product Code: LZA
Dated: June 25, 1999
Received: June 28, 1999

Dear Mr. Jap:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 - Mr. Jap

Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control,
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

PT. BERLIAN GLOVINDO

JLN. SUTOMO UJUNG 96
MEDAN 20235 - INDONESIA

Tel : 62-61-616944
Tel : 62-4-92777
Fax : 62-61-8459180

ANNEXURE II

INDICATION FOR USE

Applicant : Mr. Tony Jap
Device Name : Nitrile Patient Examination Gloves Pre-powdered
Indication for use :
K992177

A patient examination glove is a disposable device intended for medical purpose that is worn on the examiner's hand or finger to prevent contamination between patient and examiners.


(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
E-10(k) Number K 992177