

Encore Orthopaedic Powder Free Surgical Gloves (Protein Labeling Claim)
 Ansell Perry
 1875 Harsh Avenue SE
 Massillon, Ohio 44646
 Telephone: 330-833-2811
 Fax: 330-833-6213

Checklist
 Section 21.0

[1] 510 (k) Summary

[2] Ansell Perry Inc.
 1875 Harsh Avenue SE
 Massillon, Ohio 44646

Telephone: 330-833-2811
 Fax: 330-833-6213

Contact: James R. Chatterton
 Telephone: 330-833-2811
 Fax: 330-833-6213

August 31, 1999

[3] Trade Name: Encore Orthopaedic Powder Free Surgical Gloves (Protein Labeling Claim)
 Common Name: Surgical Gloves
 Classification Name: Surgeon's Glove

[4] Encore Orthopaedic Powder Free Surgical Gloves (Protein Labeling Claim), meet all of the requirements of ASTM D 3577, Type 1.

[5] Encore Orthopaedic Powder Free Surgical Gloves (Protein Labeling Claim) meet all the current specifications for ASTM D 3577 Rubber Surgical Gloves.

[6] Encore Orthopaedic Powder Free Surgical Gloves (Protein Labeling Claim) are sterile disposable devices intended to be worn by operating room personnel to protect a surgical wound from contamination.

[7] Encore Orthopaedic Powder Free Surgical Gloves (Protein Labeling Claim) are summarized with the following technological characteristics compared to ASTM or equivalent standards.

Characteristics	Standard
Dimensions	Meets ASTM D 3577
Physical Properties	Meets ASTM D 3577, Type 1

992953

Encore Orthopaedic Powder Free Surgical Gloves (Protein Labeling Claim)
Ansell Perry
1875 Harsh Avenue SE
Massillon, Ohio 44646
Telephone: 330-833-2811
Fax: 330-833-6213

Freedom from holes	Meets ASTM D 3577 Meets ASTM D 5151
Powder-Free Meets described test in Attachment VI	Meets ASTM D 6124 Not more than 2 mg residue by mass.
Protein Label Claim	This latex glove contains 50 micrograms or less of total water extractable protein per gram.
Biocompatibility Primary Skin Irritation in Rabbits Guinea Pig Sensitization	Passes Passes

- [8] The performance test data of the non clinical tests are the same as mentioned immediately above.
- [9] Clinical data is not needed for medical gloves or for most devices cleared by the 510(k) process.
- [10] It is concluded that Encore Orthopaedic Powder Free Surgical Gloves (Protein Labeling Claim) are as safe, as effective, and perform as well as the glove performance standards referenced in Section 7 above and therefore meet:

ASTM listed standards,
FDA hole requirements, and
labeling claims for the product.

- [11] This summary will include any other information reasonably deemed necessary by The FDA.



OCT 29 1999

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. James R. Chatterton
Vice President Regulatory
Ansell Perry
Ansell Healthcare Products, Incorporated
1875 Harsh Avenue S.E.
Massillon, Ohio 44646

Re: K992953
Trade Name: Encore Orthopaedic Powder Free Surgical
Gloves (Protein Labeling Claim), 50 Micrograms or Less
Regulatory Class: I
Product Code: KGO
Dated: August 31, 1999
Received: September 1, 1999

Dear Mr. Chatterton:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Pre-market Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in

the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

3.0 Indications for Use Statement:

INDICATIONS FOR USE

Applicant: Ansell Perry

510(K) Number (if known): K992953 *

Device Name: Surgeons Glove, brown latex polymer coated, powder free with protein label claim

Indications For Use: (50 MICROGRAMS OR LESS)

A device made of natural rubber intended to be worn by operating room personnel to protect a surgical wound from contamination.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH Office of Device Evaluation (ODE)

Prescription Use _____
Per 21 CFR 801.109

OR

Over-The-Counter X

Chin S. Lim
(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number K992953 3

(Optional Format 1-2-96)