

K111094 #1 | 2

**510(K) SUMMARY**  
ARTHROCare CORPORATION  
SPEEDLOCK KNOTLESS FIXATION SYSTEM

AUG - 9 2011

**General Information**

**Submitter Name/Address:** ArthroCare Corporation  
680 Vaqueros Avenue  
Sunnyvale, CA 94085-3523

**Establishment Registration No.:** 2951580

**Contact Person:** Laura N. Kasperowicz  
Sr. Manager, Regulatory Affairs

**Date Prepared:** April 14, 2011

**Device Description**

**Model Name:** SpeedLock Knotless Fixation Device  
**Generic/Common Name:** Bone Anchor, Fastener, Fixation, Soft Tissue  
**Classification Name:** Fastener, Fixation, Nondegradeable, Soft Tissue  
**Device Classification:** Class II per 21 CFR 888.3040, Product code: MBI

**Model Name:** Drill, 3.0 mm  
**Generic/Common Name:** Bone Drill

**Model Name:** PathFinder® Obturator  
**Generic/Common Name:** Bone Hole Locator

**Model Name:** Sharp-Tipped Obturator  
**Generic/Common Name:** Bone Hole Locator

**Model Name:** Drill Guide, 3.5mm High Visibility  
**Generic/Common Name:** Drill Guide

**Model Name:** Drill Guide, 3.5mm Low Profile  
**Generic/Common Name:** Drill Guide

**Predicate Devices**

Opus® SpeedLock®  
Knotless Fixation Device K090615 (June 3, 2010)

**Product Description**

The SpeedLock Knotless Fixation Device is a bone anchor system with inserter handle designed for specific indications in arthroscopic and orthopedic procedures.

K111044 02/2

## 510(k) SUMMARY

### **Indications For Use**

The SpeedLock Knotless Fixation Device with inserter is indicated for use in fixation of soft tissue to bone. Examples of such procedures include:

**Shoulder:** Bankart Repair, SLAP lesion repair, acromio-clavicular separation, rotator cuff repair, capsule shift/capsule-labral reconstruction, biceps tenodesis and deltoid repair

**Ankle:** Lateral instability, medial instability, Achilles tendon repair/reconstruction and midfoot reconstruction

**Foot:** Hallux valgus reconstruction

**Elbow:** Tennis elbow repair, biceps tendon attachment

**Knee:** Extra-capsular repairs; reattachment of: medial collateral ligament, posterior oblique ligament or joint capsule closure to anterior proximal tibia; extra capsular reconstruction, ITB tenodesis; patellar ligament and tendon avulsions

### **Substantial Equivalence**

By definition, substantial equivalence means that a device has the same intended use and technical characteristics as the predicate device, or has the same intended use and different technological characteristics, but can be demonstrated to be as safe and effective as the predicate device. The SpeedLock Knotless Fixation device design and technology is substantially equivalent to the existing SpeedLock Knotless Fixation Device [K090615]. Side by side comparison bench testing was performed on the proposed and predicate device per the US FDA Guidance Document for Testing Bone Anchors. The differences between the SpeedLock and the predicate device do not raise any questions regarding the safety and effectiveness of the implant. Furthermore, the materials are well characterized and have been used in predicate devices with similar indications. The proposed device, as designed, is as safe and effective as predicate devices.

### **Summary and Reason for 510k Notification**

The purpose of this 510(k) is to notify the Food and Drug Administration of a proposed modification to an existing product. The proposed device, the SpeedLock Knotless Fixation Device is substantially equivalent to the SpeedLock Knotless Fixation Device originally cleared under K090615.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

Arthrocare Corporation  
% Ms. Laura Kasperowicz  
Senior Manager, Regulatory Affairs  
680 Vaqueros Avenue  
Sunnyvale, California 94085-3523

AUG - 9 2011

Re: K111044

Trade/Device Name: Speedlock® Knotless Fixation Device  
Regulation Number: 21 CFR 888.3040  
Regulation Name: Smooth or threaded metallic bone fixation fastener  
Regulatory Class: Class II  
Product Code: MBI  
Dated: June 30, 2011  
Received: July 05, 2011

Dear Ms. Kasperowicz:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic  
and Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## INDICATIONS FOR USE STATEMENT

510(k) Number: K 111044

Device Name: SpeedLock® Knotless Fixation System

Indications for Use:

The SpeedLock Knotless Fixation Device is indicated for use in fixation of soft tissue to bone. Examples of such procedures include:

**Shoulder:** Bankart Repair, SLAP lesion repair, acromio-clavicular separation, rotator cuff repair, capsule shift/capsule-labral reconstruction, biceps tenodesis and deltoid repair

**Ankle:** Lateral instability, medial instability, Achilles tendon repair/reconstruction and midfoot reconstruction

**Foot:** Hallux valgus reconstruction

**Elbow:** Tennis elbow repair, biceps tendon attachment

**Knee:** Extra-capsular repairs; reattachment of: medial collateral ligament, posterior oblique ligament or joint capsule closure to anterior proximal tibia; extra capsular reconstruction, ITB tenodesis; patellar ligament and tendon avulsions

|   |              |        |  |               |
|---|--------------|--------|--|---------------|
| Prescription Use<br>(Part 21 CFR 801 Subpart D) | <u>  X  </u> | AND/OR | Over-The-Counter Use<br>(21 CFR 801 Subpart C) | <u>  NO  </u> |
|---|--------------|--------|--|---------------|

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)  
Division of Surgical, Orthopedic,  
and Restorative Devices

510(k) Number 111044