



Food and Drug Administration 10903 New Hampshire Avenue Document Control Room – WO66-G609 Silver Spring, MD 20993-0002

Matang Manufacturing Sdn. Bhd Ms. Leslie Leong Operation Manager Lot 10, Kawasan Perindustrian Serkam Mukim Serkam, Merlimau Melaka, Malaysia 77300

OCT 1 7 2012

Re: K122131

Trade/Device Name: Powder Free Nitrile Examination Glove - Blue (Chemotherapy)

Regulation Number: 21 CFR 880.6250

Regulation Name: Patient Examination Glove

Regulatory Class: I

Product Code: LZA, LZC Dated: October 4, 2012 Received: October 12, 2012

Dear Ms. Leong:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to

http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm.

Sincerely yours,

Anthony D. Watson, B.S., M.S., M.B.A.

Director

Division of Anesthesiology, General Hospital, Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

Premarket Notification For Powder Free Nitrile Examination Glove - Blue (Chemotherapy)

3.0	Indication	For Use	Statement
J.U	HUUCAUUH	LOI CSC	Statement

510(K) Number (if known): <u>K12213</u>

Device Name: Powder Free Nitrile Examination Glove - Blue (Chemotherapy)

Indications For Use : A Powder Free Nitrile Examination Glove – Blue (Chemotherapy) is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between the patient and examiner.

In addition, these Chemotherapy gloves were tested for use with the

following drug concentrations per ASTM D 6978-05

Permeation Test Results

TEST CHEMOTHERAPY DRUG AND CONCENTRATION	MINIMUM BREAKTHROUGH DETECTION TIME (Specimen 1/2/3)(Minutes)
Carmustine (BCNU), 3.3 mg/ml (3,300 ppm)	60.3 (60.3, 60.5, 60.6)
Cisplatin, 1.0 mg/ml (1,000 ppm)	No breakthrough up to 240 min.
Cyclophosphamide (Cytoxan), 20.0 mg/ml (20,000ppm)	No breakthrough up to 240 min.
Dacarbazine (DTIC), 10.0 mg/ml (10,000 ppm)	No breakthrough up to 240 min.
Doxorubicin Hydrochloride, 2.0 mg/ml (2,000 ppm)	No breakthrough up to 240 min.
Etoposide (Toposar), 20.0 mg/ml (20,000 pm)	No breakthrough up to 240 min.
Fluorouracil, 50.0 mg/ml (50,000 pm)	No breakthrough up to 240 min.
Ifosfamide, 50.0 mg/ml (50,000 ppm)	No breakthrough up to 240 min.
Methotrexate, 25 mg/ml (25,000 ppm)	No breakthrough up to 240 min.
Mitomycin C, 0.5 mg/ml (500 ppm)	No breakthrough up to 240 min.
Mitoxantrone, 2 mg/ml (2,000 ppm)	No breakthrough up to 240 min.
Paclitaxel (Taxol), 6.0 mg/ml (6,000 ppm)	No breakthrough up to 240 min.
Thiotepa, 10.0 mg/ml (10,000 ppm)	156.0 (163.8, 156.0, 165.0)
Vincristine Sulfate, 1.0 mg/ml (1,000 ppm)	No breakthrough up to 240 min.

Caution: Testing showed the following drugs have low average breakthrough detection time:

* Carmustine - 60.3 min

* Thiotepa – 156.0 min

Prescription Use	AND/OR	Over-The-Counter Use	
(Part 21 CFR 801 Subpart D)		(21 CFR 801 Subpart C)	

AGE OF NEEDED

(Part 21 CFR 801 Subpart D)	(21 CFR 801 Subj
 PLEASE DO NOT WRITE BELOW THIS LINE-COM	VTINUE ON ANOTHER PA
Concurrence of CDRH, Off (Division Sign-Off) Division of Anesthesiology, General Hospital Infection Control, Dental Devices 510(k) Number: K12213	ice of Device Evaluation (ODE)