

Section 5. 510(k) Summary

MAR 14 2013

510(k) Owner

Pivot Medical Inc.
247 Humboldt Court
Sunnyvale CA 94089
Phone: 408-774-1452
Fax: 408-739-4199
Jon Cook
Director Regulatory Affairs and Quality Assurance

FDA Contact

Jon Cook
Director Regulatory Affairs and Quality Assurance
Pivot Medical
247 Humboldt Court
Sunnyvale CA 94089
Telephone: (408) 774-1452
Facsimile: (408) 739-4199
Email: jcook@pivotmedical.com

Date Summary Prepared: February 8, 2013

Device Name

Trade Name: NanoTack Suture Anchor 1.4mm with Flex Inserter
Common Name: Bone Anchor
Classification Name: Smooth or threaded metallic bone fixation fastener
Regulation number: 21 CFR888.3040
Product Code: MBI

Predicate Device

NanoTack Suture Anchor 1.4mm – K110473

Device Description

The Pivot NanoTack Suture Anchor 1.4mm is a non-degradable suture anchor manufactured from PEEK-OPTIMA® polymer attached to a non-degradable ultra high molecular weight polyethylene (UHMWPE) blue co-braid #1 suture that is pre-assembled to a Flex stainless steel Inserter. The NanoTack Suture Anchor 1.4mm with Flex Inserter is provided as a single use sterile device.

Intended Use / Indication for Use

The NanoTack Suture Anchor is intended for the fixation of soft tissue to bone in the hip, and is indicated for the reattachment of hip labrum to the acetabulum.

Summary of Technological Characteristics

The NanoTack Suture Anchor 1.4mm with Flex Inserter is substantially equivalent in design, materials, and function to the NanoTack Suture Anchor 1.4mm with Inserter (K110473). The only difference between the NanoTack Suture Anchor 1.4mm with Flex Inserter and the predicate is the laser cut pattern on the distal end of the Inserter shaft which allows the Inserter to flex through a curved guide.

Summary of Performance Testing

The performance testing conducted demonstrates that the insertion and fixation properties of the NanoTack Suture Anchor 1.4mm with Flex Inserter are substantially equivalent to the predicate device. Pre-clinical testing includes insertion strength and anchor strength.

Summary of Substantial Equivalence

The NanoTack Suture Anchor 1.4mm with Flex Inserter design, materials, and technology is substantially equivalent to the NanoTack Suture Anchor 1.4mm cleared in K110473. The modification to the inserter shaft does not raise new questions regarding the safety and effectiveness of the NanoTack Suture Anchor 1.4mm. Based upon the indications for use, technological characteristics, and comparison to the predicate device, the Pivot NanoTack Suture Anchor 1.4mm with Flex Inserter is substantially equivalent to the predicate device.



March 14, 2013

Pivot Medical Incorporated
% Mr. Jon Cook
Director Regulatory Affairs and Quality Assurance
247 Humboldt Court
Sunnyvale, California 94089

Re: K130351

Trade/Device Name: NanoTack Suture Anchor
Regulation Number: 21 CFR 888.3040
Regulation Name: Smooth or threaded metallic bone fixation fastener
Regulatory Class: Class II
Product Code: MBI
Dated: February 8, 2013
Received: February 12, 2013

Dear Mr. Cook:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA).

You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

Page 2 – Mr. Jon Cook

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Erin  Keith

Mark N. Melkerson
Director
Division of Orthopedic Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Section 4. **Indications for Use**

510(k) Number (if known): **K130351**

Device Name: NanoTack Suture Anchor

Indications for Use:

The NanoTack Suture Anchor is indicated for the reattachment of hip labrum to the acetabulum.

Prescription Use **X**
(Part 21 CFR 801 Subpart D)

AND / OR

Over-The-Counter-Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Casey  Ph.D.
Division of Orthopaedic Devices