

K132171
Page 1 of 5

510(k) Summary of Safety and Effectiveness

This summary of safety and effectiveness information is submitted in accordance with 21 CFR §807.92.

1) Submitter's name, address, telephone number, contact person

Submitted by:

SuperSonic Imagine, S.A.
Les Jardins de la Duranne – Bât. E & F
510, rue René Descartes
13857 Aix-en-Provence Cedex
France
Telephone: +33 442 99 24 24

Distributed by:

SuperSonic Imagine, Inc.
11714 North Creek Parkway N
Suite 150
Bothell, WA 98011
North America
Telephone: +1(425) 686 6380

SEP 24 2013

Corresponding Official:

Jacques Souquet
Chief Executive Officer
Telephone: +33 442 99 24 35

Date: 2013/07/12

2) Name of the device, including the trade or proprietary name if applicable, the common or usual name, and the classification name, if known:

Common/Usual Name: Diagnostic Ultrasound System with Accessories
Proprietary Name: Aixplorer®

Classification:

Regulatory Class: II

| Classification Name: | 21 CFR Section | Product Code |
|--|----------------|--------------|
| Ultrasonic Pulsed Doppler Imaging System | 892.1550 | 90-IYN |
| Ultrasonic Pulsed Echo Imaging System | 892.1560 | 90-IYO |
| Diagnostic Ultrasound Transducer | 892.1570 | 90-ITX |

3) Substantially Equivalent/Predicate Devices

AIXPLORER® Ultrasound Imaging System (K121329), cleared on 08/24/2012
AIXPLORER® Ultrasound Imaging System (K112255), cleared on 08/28/2012
Siemens Acuson S2000TM Diagnostic Ultrasound System (K072786), cleared on 11/13/2007
Philips iU22 Ultrasound System (K093563), cleared on 02/01/2010

4) Description of Device

The SuperSonic Imagine AIXPLORER® system is a cart based ultrasound imaging system used to perform non-invasive diagnostic general purpose ultrasound imaging studies. The system contains a scan converter and can be coupled to a variety of linear, curved, micro-convex, and motorized linear array transducers to produce images, which are displayed on a LCD monitor. An adjustable control panel with integrated touch screen allows the user to perform an ultrasound exam quickly and efficiently in accordance with ALARA principles. The system also allows the user to perform measurements, capture images to digital memory or to an external device (such as a printer), and review diagnostic studies in the form of a report. The system functions in a manner identical to the predicate devices and transducers for the imaging modes: B-Mode, M-mode, Color Flow, Pulsed Wave Doppler, Harmonic Imaging, Amplitude Doppler, 3D imaging and for ShearWave™ elastography.

5) Intended Use

The SuperSonic Imagine AIXPLORER® ultrasound system and transducer are intended for general purpose pulse echo ultrasound imaging and Doppler fluid flow analysis of the human body.

The SuperSonic Imagine AIXPLORER® ultrasound system is indicated for use in the following applications: Abdominal, Small Organs, Musculoskeletal, Superficial Musculoskeletal, Vascular, Peripheral Vascular, OB-GYN, Pelvic, Pediatric, Urology, Trans-rectal, Trans-vaginal and Neonatal Cephalic.

The system also provides the ability to measure anatomical structures (Abdominal, Small Organs, Musculoskeletal, Superficial Musculoskeletal, Peripheral Vascular, GYN, Pelvic, Pediatric, Urology, Trans-rectal, Trans-vaginal, Neonatal Cephalic, Fetal/Obstetrics).

6) Summary of Technological Characteristics – New Device compared to Predicates

| | Philips iu22 (predicate K093563) | Siemens Acuson S2000™ (predicate K072786) | SuperSonic Imagine AIXPLORER® (predicate K121329) | SuperSonic Imagine AIXPLORER® (predicate K112255) | SuperSonic Imagine AIXPLORER® (submission device) |
|-----------------------------------|--|--|---|---|---|
| Clinical Applications | ---** | General Radiology | ---** | ---** | ---** |
| | Abdominal, | Identical | Identical | Identical | Identical |
| | Small Organs* | Identical | Identical | Identical | Identical |
| | Musculoskeletal | Identical | Identical | Identical | Identical |
| | --- | Superficial Musculoskeletal | Identical | Identical | Identical |
| | Fetal | Identical | --- | --- | Identical |
| | --- | Transcranial | --- | --- | --- |
| | --- | OB | --- | --- | Identical |
| | Cardiac | Identical | Identical | Identical | Identical |
| | --- | Pelvic | Identical | Identical | Identical |
| Adult and neonatal cephalic | Identical | Identical (for neonatal cephalic) | --- | Identical (for neonatal cephalic) | |
| Pediatric | Identical | Identical | Identical | Identical | |

| | | | | | |
|-------------------------|-----------------------|---------------------|-----------|-----------|-----------|
| | Urology | Identical | Identical | Identical | Identical |
| | — | Vascular | Identical | Identical | Identical |
| | Peripheral Vascular | Identical | Identical | Identical | Identical |
| | Ophthalmic | --- | --- | --- | --- |
| | Intra-operative | --- | --- | --- | --- |
| | Laparoscopic | --- | --- | --- | --- |
| | Trans-rectal | --- | Identical | Identical | Identical |
| | Trans-vaginal | --- | Identical | Identical | Identical |
| | Fetal echo | --- | --- | --- | Identical |
| | | | | | |
| Imaging Modes | | | | | |
| Conventional | B-mode, | Identical | Identical | Identical | Identical |
| | M-mode, | Identical | --- | --- | Identical |
| | PW, | Identical | Identical | Identical | Identical |
| | CW (continuous Wave), | Identical | --- | --- | --- |
| | Color Doppler, | Identical | Identical | Identical | Identical |
| | Amplitude Doppler | Identical | Identical | Identical | Identical |
| Other | Harmonic imaging, | Identical | Identical | Identical | Identical |
| | Spatial Compounding, | Identical | Identical | --- | Identical |
| | Panoramic, | Identical | --- | --- | --- |
| | Contrast | Identical | Identical | --- | --- |
| | --- | Elastography | Identical | Identical | Identical |
| Combination | B-mode+Color, | Identical, | Identical | Identical | Identical |
| | B-mode+Color+PW | Identical | Identical | --- | Identical |
| | B-mode +PW | Identical | Identical | Identical | Identical |
| | B-mode+M-mode | Identical | --- | --- | Identical |
| | --- | B-mode+Elastography | Identical | Identical | Identical |
| | | | | | |
| Transducers | | | | | |
| Transducer types | Linear Array | Identical | Identical | Identical | Identical |
| | Curved Array | Identical | Identical | Identical | Identical |
| | Phased Array | --- | --- | --- | --- |
| | Laparoscopic probe | --- | --- | --- | --- |
| | | | | | |

| | | | | | |
|---------------------------------------|---|------------------------|------------------------|------------------------|------------------------|
| | Motorized Linear Probe Microconvex probe | Identical Identical | Identical Identical | Identical Identical | Identical Identical |
| Track | Track 3 (Acoustic Output Display) | Identical | Identical | Identical | Identical |
| Patient Contact Materials | Yes, per ISO-10993-1 | Identical | Identical | Identical | Identical |
| Acoustic Output within FDA guidelines | Yes, as per NEMA UD-3 | Identical | Identical | Identical | Identical |
| General Safety | Conforms to IEC 60601-1, IEC 60601-2 | Identical | Identical | Identical | Identical |

Note:

*: Breast, Thyroid, Testicle, etc

** : --- means not applicable

7) A brief discussion of the clinical tests submitted, referenced, or relied on in the premarket notification submission for a determination of substantial equivalence

Non-clinical testing was conducted per the following standards to support a determination of substantial equivalence to the predicate devices.

| Reference Standard | Tests Performed |
|--|---|
| IEC 60601-1 3 rd Edition | All applicable electrical, basic safety and essential performance tests. |
| UL 60601-1 1 st Edition | All applicable electrical, basic safety and essential performance tests specific to the U.S.A. |
| IEC 60601-1-1 2 nd Edition | All applicable tests pertaining to Medical Electrical Systems. |
| IEC 60601-1-2 3 rd Edition | All applicable testing pertaining to electromagnetic compatibility. |
| IEC 60601-2-37 2 nd Edition | All applicable testing pertaining to the particular requirements for the safety of ultrasonic medical diagnostic and monitoring equipment. |
| NEMA UD 2 (Rev. 3) | All tests applicable in order to demonstrate compliance with the "Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment". |
| NEMA UD 3 (Rev. 2) | All tests applicable in order to demonstrate compliance with the "Standard For Real Time Display Of Thermal And Mechanical Acoustic Output Indices On Diagnostic Ultrasound Equipment". |
| ISO 10993-1 | Applicable biocompatibility tests per FDA 510(k) Memorandum - #G95-1 – per the appropriate |

K132171
Page 5 of 5

| | |
|--|------------------|
| | device category. |
|--|------------------|

In addition to the referenced standards testing, performance tests were conducted with respect to Fetal/Obstetrics features.

The above testing confirmed that the Aixplorer System performs according to the stated intended use. All data fell within pre-determined product specifications and external standard requirements. Results of non-clinical testing confirmed the substantial equivalence of the Aixplorer System to the predicate device(s).

8) A brief discussion of the clinical tests submitted, referenced, or relied on in the premarket notification submission for a determination of substantial equivalence

Clinical data is not required as the Aixplorer System uses the same technology and principles as predicate devices.

9) Conclusion

The manufacturer and the design and development of the submission device comply with 21 CFR Part 820 and ISO 13485 (2003) Quality Standards. The submission device, designed to comply with applicable safety standards, is tested during manufacturing process to ensure compliance with these standards. Consequently, according tests performed, the opinion of SuperSonic Imagine is the submission device is as safe and effective as the predicate devices cited in item 3.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - W066-G609
Silver Spring, MD 20993-0002

September 24, 2013

Supersonic Imagine, S.A.
% Mr. Aurelie Gruener
Les Jardins de la Duranne
510 Rue René Descartes – Bât. E et F
Aix -en-Provence Cedex 13 857
FRANCE

Re: K132171

Trade/Device Name: Aixplorer®
Regulation Number: 21 CFR 892.1550
Regulation Name: Ultrasonic pulsed doppler imaging system
Regulatory Class: II
Product Code: IYN, IYO, ITX
Dated: September 3, 2013
Received: September 3, 2013

Dear Mr. Gruener:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

This determination of substantial equivalence applies to the following transducers intended for use with the Aixplorer®, as described in your premarket notification:

Transducer Model Number

SL15-4
SC6-1

SE12-3
SLV16-5

SL10-2
SMC12-3

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



for

Janine M. Morris
Director, Division of Radiological Health
Office of In Vitro Diagnostics
and Radiological Health
Center for Devices and Radiological Health

Enclosure

Diagnostic Ultrasound Indications for Use

510(k) number (if known): XXXXXXXXXX K132171

Device Name: AIXPLORER® Ultrasound System

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Indications for Use:

The SuperSonic Imagine AIXPLORER® ultrasound system is indicated for use in the following applications: Abdominal, Small Organs, Musculoskeletal, Superficial Musculoskeletal, Vascular, Peripheral Vascular, OB-GYN, Pelvic, Pediatric, Urology, Trans-rectal, Trans-vaginal and Neonatal Cephalic.

The system also provides the ability to measure anatomical structures (Abdominal, Small Organs, Musculoskeletal, Superficial Musculoskeletal, Peripheral Vascular, GYN, Pelvic, Pediatric, Urology, Trans-rectal, Trans-vaginal, Neonatal Cephalic, Fetal/Obstetrics).


Prescription Use (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH; Office of In Vitro Diagnostics and Radiological Health (OIR)



(Division Sign-Off)
Division of Radiological Health
Office of In Vitro Diagnostics and Radiological Health
510(k) _____ K132171 _____

Diagnostic Ultrasound Indications for Use

510(k) number (if known):

Device Name: AIXPLORER® Ultrasound System

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

| Clinical Application | | Mode of Operation | | | | | | |
|------------------------|--|-------------------|---|-----|-----|---------------|--------------------|----------------------|
| General (Track 1 Only) | Specific (Tracks 1 & 3) | B | M | PWD | CWD | Color Doppler | Combined (Specify) | Other* (Specify) |
| Ophthalmic | Ophthalmic | | | | | | | |
| Fetal Imaging & Other | Fetal | N | N | N | | N | N, 1, 3, 4, 11 | N, 5, 6 |
| | Abdominal (including urology) | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 10, 9 |
| | Intra-operative (Specify) | | | | | | | |
| | Intra-operative (Neuro) | | | | | | | |
| | Laparoscopic | | | | | | | |
| | Pediatric | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 10, 9 |
| | Small Organ (Breast, Thyroid, Testicle, Prostate, penis, etc...) | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 10, 9 |
| | Neonatal Cephalic | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 9 |
| | Adult Cephalic | | | | | | | |
| | Trans-rectal | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8 |
| | Trans-vaginal | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8 |
| | Trans-urethral | | | | | | | |
| | Trans-esoph. (non-Card.) | | | | | | | |
| | Musculo-skeletal (Conventional) | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 10, 9 |
| | Musculo-skeletal (Superficial) | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 10, 9 |
| | Intravascular | | | | | | | |
| GYN | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 10 | |
| Pelvic | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 10 | |
| Other (Specify) | | | | | | | | |
| Cardiac | Cardiac Adult | | | | | | | |
| | Cardiac Pediatric | | | | | | | |
| | Intravascular (Cardiac) | | | | | | | |
| | Trans-esoph. (Cardiac) | | | | | | | |
| | Intra-cardiac | | | | | | | |
| | Other (Specify) | | | | | | | |
| Peripheral Vessel | Peripheral vessel | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 10, 9 |
| | Other (Specify) | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |

N = new indication; P = previously cleared by FDA (K121329)

Additional Comments:

1: Combined modes include: B+ Color Flow

2: Combined modes include: B+ ShearWave™ Elastography

3: Combined modes include: B+ Pulsed Wave

4: Combined modes include: B+ Pulsed Wave + Color Flow

5: Harmonic Imaging

6: Spatial Compounding

7: ShearWave™ Elastography

8: Imaging Guidance for Biopsies

9: Panoramic Imaging

10: 3D Imaging

11: Combined modes include: B+ M modes

Prescription Use (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH: Office of In Vitro Diagnostics and Radiological Health (OIR)

Diagnostic Ultrasound Indications for Use

510(k) Number (if known):

Device Name: SL15-4 transducer (1D Linear Array Transducer)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

| Clinical Application | | Mode of Operation | | | | | | |
|------------------------|--|-------------------|---|-----|-----|---------------|--------------------|------------------|
| General (Track 1 Only) | Specific (Tracks 1 & 3) | B | M | PWD | CWD | Color Doppler | Combined (Specify) | Other* (Specify) |
| Ophthalmic | Ophthalmic | | | | | | | |
| Fetal Imaging & Other | Fetal | | | | | | | |
| | Abdominal | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 9 |
| | Intra-operative (Specify) | | | | | | | |
| | Intra-operative (Neuro) | | | | | | | |
| | Laparoscopic | | | | | | | |
| | Pediatric | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 9 |
| | Small Organ (for example Breast, Thyroid, Testicle, Prostate, Penis) | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 9 |
| | Neonatal Cephalic | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 9 |
| | Adult Cephalic | | | | | | | |
| | Trans-rectal | | | | | | | |
| | Trans-vaginal | | | | | | | |
| | Trans-urethral | | | | | | | |
| | Trans-esoph. (non-Card.) | | | | | | | |
| | Musculo-skeletal (Conventional) | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 9 |
| | Musculo-skeletal (Superficial) | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 9 |
| | Intravascular | | | | | | | |
| | GYN | | | | | | | |
| Pelvic | | | | | | | | |
| Other (Specify) | | | | | | | | |
| Cardiac | Cardiac Adult | | | | | | | |
| | Cardiac Pediatric | | | | | | | |
| | Intravascular (Cardiac) | | | | | | | |
| | Trans-esoph. (Cardiac) | | | | | | | |
| | Intra-cardiac | | | | | | | |
| Other (Specify) | | | | | | | | |
| Peripheral | Peripheral vessel | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 9 |
| Vessel | Other (Specify) | | | | | | | |

N = new indication; P = previously cleared by FDA (K121329)

Additional Comments:

- 1: Combined modes include: B+ Color Flow
- 2: Combined modes include: B+ ShearWave™ Elastography
- 3: Combined modes include: B+ Pulsed Wave
- 4: Combined modes include: B+ Pulsed Wave + Color Flow
- 5: Harmonic Imaging

- 6: Spatial Compounding
- 7: ShearWave™ Elastography
- 8: Imaging Guidance for Biopsies
- 9: Panoramic Imaging
- 10: 3D Imaging
- 11: Combined modes include: B+ M modes

Prescription Use (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____ (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH: Office of In Vitro Diagnostics and Radiological Health (OIR)

Diagnostic Ultrasound Indications for Use

510(k) Number (if known):

Device Name: SC6-1 transducer (curved array transducer)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

| General (Track 1 Only) | Clinical Application Specific (Tracks 1 & 3) | Mode of Operation | | | | | | |
|---------------------------|---|-------------------|---|-----|-----|------------------|-----------------------|------------------|
| | | B | M | PWD | CWD | Color Doppler | Combined (Specify) | Other* (Specify) |
| Ophthalmic | Ophthalmic | | | | | | | |
| Fetal Imaging & Other | Fetal | N | N | N | | N | N, 1, 3, 4, 11 | N, 5, 6 |
| | Abdominal (including urology) | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 9 |
| | Intra-operative (Specify) | | | | | | | |
| | Intra-operative (Neuro) | | | | | | | |
| | Laparoscopic | | | | | | | |
| | Pediatric | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 9 |
| | Small Organ (Breast, Thyroid, Testicle, Prostate, penis, etc...) | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8 |
| | Neonatal Cephalic | | | | | | | |
| | Adult Cephalic | | | | | | | |
| | Trans-rectal | | | | | | | |
| | Trans-vaginal | | | | | | | |
| | Trans-urethral | | | | | | | |
| | Trans-esoph. (non-Card.) | | | | | | | |
| | Musculo-skeletal (Conventional) | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 9 |
| | Musculo-skeletal (Superficial) | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 9 |
| | Intravascular | | | | | | | |
| | GYN | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8 |
| Pelvic | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8 | |
| Other (Specify) | | | | | | | | |
| Cardiac | Cardiac Adult | | | | | | | |
| | Cardiac Pediatric | | | | | | | |
| | Intravascular (Cardiac) | | | | | | | |
| | Trans-esoph. (Cardiac) | | | | | | | |
| | Intra-cardiac | | | | | | | |
| Other (Specify) | | | | | | | | |
| Peripheral | Peripheral vessel | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 9 |
| Vessel | Other (Specify) | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 9 |

N = new indication; P = previously cleared by FDA (K121329)

Additional Comments:

- 1: Combined modes include: B+ Color Flow
- 2: Combined modes include: B+ ShearWave™ Elastography
- 3: Combined modes include: B+ Pulsed Wave
- 4: Combined modes include: B+ Pulsed Wave + Color Flow
- 5: Harmonic Imaging

- 6: Spatial Compounding
- 7: ShearWave™ Elastography
- 8: Imaging Guidance for Biopsies
- 9: Panoramic Imaging
- 10: 3D Imaging
- 11: Combined modes include: B+ M modes

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH: Office of In Vitro Diagnostics and Radiological Health (OIR)

Diagnostic Ultrasound Indications for Use

510(k) Number (if known):

Device Name: SE12-3 transducer (endocavitary transducer)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

| General (Track 1 Only) | Clinical Application Specific (Tracks 1 & 3) | Mode of Operation | | | | | | |
|---------------------------|--|-------------------|---|-----|-----|---------------|--------------------|------------------|
| | | B | M | PWD | CWD | Color Doppler | Combined (Specify) | Other* (Specify) |
| Ophthalmic | Ophthalmic | | | | | | | |
| Fetal Imaging & Other | Fetal | N | N | | | N | N, 1, 3, 4, 11 | N, 5, 6 |
| | Abdominal | | | | | | | |
| | Intra-operative (Specify) | | | | | | | |
| | Intra-operative (Neuro) | | | | | | | |
| | Laparoscopic | | | | | | | |
| | Pediatric | | | | | | | |
| | Small Organ (for example Breast, Thyroid, Testicle, Prostate, penis, etc...) | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8 |
| | Neonatal Cephalic | | | | | | | |
| | Adult Cephalic | | | | | | | |
| | Trans-rectal | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8 |
| | Trans-vaginal | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8 |
| | Trans-urethral | | | | | | | |
| | Trans-esoph. (non-Card.) | | | | | | | |
| | Musculo-skeletal (Conventional) | | | | | | | |
| | Musculo-skeletal (Superficial) | | | | | | | |
| | Intravascular | | | | | | | |
| | GYN | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8 |
| | Pelvic | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8 |
| | Other (Specify) | | | | | | | |
| Cardiac | Cardiac Adult | | | | | | | |
| | Cardiac Pediatric | | | | | | | |
| | Intravascular (Cardiac) | | | | | | | |
| | Trans-esoph. (Cardiac) | | | | | | | |
| | Intra-cardiac | | | | | | | |
| | Other (Specify) | | | | | | | |
| Peripheral Vessel | Peripheral vessel | | | | | | | |
| | Other (Specify) | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8 |

N = new indication; P = previously cleared by FDA (K121329)

Additional Comments:

- 1: Combined modes include: B+ Color Flow
- 2: Combined modes include: B+ ShearWave™ Elastography
- 3: Combined modes include: B+ Pulsed Wave
- 4: Combined modes include: B+ Pulsed Wave + Color Flow
- 5: Harmonic Imaging

- 6: Spatial Compounding
- 7: ShearWave™ Elastography
- 8: Imaging Guidance for Biopsies
- 9: Panoramic Imaging
- 10: 3D Imaging
- 11: Combined modes include: B+ M modes

Prescription Use (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____ (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH; Office of In Vitro Diagnostics and Radiological Health (OIR)

Diagnostic Ultrasound Indications for Use

510(k) Number (if known):

Device Name: SLV16-5 transducer (motorized linear transducer)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

| Clinical Application | | Mode of Operation | | | | | | | |
|------------------------|--|-------------------|---|-----|-----|---------------|--------------------|---------------------|---------------------|
| General (Track 1 Only) | Specific (Tracks 1 & 3) | B | M | PWD | CWD | Color Doppler | Combined (Specify) | Other* (Specify) | |
| Ophthalmic | Ophthalmic | | | | | | | | |
| Fetal Imaging & Other | Fetal | | | | | | | | |
| | Abdominal | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 10, 9 | |
| | Intra-operative (Specify) | | | | | | | | |
| | Intra-operative (Neuro) | | | | | | | | |
| | Laparoscopic | | | | | | | | |
| | Pediatric | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 10, 9 | |
| | Small Organ (for example Breast, Thyroid, Testicle, Prostate, penis, etc...) | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 10, 9 | |
| | Neonatal Cephalic | | | | | | | | |
| | Adult Cephalic | | | | | | | | |
| | Trans-rectal | | | | | | | | |
| | Trans-vaginal | | | | | | | | |
| | Trans-urethral | | | | | | | | |
| | Trans-esoph. (non-Card.) | | | | | | | | |
| | Musculo-skeletal (Conventional) | P | | P | | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 10, 9 |
| | Musculo-skeletal (Superficial) | P | | P | | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 10, 9 |
| | Intravascular | | | | | | | | |
| | GYN | | | | | | | | |
| Pelvic | | | | | | | | | |
| Other (Specify) | | | | | | | | | |
| Cardiac | Cardiac Adult | | | | | | | | |
| | Cardiac Pediatric | | | | | | | | |
| | Intravascular (Cardiac) | | | | | | | | |
| | Trans-esoph. (Cardiac) | | | | | | | | |
| | Intra-cardiac | | | | | | | | |
| Other (Specify) | | | | | | | | | |
| Peripheral | Peripheral vessel | P | | P | | P | P 1, 2, 3, 4 | P 5, 6, 7, 8, 10, 9 | |
| Vessel | Other (Specify) | | | | | | | | |

N = new indication; P = previously cleared by FDA (K121329)

Additional Comments:

- 1: Combined modes include: B+ Color Flow
- 2: Combined modes include: B+ ShearWave™ Elastography
- 3: Combined modes include: B+ Pulsed Wave
- 4: Combined modes include: B+ Pulsed Wave + Color Flow
- 5: Harmonic Imaging

- 6: Spatial Compounding
- 7: ShearWave™ Elastography
- 8: Imaging Guidance for Biopsies
- 9: Panoramic Imaging
- 10: 3D Imaging
- 11: Combined modes include: B+ M modes

Prescription Use (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____ (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH: Office of In Vitro Diagnostics and Radiological Health (OIR)

Diagnostic Ultrasound Indications for Use

510(k) Number (if known):

Device Name: SL10-2 transducer (linear transducer)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

| Clinical Application | | Mode of Operation | | | | | | |
|------------------------|--|-------------------|---|-----|-----|---------------|--------------------|------------------|
| General (Track 1 Only) | Specific (Tracks 1 & 3) | B | M | PWD | CWD | Color Doppler | Combined (Specify) | Other* (Specify) |
| Ophthalmic | Ophthalmic | | | | | | | |
| Fetal Imaging & Other | Fetal | | | | | | | |
| | Abdominal | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |
| | Intra-operative (Specify) | | | | | | | |
| | Intra-operative (Neuro) | | | | | | | |
| | Laparoscopic | | | | | | | |
| | Pediatric | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |
| | Small Organ (for example Breast, Thyroid, Testicle, Prostate, penis, etc...) | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |
| | Neonatal Cephalic | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 9 |
| | Adult Cephalic | | | | | | | |
| | Trans-rectal | | | | | | | |
| | Trans-vaginal | | | | | | | |
| | Trans-urethral | | | | | | | |
| | Trans-esoph. (non-Card.) | | | | | | | |
| | Musculo-skeletal (Conventional) | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |
| | Musculo-skeletal (Superficial) | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |
| Intravascular | | | | | | | | |
| GYN | | | | | | | | |
| Pelvic | | | | | | | | |
| Other (Specify) | | | | | | | | |
| Cardiac | Cardiac Adult | | | | | | | |
| | Cardiac Pediatric | | | | | | | |
| | Intravascular (Cardiac) | | | | | | | |
| | Trans-esoph. (Cardiac) | | | | | | | |
| | Intra-cardiac | | | | | | | |
| Other (Specify) | | | | | | | | |
| Peripheral | Peripheral vessel | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |
| Vessel | Other (Specify) | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |

N = new indication; P = previously cleared by FDA (K121329)

Additional Comments:

- 1: Combined modes include: B+ Color Flow
- 2: Combined modes include: B+ ShearWave™ Elastography
- 3: Combined modes include: B+ Pulsed Wave
- 4: Combined modes include: B+ Pulsed Wave + Color Flow
- 5: Harmonic Imaging

- 6: Spatial Compounding
- 7: ShearWave™ Elastography
- 8: Imaging Guidance for Biopsies
- 9: Panoramic Imaging
- 10: 3D Imaging
- 11: Combined modes include: B+ M modes

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH: Office of In Vitro Diagnostics and Radiological Health (OIR)

Diagnostic Ultrasound Indications for Use

510(k) Number (if known):

Device Name: SMC12-3 transducer (micro-curved transducer)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

| Clinical Application | | Mode of Operation | | | | | | |
|------------------------|--|-------------------|---|-----|-----|---------------|--------------------|------------------|
| General (Track 1 Only) | Specific (Tracks 1 & 3) | B | M | PWD | CWD | Color Doppler | Combined (Specify) | Other* (Specify) |
| Ophthalmic | Ophthalmic | | | | | | | |
| Fetal Imaging & Other | Fetal | | | | | | | |
| | Abdominal | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |
| | Intra-operative (Specify) | | | | | | | |
| | Intra-operative (Neuro) | | | | | | | |
| | Laparoscopic | | | | | | | |
| | Pediatric | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |
| | Small Organ (for example Breast, Thyroid, Testicle, Prostate, penis, etc...) | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |
| | Neonatal Cephalic | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 9 |
| | Adult Cephalic | | | | | | | |
| | Trans-rectal | | | | | | | |
| | Trans-vaginal | | | | | | | |
| | Trans-urethral | | | | | | | |
| | Trans-esoph. (non-Card.) | | | | | | | |
| | Musculo-skeletal (Conventional) | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |
| | Musculo-skeletal (Superficial) | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |
| | Intravascular | | | | | | | |
| GYN | | | | | | | | |
| Pelvic | | | | | | | | |
| Other (Specify) | | | | | | | | |
| Cardiac | Cardiac Adult | | | | | | | |
| | Cardiac Pediatric | | | | | | | |
| | Intravascular (Cardiac) | | | | | | | |
| | Trans-esoph. (Cardiac) | | | | | | | |
| | Intra-cardiac | | | | | | | |
| Other (Specify) | | | | | | | | |
| Peripheral | Peripheral vessel | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |
| Vessel | Other (Specify) | P | | P | | P | P, 1, 2, 3, 4 | P, 5, 6, 7, 8, 9 |

N = new indication; P = previously cleared by FDA (K121329)

Additional Comments:

- 1: Combined modes include: B+ Color Flow
- 2: Combined modes include: B+ ShearWave™ Elastography
- 3: Combined modes include: B+ Pulsed Wave
- 4: Combined modes include: B+ Pulsed Wave + Color Flow
- 5: Harmonic Imaging

- 6: Spatial Compounding
- 7: ShearWave™ Elastography
- 8: Imaging Guidance for Biopsies
- 9: Panoramic Imaging
- 10: 3D Imaging
- 11: Combined modes include: B+ M modes

Prescription Use (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____ (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH: Office of In Vitro Diagnostics and Radiological Health (OIR)