


Labeling of Patient Spoon



Patient Spoon labeling

Method: laser engraving, Size: L 7 mm x W 7 mm Data-Matrix code and serial number.

Labeling of Docking Station



Product label Docking Station



The rating plate is glued on the back of the device housing of the Docking Station. The barcode 128 with the HW-ID (individual hardware identification of the circuit board) is affixed by the manufacturer (EKM) on the rating plate in the field next to the symbol for the serial number on the rear of the device housing. The algorithms and the above representation of the SN are described in the Labelling specification.