



SEP 30 2003

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Richard Woortman
Manager Technical Services
Cavex Holland B.V.
Harmanjansweg 19
P.O. Box 852
Haarlem, Holland

Re: K032116
Trade/Device Name: Tulip Colorswitch Alginate Impression Materials
Regulation Number: 872.3660
Regulation Name: Impression Material
Regulatory Class: II
Product Code: ELW
Dated: July 7, 2003
Received: September 17, 2003

Dear Mr. Woortman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4613. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink that reads "Susan Runner". The signature is written in a cursive style with a large, looped initial "S".

Susan Runner, DDS, MA
Interim Director
Division of Anesthesiology, General Hospital
Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) NUMBER (IF KNOWN): K032116

DEVICE NAME: Tulip Colorswitch
Alginate Impression Materials

INDICATIONS FOR USE:

Tulip Colorswitch is a dental impression material based on alginate. It is used for taking impressions of the oral cavity with the purpose of constructing a gypsum cast that is a copy of the situation in the mouth. It is a general purpose impression material for making study models, first impressions for the construction of individual trays, situation models, orthodontic impressions etc.

Tulip Colorswitch has the special characteristic of being a color-changing alginate: it becomes violet upon contact with water in the mixing bowl, changes to pink indicating the end of the mixing time and then finally changes to white indicating the end of the setting time in the mouth.

MCA Che for MSR

(Division Sign-Off)
Division of Anesthesiology, General Hospital,
Infection Control, Dental Devices

510(k) Number: K032116

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Prescription Use _____
(Per 21 CFR 801.109)

OR

Over-the-Counter Use _____