

AUG 25 2005

K052080

**Special 510(k) Summary  
Linvatec Biomaterials  
Matryx™ Interference Screw**

**Submitter's Name, Address, Telephone Number, and Contact Person**

Linvatec Biomaterials Ltd.  
Tuija Annala  
Director, Quality, Regulatory Affairs and Product  
Development  
P.O.Box 3  
FIN-33721 Tampere  
Finland, Europe  
Phone: 358-3-316 5600  
Facsimile: 358-3-316 5629

**Date prepared:** June 30<sup>th</sup>, 2005

**Name of the device:**

- A. Trade or Proprietary Name: Matryx™ Interference Screw
- B. Common Name: Bioabsorbable Interference Screw
- C. Classification Name: Bone Fixation Screw
- D. Device Product Code: MAI and HWC

**Predicate Device:**

1. Linvatec Biomaterials Ltd Osteo ACL Screw (K032894)

**Intended Use:**

The Matryx™ Interference Screw is intended for use in interference fixation of bone-patellar tendon – bone and soft tissue grafts in anterior and posterior cruciate ligament reconstructions.

The Matryx™ Interference Screw is not intended for use in and is contraindicated for 1) insufficient quality and quantity of bone for attachment of graft, 2) Blood supply limitation and/or previous infections, which could retard healing, 3) Foreign body sensitivity to the implant material. Where the material is suspected a test should be made prior to implantation to rule out sensitivity, 4) Patients with active sepsis or infection, 5) Conditions, which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing and rehabilitation period, 6) ACL repairs, which would not be appropriate for fixation with metallic screws.

**Device Description:**

The device description of the Matryx™ Interference Screw is as follows.

- The implant is composed of mixture of poly-96L/4D-lactide copolymer and tricalcium phosphate. The predicate device Osteo ACL Screw is made of the very same raw material.
- Lengths of implant are 20 - 30 mm
- Diameters of implant are 7mm-9mm.

The only modifications that were made are:

- Change of trade name. This change is updated in labelling.
- Minor modifications and further definition of tolerances of screw design.

**Substantial Equivalence:**

Linvatec Biomaterials Ltd Matryx Interference Screw is substantially equivalent to the cleared predicate device. The applied modifications do not raise any new concerns of safety and efficacy of the implant.



AUG 25 2005

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Tuija Annala  
Director, Quality, Regulatory Affairs  
and Product Development  
Linvatec Biomaterials Ltd.  
P.O. Box 3  
FIN-33721 Tampere  
Finland, Europe

Re: K052080  
Trade/Device Name: Matryx™ Interference Screw  
Regulation Number: 21 CFR 888.3040  
Regulation Name: Smooth or threaded metallic bone fixation fastener  
Regulatory Class: II  
Product Code: HWC, MAI  
Dated: June 30, 2005  
Received: August 02, 2005

Dear Ms. Annala:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the [Federal Register](#).


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Tuija Annala

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Mark N. Melkerson', is written over a horizontal line. The signature is stylized and cursive.

Mark N. Melkerson  
Acting Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## INDICATIONS FOR USE

510(K) Number (if known): K052080

Device Name: Matryx™ Interference Screw

### Indications for Use:

The Matryx™ Interference Screw is intended for use in interference fixation of bone-patellar tendon – bone and soft tissue grafts in anterior and posterior cruciate ligament reconstructions.

The Matryx™ Interference Screw is not intended for use in and is contraindicated for 1) insufficient quality and quantity of bone for attachment of graft, 2) Blood supply limitation and/or previous infections, which could retard healing, 3) Foreign body sensitivity to the implant material. Where the material is suspected a test should be made prior to implantation to rule out sensitivity, 4) Patients with active sepsis or infection, 5) Conditions, which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing and rehabilitation period, 6) ACL repairs, which would not be appropriate for fixation with metallic screws.

Prescription Use  AND/OR Over-The-Counter Use   
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)

**Division of General, Restorative,  
and Neurological Devices**

510(k) Number K052080