

**510(K) SUMMARY**

K060474

**Lisfranc Plates**

AUG - 3 2006

**Submitter's name and address:**

Newdeal SAS  
10, place d'Helvétie  
69006 Lyon, France  
Tel: +33 4 37 47 51 51  
Fax: +33 4 37 47 51 52

**Contact person and telephone number**

Morgane Grenier  
Regulatory and Clinical Affairs Director  
Newdeal SAS  
10, place d'Helvétie  
69006 Lyon, France  
Tel: +33 4 37 47 51 51  
Fax: +33 4 37 47 51 52

**Alternate Contacts****Authorized Agent in the United States**

Judith E. O'Grady, RN, MSN  
Sr. Vice President, Regulatory Affairs, Quality Assurance and Clinical Affairs  
Integra LifeSciences Corporation  
311 Enterprise Drive  
Plainsboro, NJ 08536, USA  
Tel: (609) 936-2311  
Fax: (609) 275-9445  
E-mail: [jogrady@integra-ls.com](mailto:jogrady@integra-ls.com)

**Date Summary was prepared:**

July 17, 2006

**Name of the device:**

Proprietary Name: Lisfranc Plate  
Common Name: Plate, Fixation, Bone  
Classification Name: Single/multiple component metallic bone fixation appliances and accessories (21CFR 888.3030)  
Device Product Code: HRS  
Classification Panel: Orthopedic

**Substantial Equivalence:**

The **Lisfranc Plate** is substantially equivalent to the Acumed Lower Extremity Congruent Bone Plate System, K033639, the Synthes Modular Foot System – 2.7 mm Module, K010321 and the Synthes Modular Foot System, K001941.

**Device Description:**

The NEWDEAL<sup>®</sup> Lisfranc Plate consists of an osteosynthesis plate designed to bridge the tarsometatarsal joints (Lisfranc joints). It is available in different sizes, and is implanted using NEWDEAL<sup>®</sup> locking system fixation screws and washers. The NEWDEAL<sup>®</sup> locking system includes as many fixation screws as there are threaded lipped sockets on the plate and as many washers as implanted screws. The NEWDEAL<sup>®</sup> locking system creates a single implant/screw unit fixed into the bone. The osteosynthesis screws must be driven into the bone through the holes in the plate. The system is locked by means of washers drilled into the threaded lipped socket at the top of each hole, thus blocking each screw head.

**Intended Use:**

The NEWDEAL Lisfranc Plates are intended to be used for fractures, fusions, osteotomies and replantations of small bones of the tarsometatarsal joints (Lisfranc joints).

**Testing and Test Results:**

Performance testing was conducted compared to a legally marketed predicate device. The results demonstrated equivalence between the Lisfranc Plates and the predicate device.

**Conclusion**

The Newdeal Lisfranc Plates are substantially equivalent to commercially marketed devices, the Acumed Lower Extremity Congruent Bone Plate System, K033639 the Synthes Modular Foot System – 2.7 mm Module, K010321 and the Synthes Modular Foot System, K001941.

The Newdeal Lisfranc Plates do not raise any new issues of scientific technology, safety or effectiveness



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

AUG - 3 2006

Newdeal SA  
c/o Ms. Judith O'Grady  
Sr. VP, Regulatory Affairs  
Integra Lifesciences Corporation  
311 Enterprise Drive  
Plainsboro, New Jersey 08536

Re: K060474  
Trade/Device Name: Lisfranc Plate  
Regulation Number: 21 CFR 888.3030  
Regulation Name: Single/multiple component metallic bone fixation appliances and accessories  
Regulatory Class: II  
Product Code: HRS  
Dated: July 17, 2006  
Received: July 18, 2006

Dear Ms. O'Grady:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must

comply with all the Act’s requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

Handwritten signature of Mark N. Melkerson in black ink, written over the typed name.

Mark N. Melkerson, M.S.

Director

Division of General, Restorative and

Neurological Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K060474

Device Name: Lisfranc Plates

### Indications For Use:

The NEWDEAL Lisfranc Plates are intended for fractures, fusions, osteotomies and replantations of small bones at the tarsometatarsal joints (Lisfranc joints)

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
Page 1 of 1  
\_\_\_\_\_  
(Division Sign-Off)  
Division of General, Restorative,  
and Neurological Devices

510(k) Number   K060474