

510(k) Summary
ConMed Linvatec Biomaterials
Matryx™ Interference Screw

MAR 02 2007

Submitter's Name, Address, Telephone Number, and Contact Person

ConMed Linvatec Biomaterials Ltd.
Päivi Mesiranta
Regulatory Specialist
P.O.Box 3
FIN-33721 Tampere
Finland, Europe
Phone: 358-3-316 5600
Facsimile: 358-3-316 5688

Date prepared: October 4th, 2006

Name of the device:

- A. Trade or Proprietary Name: Matryx™ Interference Screw
- B. Common Name: Bioabsorbable Interference Screw
- C. Classification Name: Bone Fixation Screw
- D. Device Product Code: MAI and HWC

Predicate Device:

1. ConMed Linvatec BioScrew (K933719, K952831, K973758), ConMed Linvatec Biomaterials Matryx™ Interference Screw (K052080)

Intended Use:

The Matryx™ Interference Screw is intended for use in interference fixation of bone – patellar tendon – bone and soft tissue grafts in anterior and posterior cruciate ligament reconstructions.

The Matryx™ Interference Screw is not intended for use in and is contraindicated for 1) Insufficient quality and quantity of bone for attachment of graft, 2) Blood supply limitations and/or previous infections, which could retard healing, 3) Foreign body sensitivity to the implant material. Where the material is suspected a test should be made prior to implantation to rule out sensitivity, 4) Patients with active sepsis or infection, 5) Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing and rehabilitation period, 6) ACL repairs which would not be appropriate for fixation with metallic screws.

Device Description:

The device description of the Matryx™ Interference Screw is as follows.

- The implant is composed of mixture of poly-96L/4D-lactide copolymer and tri-calcium phosphate. The predicate device previously cleared Matryx Interference Screw is made of the very same raw material.
- Lengths of implant are 25 - 38 mm
- Diameters of implant are 9 mm-11mm.

The only modifications that were made are:

- Amendment of a new screw sizes.
- Thread profile for longer 9 mm screw is similar with already cleared Matryx™ Interference Screws (K052080). The thread profile of diameter 10 mm and 11 mm screws is similar with predicate device ConMed Linvatec BioScrew (K933719, K952831, K973758)
- Amendment of Bone Taps for screws of 9 mm x 38 mm, 10 mm and 11 mm.
- Reference numbers for these new screw versions. These changes are updated in labelling.
- Revision of insert sheet covering the new screw sizes.

Substantial Equivalence:

ConMed Linvatec Biomaterials Ltd Matryx Interference Screw is substantially equivalent to the cleared predicate device. The applied modifications do not raise any new concerns of safety and efficacy of the implant.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

ConMed Linvatec Biomaterials Ltd
% Ms. Päivi Mesiranta
Regulatory Specialist
P.O. Box 3
FIN-33721 Tampere
Finland, Europe

MAR 02 2007

Re: K063588
Trade/Device Name: Matryx™ Interference Screw
Regulation Number: 21 CFR 888.3040
Regulation Name: Smooth or threaded metallic bone fixation fastener
Regulatory Class: Class II
Product Code: HWC
Dated: October 6, 2006
Received: January 16, 2007

Dear Ms. Mesiranta:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set

Page 2 – Ms. Päivi Mesiranta

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or 240-276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Mark N. Melkerson
Director
Division of General, Restorative and
Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATIONS FOR USE

510(K) Number (if known): K 063588

Device Name: Matryx™ Interference Screw

Indications for Use:

The Matryx™ Interference Screw is intended for use in interference fixation of bone – patellar tendon – bone and soft tissue grafts in anterior and posterior cruciate ligament reconstructions.

The Matryx™ Interference Screw is not intended for use in and is contraindicated for 1) Insufficient quality and quantity of bone for attachment of graft, 2) Blood supply limitations and/or previous infections, which could retard healing, 3) Foreign body sensitivity to the implant material. Where the material is suspected a test should be made prior to implantation to rule out sensitivity, 4) Patients with active sepsis or infection, 5) Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing and rehabilitation period, 6) ACL repairs which would not be appropriate for fixation with metallic screws

Prescription Use OR Over-The-Counter Use No
(Per 21 CFR 801.109)

(Please do not write below this line – continue on another page if needed)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

Division of General, Restorative,
and Neurological Devices

510(k) Number K 06 3588